

# AIR WALKER INSTRUCTIONS

## A

For first time application, loosen straps and remove liner from boot.



## B

Place foot in liner and secure with contact closure. Make sure the heel fits snugly into the posterior portion of the liner. Fasten the foot flaps on the liner first, making sure they are snug. Wrap and fasten the leg portion of the liner, again making sure that it is snug from bottom to top. Some trimming of the liner may be necessary to accomplish optimum fit.



## C

Spread the uprights using both hands and step into boot, aligning uprights with midline of the ankle.

Note: If uprights have plastic tube liner, remove them now to affix upright to liner.



## D

Secure boot straps, starting at the toes and working up the leg.



To **INFLATE** air chambers simply insert the **WHITE** tip of the air pump and squeeze pump. **DO NOT OVER INFLATE**. The air chamber should only be inflated to a comfortable level. To **DEFLATE** insert **BLACK** tip into the valve and squeeze pump until comfortable level is attained

Slight air loss may be experienced during periods of wear and air can easily be added.

Foam pads are provided to alleviate the problem of excess rubbing and pads can be placed either inside the liner or between the liner and the upright stay.

