





# Designed for life.

# Our Belief

We are committed to the highest standards of ethical conduct in all aspects of our business. We believe that honesty and integrity engender trust, which is the cornerstone of our business.

As a family owned and operated business, we take personal and professional pride in the products that we manufacture and distribute. We are proud of the fact that we have successfully cultivated a "family-friendly" working environment with generations of families contributing towards the steady growth of our company.

Our employees are respected, treated with fairness, and recognized as assets to our team. This principle philosophy allows each person to be consistently challenged in all aspects of their work, permitting them the opportunity to excel and advance professionally.

# Our Operation

Our success in the global marketplace is the result of our commitment to our manufacturing and distribution facility, which recently saw a 25% increase in space occupied. That's been our philosophy since we incorporated over twenty-five years ago, and it's why we will continue to grow in America and abroad. We thank our customers for the support and trust they've shown us, and look forward to challenging ourselves to create new products that bring value to our customers and the medical community.

# Our Intent

We will continue to make it possible for millions of people to resume every day activities; whether it is a return to work, resumption of a sports activity or simply to live a more active life. We feel growth isn't just a priority, but a necessity for our business. Corflex will achieve its goals by continued diversification, product development and greater market penetration. We are committed to not only keep pace with the health care changes and challenges, but to help influence them.

We intend to maintain the twenty-nine year tradition of nurturing our corporate culture, manufacturing a high quality product and achieving desired patient outcomes. Corflex – designed for life.







# How To Place An Order

The Corflex Customer Service Team is ready to assist you in placing your order. You may place your orders via a variety of methods:









TELEPHONE 800.426.7353, Ext. 1 603.623.3344, Ext. 1

603.623.4111

EMAIL service@corflex.com

Electronic Data Interchange For details, contact sales@corflex.com

EDI

# Terms & Conditions

Payment terms are Net 30 days. For the convenience of our customers Corflex also accepts MasterCard, Visa, Discover, American Express and wire transfers for payment of invoices. Orders are shipped F.O.B. Origin with shipping charges added to the invoice unless otherwise specified. Our minimum order is \$25; there will be a \$10 service charge on orders less than \$25.

# **Returns & Exchanges**

An authorization number must be obtained from Corflex prior to returning or exchanging any product. All returns and exchanges are subject to the following re-stocking fees: 1-45 days 15% fee; 46-90 days 25% fee; after 90 days products are non-returnable for credit or exchange. Special orders or custom sized products may not be returned for credit or exchange.

# Warranty

Corflex offers a limited lifetime warranty on the products it manufactures. Products manufactured by Corflex are warranted to the original purchaser against defects in workmanship and material for the life of the product. The foregoing warranty shall not apply to normal wear and tear; or if the product has been subjected to mis-use or application the warranty will be voided at the discretion of Corflex. Furthermore, Corflex offers no other warranty, express or implied, including but not limited to the implied warranties of merchantability and fitness of purpose.

# Symbol Key

The following symbols are used throughout the catalog to assist you in locating products to fit your needs.



# Cervical / Spine

Back pg. 3 Neck pg. 12



# **Upper Extremity**

Shoulder pg. 15 Elbow pg. 26 Wrist/Hand pg. 32



# Lower Extremity

Knee pg. 47
Foot/Ankle pg. 70
Hip pg. 80

# **Pediatric Products**

Cervical pg. 83
Shoulder pg. 84
Arm pg. 85
Wrist pg. 86
Knee pg. 87
Foot/Ankle pg. 88

### Index

Product Listing pg. 91







This symbol indicates that bariatric sizes are available for this product.




This symbol certifies that all products manufactured by Corflex are latex free.



This symbol indicates that certain precautions are needed when using this product.



This symbol identifies CoolTex products that are manufactured from a neoprene-alternative material.



# CERVICAL/ SPINE

Disc Unloader Spinal Orthosis - helping to reduce pain and get you back to the daily activities you enjoy

**SPINAL ORTHOSIS BACK SUPPORTS** BACK MISCELLANEOUS BINDERS/BELTS **CERVICAL SUPPORTS** 



# Disc Unloader Spinal Orthosis

Designed by a leading physiatrist, the Disc Unloader provides relief from chronic or acute back pain by unloading pressure off the lumbar disc and shifting the load to the abdominal region. This brace is designed with patient comfort in mind. Made from a breathable lightweight fabric, the modular design of this versatile spinal orthosis allows it to be used throughout the continuum of care to provide stability and support without interfering with daily activities.

- Unique modular design reduces intradiscal pressure through. circumferential compression, standard 6" anterior and 12" posterior panels accept anatomically designed inserts along with gel pack that provides the benefits of cold/hot therapy; available with extended 10" anterior panel for additional support of pendulous abdomen
- Indicated for chronic or acute low back pain, discogenic disc disease, post-surgical application, strains/sprains to lumbar muscles/ligaments and postural support
- Constructed of soft breathable fabric with contoured rigid A.B.S. inserts, elastic double-pulls and contact closure

### Product Number

	XS	S	M	L	XL
6" Standard	35-4000	35-4001	35-4002	35-4003	35-4004
10" Extended	35-4100	35-4101	35-4102	35-4103	35-4104
	2XL	3XL	4XL	GEL PACK	(
Standard	<b>2XL</b> 35-4005	<b>3XL</b> 35-4006	<b>4XL</b> 35-4007	<b>GEL PACK</b> 72-7006	(

### Sizing Guide

XS	S	М	L	XL
25"- 30"	30"- 35"	35"- 40"	40"- 45"	45"- 50"
2XL	3XL	4XL	GEL PAC	K
50"- 55"	55"- 60"	60"- 65"	9"x 11"	









For additional sizing, laundering or application instructions visit corflex.com/discunloader



6" Standard Anterior Panel



10" Extended Anterior Panel

# **Disc Unloader LSO**

- → Offers support from T9 to SC junction; features a modular design that helps relieve back pain by reducing intradiscal pressure through circumferential compression; 6" anterior and 15" posterior panels accept anatomically designed inserts along with a gel pack that provides the benefits of cold/hot therapy; also available with extended 10" anterior panel for additional support of pendulous abdomen
- Indicated for chronic or acute back pain, discogenic disc disease, post-surgical stabilization, spinal stenosis, herniated disc, stable non-displaced spinal fractures (L2-L5), radiculopathy, sciatica, spondylosis, spondylolisthesis, facet syndrome and general weakness of the lumboscaral spine
- Constructed of soft breathable fabric with contoured rigid A.B.S. anterior and posterior plastic inserts, elastic double pulls and contact closures

Product Number

	XS	S	М	L	XL	2XL	3XL	4XL
6" Standard	35-4200	35-4201	35-4202	35-4203	35-4204	35-4205	35-4206	35-4207
10" Extended	35-4300	35-4301	35-4302	35-4303	35-4304	35-4305	35-4306	35-4307
	GEL PAC	CK						
	72-7006							

Sizing Guide

oizing darac								
	XS	S	M	L	XL	2XL	3XL	4XL
	25"- 30"	30"- 35"	35"- 40"	40"- 45"	45"- 50"	50"- 55"	55"- 60"	60"- 65"
	<b>GEL PAC</b>	K						

9"x 11"

Measurement taken from waist circumference

See page 94 for Back Support Comparison Chart.

# ALL PRODUCTS MANUFACTURED BY CORFLEX ARE LATEX FREE.





# **E/N Lumbar Support**

- Features tapered elastic side panels with neoprene posterior back panel and pocket that accepts either a thermoplastic insert, gel pack, or pressure pad (sold separately); also features double-pull tension straps for added compression
- Indicated for moderate to severe low back pain, lumbar muscle weakness, lumbar sprains/strains, and lumbar instability
- Manufactured from premium quality elastic with contact closure



Product Number	
----------------	--

	XS	S	M	L
Beige	77-6220	77-6221	77-6222	77-6223
Black	77-6120	77-6121	77-6122	77-6123
	XL	2XL	3XL	4XL
Beige			<b>3XL</b> 77-6226	.,

# Sizing Guide

0121119 00	11010			
	XS	S	M	L
	18"- 24"	24"- 30"	30"- 36"	36"- 42"
	XL	2XL	3XL	4XL
	12"_ 18"	/8" <sub>-</sub> 52"	52"- 56"	56"- 60"

Measurement taken from hip circumference

Insert options: 37-6031 Thermoplastic Insert 77-9999 Pressure Pad 72-7006 Gel Cold/Hot Pack



# **Back-2-Form Support**

- Features tapered elastic side panel with elastic back panel and pocket that holds either a thermoplastic insert, gel pack, or pressure pad (sold separately); double-pull tension straps for added compression or web strap for static treatment
- Indicated for moderate to severe low back pain, lumbar muscle weakness, lumbar sprains/strains, and lumbar instability
- Manufactured from premium quality elastic with contact closure



### Product Number

	XS	S	М	L
	35-1840	35-1841	35-1842	35-1843
w/Strap	35-1850	35-1851	35-1852	35-1853
	XL	2XL	3XL	4XL
			<b>3XL</b> 35-1846	

### Sizing Guide

XS	5	M	L
18"- 24"	24"- 30"	30"- 36"	36"- 42"
XL	2XL	3XL	4XL

Measurement taken from hip circumference

Insert options: 37-6031 Thermoplastic Insert 77-9999 Pressure Pad 72-7006 Gel Cold/Hot Pack



# **Contour Back Support**

- Features contoured front panel to accommodate hip area and posterior pocket that accepts either a thermoplastic insert, gel pack, or pressure pad (sold separately); also features double-pull tension straps for added compression
- Indicated for moderate to severe low back pain, lumbar muscle weakness, lumbar sprains/strains, and lumbar instability
- Manufactured from premium quality elastic with contact closure



### Product Number

XS	S	M	L
35-1740	35-1741	35-1742	35-1743
XL	2XL	3XL	
25 17//	35-1745	35 17/16	

# Sizing Guide

JIZIII I G	ilde			
	XS	S	M	L
	18"- 24"	24"- 30"	30"- 36"	36"- 42"
	XL	2XL	3XL	

Measurement taken from hip circumference

Insert options: 37-6031 Thermoplastic Insert 77-9999 Pressure Pad 72-7006 Gel Cold/Hot Pack





# **Ultra Lumbo Sacral Support**

- Features soft breathable back panel with four malleable stays and double-pull tension straps for added compression
- Indicated for moderate low back pain, lumbar muscle weakness, lumbar sprains/strains
- Manufactured from premium quality elastic with malleable aluminum stays and contact closure




### Product Number

XS	S	M	L	XL
35-1500	35-1501	35-1502	35-1503	35-1504
2XL	3XL	4XL		

### Sizing Guide

XS	S	M	L	XL
18"- 24"	24"- 30"	30"- 36"	36"- 42"	42"- 48"
2XL	3XL	4XL		
48"- 52"	52"- 56"	56"- 60"		

Measurement taken from hip circumference



# **Criss-Cross Back Support**

- Features criss-cross back panel with flexible anti-roll stays and double-pull elastic tension straps for added compression
- Indicated for moderate low back pain, lumbar muscle weakness and lumbar sprains/strains
- Manufactured from premium quality elastic with contact closure



### Product Number

\/C				\/I
XS	S	М	L	XL
35-1600	35-1601	35-1602	35-1603	35-1604
2XL	3XL	4XL		
35-1605	35-1606	35-1607		

### Sizing Guide

XS	S	M	L	XL
18"- 24"	24"- 30"	30"- 36"	36"- 42"	42"- 48"
2XL	3XL	4XL		
48"- 52"	52"- 56"	56"- 60"		

Measurement taken from hip circumference



# **Target Back Wrap**

- Provides support, compression and warmth for the back, with sewn-in tapered pressure pad; universal design allows for inventory reduction
- Indicated for mild low back pain, lumbar sprains/strains or weakness
- Manufactured from premium latex-free neoprene



### Product Number

UNIV	2XL/3XL
88-1800	88-1801

### Sizing Guide

UNIV	2XL/3XL	
24"- 46"	46"- 56"	

Measurement taken from hip circumference

See page 94 for Back Support Comparison Chart

# Thoracolumbar Support

- Features convenient wrap-around design with suspenders for ease of application and double-pull tension straps for added compression, posterior pocket accepts thermoplastic insert (sold separately)
- Indicated for pain in the thoracolumbar region and compression
- Manufactured from premium quality elastic with adjustable suspenders and contact closure

### Product Number

XS	S	M	L
35-1700	35-1701	35-1702	35-1703
XL	2XL		

### Sizina Guide

0				
	XS	S	M	L
	18"- 24"	24"- 30"	30"- 36"	36"- 42"
	XL	2XL		
	42"- 48"	48"- 52"		
			,	

Measurement taken from hip circumference Insert option:

37-6032 TLSO Thermoplastic Insert



# V-Lock Back Support

- Features double elastic v-shaped pull for easy application with contoured front for optimal patient comfort; also includes plastic and optional aluminum stays for flexible or firm support
- Indicated for moderate low back pain, lumbar muscle weakness, lumbar sprains/strains
- Manufactured from premium quality elastic with contact closure; includes both plastic and aluminum stays



### Product Number

XS	S	M	L
35-1760	35-1761	35-1762	35-1763
XL	2XL	3XL	4XL

Sizing Gu	ıide			
	XS	S	M	L
	18"- 24"	24"- 30"	30"- 36"	36"- 42"
	XL	2XL	3XL	4XL
	42"- 48"	48"- 52"	52"- 56"	56"- 60"

Measurement taken from hip circumference



# V-Lock Lumbar Orthosis w/Rigid Insert

- Features pocket that accepts a preformed posterior insert and double elastic v-shaped pull for easy application with contoured front for optimal patient comfort; also includes plastic and optional aluminum stays for flexible or firm support
- Indicated for moderate to severe low back pain, lumbar muscle weakness, lumbar sprains/strains, and lumbar instability
- Manufactured from premium quality elastic with contact closure; includes rigid insert and both plastic and aluminum stays

### Product Number

M	L	XL	2XL
35-1792	35-1793	35-1794	35-1795
3XL	4XL		

# Sizing Guide

М	L	XL	2XL
30"- 36"	36"- 42"	42"- 48"	48"- 52"
3XL	4XL		
52"- 56"	56"- 60"		

Measurement taken from hip circumference







# **Industrial Back Support**

- Criss-cross design limits stress on the back and encourages proper lifting techniques; suspenders optional
- Indicated for lumbar muscle weakness, lumbar sprains/strains, moderate low back pain or prevention of back strain
- Manufactured from premium quality elastic with contact closure



### Product Number

	XS	S	M	L
Straps	77-8140	77-8141	77-8142	77-8143
No Straps	77-9150	77-9151	77-9152	77-9153
	XL	2XL	3XL	4XL
Straps	77-8144	77-8145	77-8146	77-8147
No Straps	77-9154	77-9155	77-9156	77-9157
Note: Also available in Beige and White				

### Sizing Guide

XS	S	M	L
18"- 24"	24"- 30"	30"- 36"	36"- 42"
XL	2XL	3XL	4XL

Measurement taken from hip circumference



# Low Profile Industrial Back Support

- Comfortable low profile design limits stress on back, encourages proper lifting techniques; suspenders optional
- Indicated for lumbar muscle weakness, lumbar sprains/strains, moderate low back pain or prevention of back strain
- Manufactured from premium quality elastic w/contact closure



### Product Number

	XS	S	M	L
Straps	77-8840	77-8841	77-8842	77-8843
No Straps	77-9840	77-9841	77-9842	77-9843
	XL	2XL	3XL	4XL
Straps	77-8844	77-8845	77-8846	77-8847
No Straps	77-9844	77-9845	77-9846	77-9847
Note: Also	available in	White		

### Sizing Guide

XS	S	M	L
18"- 24"	24"- 30"	30"- 36"	36"- 42"
XL	2XL	3XL	4XL

Measurement taken from hip circumference



# **Cryotherm Back Wrap**

- Economical wrap designed to provide compression and concentrated cold/ hot therapy for the back; available with optional stays
- Indicated for post-op use, arthritis, sprains/strains, and for integration into rehab program
- Manufactured from a soft breathable foam laminate with contact closure and pockets to accommodate reusable gel packs



### Product Number

	0 GELS	2 GELS	4 GELS
	72-7040	72-7042	72-7044
w/Stays	72-7040-300	72-7042-300	72-7044-300

Sizing Guide

Fits up to 48" waist circumference

Additional gels sold separately: 72-7001 - Single 5 1/2"x 11" Gel 72-7000 - Case/10 5 1/2"x 11" Gel



# Medic-Air Back Pillō

- Easily attaches to any chair to convert it into ergonomically correct seating; easily inflates to desired comfort level or deflates for portability
- Provides support for the back and relief from lower back pain
- Manufactured from medical grade nylon-coated vinyl with open-valve system

Product Number

	SINGLE	CASE/12
Blue	84-2051-100	84-2051-300
Grey	84-2054-100	84-2054-300

Sizing Guide

Product measures 18" x 15"



# Medic-Air Lumbar Röll

- Easily attaches to any chair to provide support for the lower back; easily inflates to desired comfort level or deflates for portability
- Supports posture and maintains lordotic curve while seated
- Manufactured from medical grade nylon-coated vinyl with open-valve system

Product Number

	SINGLE	CASE/12
Blue	84-3051-100	84-3051-300
Grey	84-3054-100	84-3054-300

Sizing Guide

Product measures 20" x 8"



# **Medic-Air Seat Cushion**

- Provides comfortable, stress-free sitting; easily inflates or deflates for portability
- Helps eliminate pressure points, suitable for sporting events and short-term wheelchair use
- Manufactured from medical grade nylon-coated vinyl with open-valve system

Product Number

	SINGLE	CASE/12
Black	84-5052-100	84-5052-300

Sizing Guide

Product measures 18" x 16"





# **Universal Rib Belt**

- Anatomically designed (male or female) to provide optimal comfort and support; universal design allows for inventory reduction
- Indicated for fractured or bruised ribs, intercostal muscle strains and costovertebral sprains
- Manufactured from premium quality elastic with contact closure

Proc	luct	Nu	m	ber
------	------	----	---	-----

	UNIV
Male	31-1951
Female	31-1952

# Sizing Guide

9	Garac
	UNIV
	24"_ 48

Measurement taken from chest circumference



# **Elastic Rib Belt**

- Anatomically designed (male or female) to provide optimal comfort and support
- Indicated for fractured or bruised ribs, intercostal muscle strains and costovertebral sprains
- Manufactured from premium quality elastic with contact closure

### Product Number

	S	М	L	
Male	31-1901	31-1902	31-1903	
Female	31-1911	31-1912	31-1913	
	XL	2XL		
Male	<b>XL</b> 31-1904	_,		

### Sizing Guide

_			
	S	M	L
	24"- 30"	30"- 36"	36"- 42"
	XL	2XL	
	42"- 48"	48"- 54"	

Measurement taken from chest circumference



# **Trochanter Belt**

- Designed to provide support and stability over the sacral region
- Indicated for instability of the sacroiliac joints and the lumbosacral joint
- Constructed of soft pile laminated to flannel lining with elastic side pulls and contact closure

### Product Number

UNIV
35-1419

### Sizing Guide

### UNIV

Fits up to 48" hip circumference

# **Universal Abdominal Binder**

- Features solid panel elastic for uniform support of entire abdominal region; universal design allows for inventory reduction
- Indicated for abdominal muscle strain or post-op application
- Constructed of plush elastic with soft foam laminate and contact closure

Product Number

WIDTH	UNIV
9"	32-1950
12"	32-1951

Sizing Guide

UNI	V
24"-	48

Measurement taken from chest circumference



# Panel Elastic Abdominal Binder

- Features panel elastic construction with faggoted stitching between panels to prevent rolling; available in 9" 3-panel or 12" 4-panel width
- Indicated for abdominal muscle strain or post-op application
- Manufactured from premium panel elastic with contact closure



### Product Number

WIDTH	S	М	L	XL
9"	32-1961	32-1962	32-1963	32-1964
12"	32-1971	32-1972	32-1973	32-1974
WIDTH	2XL	3XL	4XL	
WIDTH 9"	_,	<b>3XL</b> 32-1966	.,,-	

## Sizing Guide

S	M	L	XL
24"- 30"	30"- 36"	36"- 42"	42"- 48"
2XL	3XL	4XL	

Measurement taken from chest circumference



# **Ultra Abdominal Binder**

- Provides even compression around the abdomen to support muscles and help minimize pain; features hook compatible material that offers easy adjustability and secure fit; soft lining for maximum patient comfort
- Indicated after abdominoplasty, cesarean section, hernia repair, laparotomy, liposuction, bladder surgery or following bariatric surgical treatments
- Manufactured from hook compatible plush elastic material

### Product Number

	WIDTH	S	M	L
	9"	32-1921	32-1922	32-1923
/	12"	32-1931	32-1932	32-1933
	WIDTH	XL	2XL	
	9"	32-1924		
	12"	32-1934	32-1935	

### Sizing Guide

S	M	L
30"- 45"	46"- 62"	60"- 75"
XL	2XL	

Measurement taken from chest circumference







# Philadelphia Cervical Collar w/Trachea

- Features two-piece design anatomically contoured to provide superior stability and comfort with large trachea opening for easy access
- Indicated for immobilization for both extrication and rehabilitation, post-operative discectomy and spinal fusion
- Constructed of durable plastic laminated to soft hypoallergenic foam with contact closure



### Product Number

	PEDI	S	M	L
2 1/4"	19-4100	19-4121	19-4122	19-4123
3 1/4"		19-4131	19-4132	19-4133
4 1/4"		19-4141	19-4142	19-4143
5 1/4"		19-4151	19-4152	19-4153

### Sizing Guide

PEDI	S	М	L
8"- 11"	10"- 13"	13"- 16"	16"- 19"

Measurement taken from circumference of the neck at "Adam's apple"

Height measurement taken from posterior portion of mandible to the clavicle directly below



# **Ambu Perfit Collar**

- Features one-piece flat design with flip out chin piece and trachea opening for easy access
- Indicated for immobilization for extrication and rehabilitation
- Constructed of durable plastic laminated to foam with contact closure





Ace Adjustable

### Product Number

INFANT	PEDI	NECKLESS	SHORT
19-9001	19-9002	19-9003	19-9004
<b>REGULAR</b>	TALL	ACE ADJ	MINI ADJ
19-9005	19-9006	19-9007	19-9008

### Sizing Guide

INFANT	PEDI	NECKLESS	SHORT
Purple	Aqua	Red	Yellow
REGULAR	TALL	ACE ADJ	MINI ADJ
Blue	Pink	Univ Adult	Univ Pediatric



# **Ultra Cervical Collar**

- Features anatomically contoured foam with soft cotton stockinette covering for maximum patient comfort
- Indicated for cervical facet syndrome, discogenic pain, degenerated disc disease and muscle sprains or strains of cervical region
- Manufactured from medium firm density foam with 2" contact closure



# Product Number

	2XS	XS			
2"	15-5000	15-5001			
	UNV	S	М	L	XL
$2^{1}/_{2}$ "	15-5020	15-5021	15-5022	15-5023	15-5024
$3^{1}/_{2}"$	15-5030	15-5031	15-5032	15-5033	15-5034
$4^{1}/_{2}"$	15-5040	15-5041	15-5042	15-5043	15-5044

### Sizing Guide

0				
2XS	XS			
10"- 12"	14"- 16"			
UNV	S	M	L	XL
16"- 24"	16"- 20"	18"- 22"	20"- 24"	22"- 26"

Measurement taken from circumference of the neck at "Adam's apple"

Height measurement taken from posterior portion of mandible to the clavicle directly below

# **Cryotherm Neck Wrap**

- Economical wrap designed to provide concentrated cold/hot therapy for the cervical spine
- Indicated for post-op use, arthritis, strains or sprains and for integration into rehab program
- Manufactured from a soft breathable foam laminate with contact closure and pocket to hold reusable gel pack



Product Number

	0 GELS	1 GEL	2 GELS
31/2"	72-7120	72-7121	72-7122
41/2"	72-7130	72-7131	72-7132

Sizing Guide

Fits up to a 22" neck circumference

Additional gels sold separately: 72-7003 - Single 3 <sup>1</sup>/<sub>2</sub>"x 11" Gel 72-7002 - Case/10 3 <sup>1</sup>/<sub>2</sub>"x 11" Gel



# Medic-Air Cervical Sleep Pillō

- Positions between pillow and pillowcase to provide the gentle and adjustable support of air for the cervical spine; easily inflates to desired comfort level and deflates for portability
- Provides cervical support while sleeping
- Manufactured from medical grade durable nylon-coated vinyl with open-valve system

Product Number

	SINGLE	CASE/12
Blue	84-6051-100	84-6051-300

Sizing Guide

Product measures 21"x 4"



# Medic-Air Snooze Pillō

- Provides the gentle and adjustable support of air for the head and cervical spine; easily inflates to desired comfort level and deflates for portability
- Permits comfortable sleeping while in a sitting position
- Manufactured from medical grade durable nylon coated vinyl with open valve system

Product Number

	SINGLE	CASE/12
Blue	84-4051-100	84-4051-300

Sizing Guide

Product measures 14"x 9"x 4"





# UPPER EXTREMITY

Tennis Elbow Sleeve providing the support you need for your active lifestyle

SHOULDER IMMOBILIZERS

SHOULDER MISCELLANEOUS

**HUMERAL SPLINTS** 

**ARM SLINGS** 

**ELBOW SPLINTS** 

**ELBOW SLEEVES** 

TENNIS ELBOW SUPPORTS

**WRIST SPLINTS** 

**THUMB SPLINTS** 

**FINGER SPLINTS** 

HAND/FOREARM SPLINTS

15



# Shoulder Abduction Pillow w/Firm Fit Sling

Developed in association with sports medicine leaders, our Shoulder Abduction Pillow with Firm Fit Sling achieves three objectives: secure immobilization for post operative use, ease of application for both practitioner and patient, maximum comfort and wearability for long or short term use. The contoured sling features moisture wicking terry cloth material and the pillow allows for 15° of abduction.

- Features easy to apply design with terry cloth sling for optimal patient comfort
- Indicated for treatment after Bankart repairs, rotator cuff repairs, shoulder arthroscopy, separations and dislocations
- Constructed of firm density foam covered in tricot material with terry cloth sling and exercise ball to stimulate circulation


Product Number

S	М	L	XL
23-1931	23-1932	23-1933	23-1934

Sizing Guide

	S	M	L	XL
Hip	24"- 30"	30"- 36"	36"- 42"	42"- 48"
Forearm	9"- 12"	12"- 15"	15"- 18"	18"- 20"

Measurement taken from hip circumference and length of forearm from wrist to elbow





For additional sizing, laundering or application instructions visit corflex.com/shoulderabpillow



# Ultra Shoulder Abduction Pillow w/Firm Fit Sling

- → Features extended pillow design with malleable stay for proper positioning; V-lock strapping system helps prevent roll / shifting; padded neck strap for optimal patient comfort
- Indicated for treatment after Bankart repairs, rotator cuff repairs, shoulder arthroscopy, separations and dislocations
- Constructed of firm density foam covered in tricot material with terry cloth sling and exercise ball to stimulate circulation

Product Number

S	M	L	XL
23-1981	23-1982	23-1983	23-1984

Sizing Guide

	S	М	L	XL
Hip	24"- 30"	30"- 36"	36"- 42"	42"- 48"
Forearm	9"- 12"	12"- 15"	15"- 18"	18"- 20"

Measurement taken from hip circumference and length of forearm from wrist to elbow



# Ultra Shoulder Abduction Pillow w/Sling

- Features extended pillow design with malleable stay for proper positioning, V-lock strapping system to prevent roll / shifting and padded neck strap for optimal patient comfort
- Indicated for treatment after Bankart repairs, rotator cuff repairs, shoulder arthroscopy, separations and dislocations
- Constructed of firm density foam covered in tricot material with tricot sling and exercise ball to stimulate circulation

Product Number

S	M	L	XL
23-1901	23-1902	23-1903	23-1904

Sizing Guide

	S	M	L	XL
Hip	24"- 30"	30"- 36"	36"- 42"	42"- 48"
Forearm	9"- 12"	12"- 15"	15"- 18"	18"- 20"

Measurement taken from hip circumference and length of forearm from wrist to elbow



# Shoulder Abduction Pillow w/Sling

- Features pillow with malleable stay for proper positioning and padded neck strap for optimal patient comfort
- Indicated for treatment after Bankart repairs, rotator cuff repairs, shoulder arthroscopy, separations and dislocations
- Constructed of firm density foam covered in tricot material with tricot sling and exercise ball to stimulate circulation

Product Number

S	M	L	XL
23-1911	23-1912	23-1913	23-1914

Sizing Guide

	S	M	L	XL
Hip	24"- 30"	30"- 36"	36"- 42"	42"- 48"
Forearm	9"- 12"	12"- 15"	15"- 18"	18"- 20"

Measurement taken from hip circumference and length of forearm from wrist to elbow

# **ER Shoulder Abduction Pillow w/Sling**

- Features removable wedge that provides 30° of external rotation, V-lock strapping system to prevent roll / shifting and padded neck strap for optimal patient comfort
- Indicated for protocols involving non-surgical immobilization in external rotation to help reduce the incidence of recurrent dislocation
- Constructed of firm density foam covered in tricot material with tricot sling and exercise ball to stimulate circulation

Product Number

S	M	L	XL
23-1941	23-1942	23-1943	23-1944

Sizing Guide

	S	M	L	XL
Hip	24"- 30"	30"- 36"	36"- 42"	42"- 48"
Forearm	9"- 12"	12"- 15"	15"- 18"	18"- 20"

Measurement taken from hip circumference and length of forearm from wrist to elbow



# Neutral Shoulder Pillow w/Sling

- Features contoured pillow that maintains the shoulder in a natural state of external rotation and padded neck strap for optimal patient comfort
- Indicated for treatment after SLAP lesion repairs, Bankart repairs, rotator cuff repairs, separations and dislocations
- Constructed of firm density foam covered in tricot material with 3-D spacer fabric sling and exercise ball to stimulate circulation

Product Number

S	M	L	XL
23-1921	23-1922	23-1923	23-1924

Sizing Guide

	S	M	L	XL
Hip	24"- 30"	30"- 36"	36"- 42"	42"- 48"
Forearm	9"- 12"	12"- 15"	15"- 18"	18"- 20"

Measurement taken from hip circumference and length of forearm from wrist to elbow





# **Elastic Shoulder Immobilizer**

- Features 4" wide elastic chest band and humerus and forearm cuff with contact closure
- Indicated for immobilization following rotator cuff repairs, S/P glenohumeral acromioplasty surgery or other shoulder-related injuries
- Constructed of premium quality elastic with vinyl cuffs



### Product Number

XS	S	M	L
23-5240	23-5241	23-5242	23-5243
XL	2XL	3XL	

### Sizing Guide

XS	S	M	L
18"- 24"	24"- 30"	30"- 36"	36"- 42"
XL	2XL	3XL	
42"- 48"	48"- 54"	54"- 60"	

Measurements taken from chest circumference



# **Universal Comfort Shoulder Immobilizer**

- Features 6" wide tricot chest band and padded humerus and forearm cuff for optimal patient comfort; universal design allows for inventory reduction
- Indicated for immobilization following rotator cuff repairs, S/P glenohumeral acromioplasty surgery or other shoulder-related injuries
- Manufactured from tricot material with polyester kodel pile liner on cuffs with contact closure

Product Number

**UNIV** 23-4100

Sizing Guide

### UNIV

Fits 22" - 45" chest circumference



# **Universal Sling & Swathe**

- Features cradle style sling and swathe with contact closure
- Indicated for immobilization following rotator cuff repairs, S/P glenohumeral acromioplasty surgery or other shoulder-related injuries
- Manufactured from seamless, lightweight, durable material

Product Number

UNIV 23-4204

Sizing Guide

**UNIV** 24"- 48"

Measurement taken from chest circumference

# **Ultra Sling & Swathe**

- > Features high profile sling designed to encompass the elbow and swathe with contact closure
- Indicated for immobilization following rotator cuff repairs, S/P glenohumeral acromioplasty surgery, or other shoulder related injuries
- Manufactured from seamless, lightweight, durable material

Product Number

S	М	L	XL
23-4401	23-4402	23-4403	23-4404

Sizing Guide

	S	М	L	XL
Sling Pouch	9"- 10"	13"- 10"	16"- 11"	19"- 12"
Chest Cir	24"- 36"	30"- 42"	36"- 48"	42"- 54"



# Sling & Swathe

- Features cradle-style sling and swathe with contact closure
- Indicated for immobilization following rotator cuff repairs, S/P glenohumeral acromioplasty surgery or other shoulder-related injuries
- Manufactured from seamless durable material

Product Number

S

	•		_	/ -
	23-1746	23-1747	23-1748	23-1749

Sizing Guide

	S	M	L	XL
Sling Pouch	9"- 10"	13"- 10"	16"- 11"	19"- 12"
Chest Cir	24"- 36"	30"- 42"	36"- 48"	42"- 54"





This symbol indicates that pediatric sizes are available for this product.



This symbol indicates that bariatric sizes are available for this product.



This symbol certifies that all products manufactured by Corflex are latex free.



This symbol indicates that certain precautions are needed when using this product.



This symbol identifies CoolTex products that are manufactured from a neoprene-alternative material.



# Ultra Shoulder Immobilizer

- > Features 1½" wide waist strap for secure immobilization and thumb loop to help prevent wrist drop
- Indicated for immobilization following rotator cuff repairs, S/P glenohumeral acromioplasty surgery or other shoulder-related injuries
- Manufactured from seamless, lightweight, durable material with added pad for shoulder strap



# Product Number

	2XS	XS	S	M
	23-1869	23-1870	23-1871	23-1872
	L	XL		
	23-1873	23-1874		
Sizing Guide				
	avc	VC	c	
	2XS	XS	S	M
Width/Height			•	•••
Width/Height			•	•••

Measurements are sling pouch width x height



# **Economizer Shoulder Immobilizer**

- > Features 1" shoulder and waist straps with slide buckle
- Indicated for immobilization following rotator cuff repairs, S/P glenohumeral acromioplasty surgery or other shoulder-related injuries
- Manufactured from seamless lightweight muslin material



### Product Number

	XS	S	M
	23-1740	23-1741	23-1742
	L	XL	
	23-1743	23-1744	
Sizing Guide			
	XS	S	M
Width/Height	12"x 5"	15"x 6"	17"x 7"
	L	XL	
Width/Height	19"x 8"	21"x 9"	

Measurements are sling pouch width x height



# **Universal Shoulder Immobilizer**

- > Features 8" wide chest band with anti-migration strap and 21/2" wide moveable humerus and wrist cuffs for ease of application; universal design allows for inventory reduction
- Indicated for immobilization following rotator cuff repairs, S/P glenohumeral acromioplasty surgery or other shoulder-related injuries
- Constructed of foam laminate to hook sensitive fabric

### Product Number

# UNIV 23-1750 Sizing Guide

Fits 24" - 48" chest circumference

# **Cryo Pneumatic Shoulder**

- Anatomical shoulder wrap designed to combine the benefits of pneumatic compression with cold therapy; features pump bulb that detaches once desired compression is achieved
- Indicated for post-op application and as part of a physical therapy or athletic training program; helps diminish swelling/edema, decrease pain and reduce muscle spasms
- Manufactured from plush foam laminate with contact closure straps, removable gel pack with tricot lining and non-toxic gel; additional universal extension straps available



	0 GELS	1 GEL	2 GELS	
	72-7180	72-7181	72-7182	
			UNIV	
Shoulder	Gel		72-7189	
Extension	n Straps (P	kg/10)	72-7190	

Sizing Guide

UNIV

Fits up to 48" chest circumference





# **Cryotherm Shoulder Wrap**

- Designed to contour to shoulder with pockets that hold gel packs to provide benefits of cold/hot therapy
- Indicated for post-op use, treatment of arthritis or as part of a physical therapy or rehab program
- Manufactured from soft breathable foam laminate with contact closure and pockets to hold reusable gel packs



Product Number

# **0 GELS 2 GELS 4 GELS** 72-7030 72-7032 72-7034

Sizing Guide

### UNIV

Fits up to 48" chest circumference

Additional gels sold separately: 72-7001 - Single 5 <sup>1</sup>/<sub>2</sub>"x 11" Gel 72-7000 - Case/10 5 <sup>1</sup>/<sub>2</sub>"x 11" Gel



# **Tricot Shoulder Immobilizer**

- Features waist strap for secure fit; thumb loop helps prevent wrist drop;
   Velcro closure at top of sling envelope
- Indicated for immobilization following rotator cuff repairs, S/P glenohumeral acromioplasty surgery or other shoulder-related injuries
- Manufactured from soft, durable tricot material



### Product Number

2XS

	23-5059	23-5060	23-5061	23-5062
	L	XL		
	23-5063	23-5064		
Sizing Guide				
	2XS	XS	S	M
		<b>XS</b> 12"x 5"	•	•••
			•	•••

XS







# Mid Length Humeral Splint

We designed this brace to increase patient comfort and save the practitioner valuable application time. Based on comments about frequently having to trim the distal end of humeral braces, the mid length was born. Corflex also offers extended and standard length braces. All three are designed to control humeral segments through soft tissue compression.

- Controls humeral segments through circumferential soft tissue compression and features deltoid cap extension to prevent migration; shorter length eliminates need to trim distal end
- Indicated for treatment of humeral diaphyseal fractures or deep humeral bone bruises
- Constructed of rigid polyethylene laminated with closed cell foam lining and contact closure; includes two stockinette undersleeves and one arm sling

Product Number

	XS	S	М	L	XL	2XL
	37-2180	37-2181	37-2182	37-2183	37-2184	37-2185
Sizing Gu	uide					

11"- 13" 13"- 15" 15"- 17" 17"- 19" Proximal 7"-9" 9"- 11" 9"- 11" 11"- 13" 13"- 15" 15"- 17"

Measurement taken at widest part of bicep and above elbow





For additional sizing, laundering or application instructions visit corflex.com/midlengthhumeral

# **Extended Length Humeral Splint**

- Controls humeral segments through circumferential soft tissue compression and features deltoid cap extension to prevent migration
- Indicated for treatment of humeral diaphyseal fractures or deep humeral bone bruises
- Constructed of rigid polyethylene laminated with closed cell foam lining and contact closure; includes two stockinette undersleeves and one arm sling

Product N	Jumber
-----------	--------

XS	S	M	L
37-2160	37-2161	37-2162	37-2163
XL	2XL		
ΛL	ZAL		

### Sizing Guide

	XS	S	M	L
Proximal	7"- 9"	9"- 11"	11"- 13"	13"- 15"
Distal	6"-7"	7"- 9"	9"- 11"	11"- 13"
	XL	2XL		
Proximal	<b>XL</b> 15"- 17"			

Measurement taken at widest part of bicep and above elbow



# **Humeral Splint**

- Provides circumferential soft tissue compression of the humerus while allowing full range of motion at shoulder and elbow
- Indicated for treatment of distal or mid shaft humeral fractures
- Constructed of rigid polyethylene laminated with closed cell foam lining and contact closure

### Product Number

5	M	L	XL
37-2061	37-2062	37-2063	37-2064

### Sizina Guide

	S	М	L	XL
Proximal	9"- 11"	11"- 13"	13"- 15"	15"- 17"
Distal	7"- 9"	9"- 11"	11"- 13"	13"- 15"

Measurement taken at widest part of bicep and above elbow











# **Ultra Arm Sling**

- > Features double d-ring hook and loop closure for ease of application and thumb loop to help prevent wrist drop
- Indicated for injuries to the arm, shoulder, wrist or hand and cast support to maintain arm in neutral position
- Manufactured from seamless, lightweight, durable material with added shoulder strap pad



Allows for slide buckle or contact closure for ease of application

**Economizer Arm Sling** 

- Indicated for injuries to the arm, shoulder, wrist or hand and cast support to maintain arm in neutral position
- Manufactured from seamless lightweight muslin material



# Mesh Arm Sling

- Allows for slide buckle or contact closure
- Indicated for injuries to the arm, shoulder, wrist or hand and cast support to maintain arm in neutral position
- Manufactured from seamless breathable mesh material



# Product Number

XS	S	М	L
21-1870	21-1871	21-1872	21-1873
XL	2XL		

# Sizing Guide

	XS	S	М	L	
Width/Height	12"x 5"	15"x 6"	17"×7"	19"x 8"	
	XL	2XL			
Width/Height	21"x 9"	23"x 10"			
Measurements are sling pouch width x height					

Product Number

2XS	XS	S	M
21-1969	21-1970	21-1971	21-1972
L	XL		
21-1973	21-1974		

# Sizina Guide

	2XS	XS	S	M	
Width/Height	10"x 4"	12"x 5"	15"x 6"	17"×7"	
	L	XL			
Width/Height	19"x 8"	21"x 9"			
Measurements are sling pouch width x height					

# Product Number

XS	S	M	L
21-1730	21-1731	21-1732	21-1733
XL			
21-1734			

Sizing Guide				
	XS	S	M	L
Width/Height	12"x 5"	15"x 6"	17"x 7"	19"x 8"
	XL			
Width/Height	21"x 9"			

Measurements are sling pouch width x height

# Tricot Arm Sling

- > Features thumb loop to prevent wrist drop and Velcro closure at top of sling envelope to provide a secure fit
- Indicated for injuries to the arm, shoulder, wrist or hand and cast support to maintain arm in neutral position
- Manufactured from soft yet durable tricot material



### Product Number

2XS	XS	S	M
21-5059	21-5060	21-5061	21-5062
1	XL		
_	/\_		

### Sizina Guide

Sizing Garac				
	2XS	XS	S	M
Width/Height	10"x 4"	12"x 5"	15"x 6"	17"×7"
	L	XL		
Width/Height	19"x 8"	21"x 9"		

Measurements are sling pouch width x height

# **Cradle Arm Sling**

- Open-end design allows arm to easily slide into sling
- Indicated for injuries to the arm, shoulder, wrist or hand and cast support to maintain arm in neutral position
- Manufactured from seamless durable material

Product Number

**S M L** 21-1976 21-1977 21-1978

Sizing Guide

5 IVI L

Width/Height 11"x 9" 12"x 10" 13"x 11"

Measurements are sling pouch width x height



# Collar & Cuff Sling

- Features tubular neck pad and flannel lined cuff for optimal patient comfort
- Indicated for injuries to the arm, shoulder, wrist or hand
- Constructed of 1" durable strap with contact closure wrist cuff

Product Number

UNIV

21-2052

Sizing Guide

Fits 6"-10" wrist circumference



# **Strap Arm Sling**

- Features tubular neck pad for added comfort and slide buckle for maximum adjustability
- Indicated for injuries to the arm, shoulder, wrist or hand
- Manufactured from 1½" heavy duty web strap material

Product Number

5	IVI	L	XL
21-5051	21-5052	21-5053	21-5054

Sizing Guide

S	M	L	XL
20" 42"	40" 52"	50" 62"	40" 72

Measurements are minimum / maximum adjustable strap length



# Hanging Cast Sling

- Features envelope style sling with contact closure
- Indicated for control and reduction of edema in the hand and forearm, treatment of acromioclavicular sprains, anterior glenohumeral dislocation and elbow and wrist injuries
- Manufactured from durable muslin material

Product Number

UNIV

54-2170

Sizing Guide

UNIV

Width/Height 6"x 15"

Measurements are sling pouch width x height





# **XR Cubital Tunnel Splint**

- Innovative support that immobilizes and positions the elbow at 122° to help reduce pain and numbness associated with Cubital Tunnel Syndrome
- Indicated for Cubital Tunnel Syndrome, post-surgical application or post-cast removal
- Constructed of rigid A.B.S. material with medial contouring on shell and plush cloth liner (replacement liners also available)

# Product Number

S/M	M/L	L/XL	
37-2191	37-2192	37-2193	

# Sizing Guide

S/M	Fits most females & smaller framed males
M/L	Fits most males & larger framed females
L/XL	Fits large framed males & XL framed females



# **Ventral Cubital Tunnel Splint**

- Unique anterior design immobilizes and positions elbow at 122° while allowing easy access to dressing
- Indicated for Cubital Tunnel Syndrome or post-op application
- Constructed of high density polyethylene shell with plush cloth liner

### Product Number

S/M	M/L	L/XL
37-2291	37-2292	37-2293

### Sizing Guide

M/I F:	es
M/L Fits most males & larger framed female:	S
L/XL Fits large framed males & XL framed fema	ale



# **Ultra Clavicle Strap**

- Features heavy-duty felt center point pad with hook and loop closure for optimal comfort
- Indicated for clavicle fractures, Grade II/III acromiclavicular separations, sternoclavicular dislocations, postural control
- Constructed of breathable, open cell stockinette covered foam with contact closure



### Product Number

2XS	XS	S	М
25-1729	25-1730	25-1731	25-1732
L	XL	2XL	

### Sizing Guide

2XS	XS	S	M
15"- 21"	17"- 24"	25"- 31"	27"- 33"
L	XL	2XL	

Measurement taken from center back, over shoulder under arm and back to center point of back

# Elbow Orthosis w/R.O.M. Hinge

- Anatomically designed humeral and ulnar shells with low-profile R.O.M. hinge provides comfortable, secure and effective range of motion control for the elbow
- Indicated for use following elbow arthroplasty or other post-surgical or post-injury applications where limited range of motion is required
- Constructed of rigid polyethylene laminated shells with closed cell foam lining and easy-to-set R.O.M. hinge; includes stockinette undersleeves

Product Number

	S	M	L
Right	37-2201	37-2202	37-2203
Left	37-2211	37-2212	37-2213

Sizing Guide

HUMERAL	S	M	L
Proximal	9"- 11"	11"- 13"	13"- 15"
ULNAR			
	7" 0"	0" 44"	44" 46"

Measurements taken from widest part of bicep and forearm



# **Universal Elbow Immobilizer**

- Provides comfortable and secure immobilization of the elbow fixed at 90°
- Indicated for trauma stabilization, fractures of the elbow, distal humerus or proximal ulna fractures and chronic olecranon bursitis
- Constructed of foam laminated to pressure sensitive fabric with three rigid stays and bursa pad

Product Number

S/M	L/XL	
43-5071	43-5072	

### Sizing Guide

S/M	Fits most females & smaller framed males
L/XL	Fits most males & larger framed females



# Padded Elbow Wrap

- Features memory foam pad that contours to the elbow to provide optimal protection for the elbow joint; universal wrap design allows for inventory reduction
- Indicated for treatment of elbow injuries, tendonitis, bursitis, arthritis or as athletic protection
- Manufactured from premium grade soft-loop neoprene with memory foam pad and contact closure straps

Product Number

**UNIV** 88-3080

Sizing Guide

# UNIV

Fits 8"-14" elbow circumference





# **Target Elbow Sleeve**

- Provides support, warmth, protection and compression for the elbow
- Indicated for treatment of elbow injuries and athletic protection
- Manufactured from premium grade latex-free neoprene

) <i>"</i>	c	
,	3	

9"	S	M	L	XL
1/8"	88-3053	88-3054	88-3055	88-3056
<sup>3</sup> / <sub>16</sub> "	88-3153	88-3154	88-3155	88-3156
9"	2XL			
<b>9"</b> 1/8"	<b>2XL</b> 88-3057			

### Sizing Guide

S	М	L	XL
8"- 9"	9″- 10″	10"- 11"	11"- 12"
2XL			
12"- 14"			

Measurement taken at widest part of forearm



# Target Elbow Sleeve w/Pad

- Provides support, warmth and compression with additional pad over elbow for added protection
- Indicated for treatment of elbow injuries and athletic protection
- Manufactured from premium grade latex-free neoprene

### Product Number

9"	S	М	L	XL
1/8"	88-3063	88-3064	88-3065	88-3066
3/16"	88-3163	88-3164	88-3165	88-3166
	2XL			
1/8"	<b>2XL</b> 88-3067			

# Sizing Guide

S	M	L	XL
8"- 9"	9"- 10"	10"- 11"	11"- 12"
2XL			
12"- 14"			

Measurement taken at widest part of forearm



# Target Elbow Sleeve w/Strap

- Provides support, warmth and protection for the elbow with strap for added compression over affected tendons
- Indicated for treatment of medial or lateral epicondylitis, elbow injuries and athletic protection
- Manufactured from premium grade latex-free neoprene

### Product Number

9"	S	М	L	XL
1/8"	88-3073	88-3074	88-3075	88-3076
3/ <sub>16</sub> "	88-3173	88-3174	88-3175	88-3176
9"	2XL			
1/8"	88-3077			
78	00-3077			

# Sizing Guide

S	M	L	XL
8"- 9"	9"- 10"	10"- 11"	11"- 12"
2XL			
12"- 14"			

Measurement taken at widest part of forearm

# Target Tennis Elbow Sleeve w/Pad

- Features neoprene pad with circumferential strap for added compression over affected tendons
- Indicated for treatment of medial or lateral epicondylitis
- Manufactured from premium grade latex-free neoprene

Produ	ct Number			
	S	М	L	
1/8"	88-3033	88-3034	88-3035	
<sup>3</sup> / <sub>16</sub> "	88-3043	88-3044	88-3045	
	XL	2XL		
1/8"	88-3036	88-3037		
<sup>3</sup> / <sub>16</sub> "	88-3046	88-3047		
Sizing	Guide			
	S N	1 L	XL	2XL

8	3"- 9"	9"- 10"	10"- 11"	11"- 12"	12"
Massuram	ont take	n at widest	t part of for	roarm	

Measurement taken at widest part of forearn



# **Target Contoured Tennis Elbow Wrap**

- Contoured design places pressure over affected tendons while preserving complete range of motion of the elbow
- Indicated for treatment of medial or lateral epicondylitis
- Manufactured from latex-free soft loop neoprene

Product Number

	UNIV	S	M
1/8"	88-3010	88-3011	88-3012
	L	XL	

Sizing Guide

UNIV	5	IVI	L	XL
8"- 12"	8"- 9"	9"- 10"	10"- 11"	11"- 12"

Measurement taken at widest part of forearm



# Target Universal Tennis Elbow Wrap

- Features easy wrap around design and available with or without pad; universal design allows for inventory reduction
- Indicated for treatment of medial or lateral epicondylitis
- Manufactured from premium grade latex-free neoprene

Product Number

	UNIV	
1/8"	88-3020	
1/8" w/Pad	88-3021	


Sizing Guide

UNIV

8"- 12"

Measurement taken at widest part of forearm





# **Band It Elbow Strap**

- Unique design provides pressure over affected tendons while preserving complete range of motion of the elbow
- Indicated for treatment of medial or lateral epicondylitis
- Constructed of two-piece rigid padded plastic with contact closure strap

Product Number

UNIV

99-5600

Sizing Guide

UNIV

6"-15"

Measurement taken at widest part of forearm



# **Cryo Pneumatic Elbow Wrap**

- Anatomical elbow wrap designed to combine the benefits of pneumatic compression with cold therapy; features pump bulb that detaches once desired compression is achieved
- Indicated for post-op application, as part of a physical therapy or athletic training program; helps diminish swelling / edema, decrease pain and reduce muscle spasms
- Manufactured from plush foam laminate that is Velcro compatible; removable gel pack with tricot lining; non-toxic gel



## Product Number

	0 GELS	1 GEL	2 GELS	
	72-7210	72-7211	72-7212	
Elbow Gel	72-7219			

Sizing Guide

### UNIV

Fits up to 16" elbow circumference



# **Cryotherm Elbow Wrap**

- Economical wrap designed to provide compression and concentrated cold/hot therapy for the elbow
- Indicated for post surgical application, arthritis, strains or sprains, and for integration into rehab program
- Manufactured from soft breathable foam laminate with contact closure and pocket to hold reusable gel packs

Product Number

0 GELS	1 GEL	2 GELS
72-7050	72-7051	72-7052


Sizing Guide

### UNIV

Fits up to 18" bicep circumference

Additional gels sold separately: 72-7003 - Single 3<sup>1</sup>/<sub>2</sub>"x 11" Gel 72-7002 - Case/10 3<sup>1</sup>/<sub>2</sub>"x 11" Gel Measurement taken at widest part of bicep







# **Ultra Fit Wrist Splint**

Our Ultra Fit Wrist Splint provides uniform compression and immobilizes the wrist joint at a neutral angle. Made from a durable lightweight foam laminate the brace breathes well ensuring patient comfort. Preformed malleable aluminum stay helps provide anatomically correct fit. Available in 8" or 10" lengths, sized or universal.

- Features easy one-hand application that provides comfortable support of the wrist while preserving finger mobility; available in two lengths
- Indicated for Carpal Tunnel Syndrome, tendonitis, wrist strains or sprains and post-cast removal
- Manufactured from durable, lightweight material with removable malleable aluminum palmar stay, flexible dorsal stay and contact closure; also available in universal size

Product Number

	UNIV	XS	S	M	L	XL
8" Right	73-1000	73-1002	73-1003	73-1004	73-1005	73-1006
8" Left	73-1001	73-1012	73-1013	73-1014	73-1015	73-1016
	UNIV	XS	S	M	L	XL
10" Right		<b>XS</b> 73-1042	_		_	- · -

# Sizing Guide

UNIV	XS	S	M	L	XL
5"- 9"	4"- 5"	5"- 6"	6"-7"	7"- 8"	8"- 9"

Measurement taken from wrist circumference





For additional sizing, laundering or application instructions visit corflex.com/ultrafitwrist

# Wrist/Hand Orthosis (white)

- Anatomically designed to conform to the natural shape of the wrist providing comfortable and secure immobilization
- Indicated for immobilization of the wrist following post-operative rehabilitation, stabilization of minor forearm and metacarpal fractures, post-cast removal or cast substitute for long-term use
- Constructed of white rigid polyethylene laminated to soft closed cell foam lining; includes two white stockinette undersleeves

Product Number				
	XS	S	М	L
White Right	37-0600	37-0601	37-0602	37-0603
White Left	37-0610	37-0611	37-0612	37-0613
	XL			
White Right	37-0604			
White Left	37-0614			
Sizing Guide				
XS	S	М	L	XL
5" 4"	6"-7"	7"_8"	8"-9"	9"_10"

Measurement taken from wrist circumference



# Wrist/Hand Orthosis (black)

- Anatomically designed to conform to the natural shape of the wrist providing comfortable and secure immobilization
- Indicated for immobilization of the wrist following post-operative rehabilitation, stabilization of minor forearm and metacarpal fractures, post-cast removal or cast substitute for long-term use
- Constructed of black rigid polyethylene laminated to soft closed cell foam lining; includes two black stockinette undersleeves

Product Number				
	XS	S	М	L
Black Right	37-0500	37-0501	37-0502	37-0503
Black Left	37-0510	37-0511	37-0512	37-0513
	XL			
Black Right	37-0504			
Black Left	37-0514			
Sizing Guide				
XS	S	M	L	XL
F" /"	,,, ,,,,	7" 0"	0" 0"	0" 40"

5"-6" 6"-7" 7"-8" 8"-9" 9"Measurement taken from wrist circumference



# **Ultra Fit Cool Wrist Splint**

- Features easy one-hand application that provides cool and comfortable support of the wrist using breathable spacer material
- Indicated for Carpal Tunnel Syndrome, tendonitis, wrist strains or sprains and post-cast removal
- Manufactured from durable, lightweight spacer fabric with removable malleable aluminum palmar stay, flexible dorsal stay and contact closure

### Product Number

	XS	S	М	L	XL
8" Right	73-1022	73-1023	73-1024	73-1025	73-1026
8" Left	73-1032	73-1033	73-1034	73-1035	73-1036

# Sizing Guide

XS	S	M	L	XL
4"- 5"	5″- 6″	6"-7"	7"- 8"	8"-9"

Measurement taken from wrist circumference





Ultra Fit Cool Wrist Splint has the same great features as the Ultra Fit Wrist Splint



# **Signature Vinyl Wrist Splint**

- Features v-style elastic loop locks with ultra thin web strap to provide secure immobilization of wrist while preserving finger mobility
- Indicated for Carpal Tunnel Syndrome, tendonitis, wrist strains or sprains and post-cast removal
- Constructed of premium grade perforated vinyl laminated to soft flannel lining with malleable aluminum palmar and dorsal stays and contact closures



### Product Number

	S	M	L	XL
8" Right	73-5551	73-5552	73-5553	73-5554
8" Left	73-5561	73-5562	73-5563	73-5564
10" Right	73-5571	73-5572	73-5573	73-5574
10" Left	73-5581	73-5582	73-5583	73-5584
	2XS	XS	S	М
5" Right	73-5540	73-5542		
5" Left	73-5541	73-5543		
6" Right			73-5521	73-5522
6" Left			73-5531	73-5532
	L	XL		
6" Right	73-5523	73-5524		
6" Left	73-5533	73-5534		

# Sizing Guide

2XS	XS	S	M	L	XL
3"- 4"	4"- 5"	5"- 6"	6"-7"	7"- 8"	8"- 9"

Measurement taken from wrist circumference



# **Ultra Vinyl Wrist Splint**

- Features v-style loop locks with ultra thin web strap to provide secure immobilization of wrist while preserving finger mobility
- Indicated for Carpal Tunnel Syndrome, tendonitis, wrist strains or sprains and post-cast removal
- Constructed of premium grade perforated vinyl laminated to soft flannel lining with malleable aluminum palmar and dorsal stays and contact closures



### Product Number

	S	M	L	XL
8" Right	41-1521	41-1522	41-1523	41-1524
8" Left	41-1531	41-1532	41-1533	41-1534
10" Right	41-1541	41-1542	41-1543	41-1544
10" Left	41-1551	41-1552	41-1553	41-1554
	2XS	XS		
6" Right	41-1510	41-1513		
6" Left	41-1511	41-1514		

### Sizing Guide

2XS	XS	S	M	L	XL
3"- 4"	4"- 5"	5"- 6"	6"-7"	7"- 8"	8"- 9"

Measurement taken from wrist circumference



# **Suede Wrist Lacer Splint**

- Features single-pull lace closure for easy one-hand application
- Indicated for Carpal Tunnel Syndrome, tendonitis, wrist strains or sprains and post-cast removal
- Manufactured from perforated suede-like material with malleable aluminum palmar and dorsal stay

# Product Number

	XS	S	М	L	XL
8" Right	73-7550	73-7551	73-7552	73-7553	73-7554
8" Left	73-7560	73-7561	73-7562	73-7563	73-7564
	UNIV				
8" Right	<b>UNIV</b> 73-7600				

### Sizing Guide

UNIV	XS	S	М	L	XL
5"-9"	4"- 5"	5″- 6″	6"- 7"	7"- 8"	8"- 9"

Measurement taken from wrist circumference

## **Canvas Wrist Lacer Splint**

- Features cinch laces with contact closure for ease of application
- Indicated for Carpal Tunnel Syndrome, tendonitis and wrist strains or sprains
- Constructed of durable canvas with soft flannel inner lining and malleable palmar stay

#### Product Number

	XS	S	M	L	XL
6" Right	49-5030	49-5031	49-5032	49-5033	49-5034
6" Left	49-5130	49-5131	49-5132	49-5133	49-5134

#### Sizing Guide

XS	S	M	L	XL
4"- 5"	5″- 6″	6"- 7"	7″- 8″	8"- 9"

Measurement taken from wrist circumference



# Signature Foam Flannel Wrist Splint

- Provides support of the wrist while preserving finger mobility
- Indicated for Carpal Tunnel Syndrome, tendonitis, wrist strains or sprains and post-cast removal
- Constructed of cotton with flannel lining; malleable aluminum palmar stay, flexible dorsal stay and contact closure



#### Product Number

	XS	S	M	L	XL
8" Right	73-5920	73-5921	73-5922	73-5923	73-5924
8" Left	73-5930	73-5931	73-5932	73-5933	73-5934
6" Right	73-5720	73-5721	73-5722	73-5723	73-5724
6" Left	73-5730	73-5731	73-5732	73-5733	73-5734
	4XS	3XS			3XS
5" Right	73-5687	73-5688		6" Right	73-5718
5" Left	73-5697	73-5698		6" Left	73-5728
			,		

Measurement taken from wrist circumference



4XS	3XS	XS	S	М	L	XL
21/2"- 31/2"	3"-4"	4"-5"	5"-6"	6"-7"	7"-8"	8"-9"

## **Industrial Wrist Splint**

- Low profile design provides flexible support, making it ideal for the workplace
- Indicated for Carpal Tunnel Syndrome, tendonitis, wrist strains or sprains and repetitive stress injuries
- Manufactured from open-weave material, flexible stays and heavy-duty circumferential elastic wrist strap

#### Product Number

	XS	S	M	L	XL
6" Right	73-5620	73-5621	73-5622	73-5623	73-5624
6" Left	73-5630	73-5631	73-5632	73-5633	73-5634

Sizing Guide

XS	S	M	L	XL
4"- 5"	5"- 6"	6"-7"	7"- 8"	8"- 9"

Measurement taken from wrist circumference



## **Elastic Wrist Splint**

- Contoured design allows for optimal finger mobility while providing moderate wrist support
- Indicated for Carpal Tunnel Syndrome, tendonitis and wrist strains or sprains
- Manufactured from premium grade woven-cotton elastic with palmar stay and contact closure

#### Product Number

	XS	S	M	L	XL
8" Right	73-5770	73-5771	73-5772	73-5773	73-5774
8" Left	73-5780	73-5781	73-5782	73-5783	73-5784
	XS	S	M	L	XL
6" Right			<b>M</b> 73-5752	_	, ·-

## Sizing Guide

XS	S	M	L	XL
4"- 5"	5″- 6″	6"-7"	7"- 8"	8"- 9"

Measurement taken from wrist circumference







- Anatomical wrist wrap designed to combine the benefits of pneumatic compression with cold therapy; features pump bulb that detaches once desired compression is achieved
- Indicated for post-op application, severe sprains or strains, as part of a physical therapy or occupational therapy program; helps diminish swelling/edema
- Manufactured from plush foam laminate that is Velcro compatible; removable gel pack with tricot lining; non-toxic gel; movable palmar stay

Product Number

	0 GELS	1 GEL	2 GELS	
Right	72-7240	72-7241	72-7242	
Left	72-7250	72-7251	72-7252	
Wrist Gel	72-7259			

Sizing Guide

Fits up to 9" wrist circumference





## **Cryotherm Wrist Wrap**

- Economical wrap designed to provide compression and concentrated cold/hot therapy for the wrist
- Indicated for post-op use, arthritis, wrist strains or sprains and for integration into rehab program
- Manufactured from a soft breathable foam laminate with contact closure and pockets to hold reusable gel packs

Product Number

0 GELS	1 GEL	2 GELS
72-7150	72-7151	72-7152

Sizing Guide

#### UNIV

Fits up to 9" wrist circumference

Additional gels sold separately: 72-7003 - Single 3 1/2"x 11" Gel 72-7002 - Case/10 3 1/2"x 11" Gel





# **Ambo Elastic Wrist Splint**

- > Features reversible palmar stay that allows it to fit either left or right hand
- Indicated for Carpal Tunnel Syndrome, tendonitis and wrist strains or sprains
- Manufactured from premium grade elastic with contact closure

Product Number

	S	M	L	XL
6"	41-1751	41-1752	41-1753	41-1754
8"	41-1771	41-1772	41-1773	41-1774
10"	41-1791	41-1792	41-1793	41-1794

Sizing Guide

S	M	L	XL
5"- 6"	6"-7"	7″- 8″	8″- 9

Measurement taken from wrist circumference



## Wrist & Forearm Immobilizer

- Features palmar stay with movable dorsal stay/strap assembly to provide wrist immobilization; universal design allows for inventory reduction
- Indicated for Carpal Tunnel Syndrome, tendonitis and wrist strains or sprains
- Constructed of soft breathable foam laminated with contact closure



#### Product Number

	UNIV
Right	41-1796
Left	41-1797
X-Wide Right	41-1798
X-Wide Left	41-1799

Sizing Guide

UNIV	UNIV X-WIDE
5″- 10″	7″- 13″

Measurement taken from wrist circumference

## Target Wrist Stabilizer w/Stay

- Easy to apply wrap-around design provides support, warmth, compression and protection for the wrist
- Indicated for Carpal Tunnel Syndrome, tendonitis and wrist strains or sprains
- Manufactured from latex-free soft loop neoprene with removable malleable palmar stay



	UNIV R	UNIV L
7" Pedi	88-2056	88-2057
8" Standard	88-2060	88-2070
9" Extended	88-2061	88-2071

Sizing Guide

	UNIV
Standard/Extended	Fits up to 9" wrist circumference
Pedi	Fits up to 7" wrist circumference



## **Target Wrist Splint**

- Designed to provide support, warmth, compression and protection for the wrist
- Indicated for Carpal Tunnel Syndrome, tendonitis and wrist strains or sprains
- Manufactured from latex-free soft loop neoprene with removable malleable palmar stay and crossdirectional closure straps

Product Number

	UNIV
Right	88-2020
Left	88-2030

## Sizing Guide

## UNIV

Fits up to 9" wrist circumference



## **Target Wrist Stabilizer**

- Easy to apply wrap-around design provides warmth, compression and protection for the wrist
- Indicated for Carpal Tunnel Syndrome, tendonitis and wrist strains or sprains
- Manufactured from latex-free soft loop neoprene

Product Number

**UNIV** 88-2050

Sizing Guide

#### UNIV

Fits up to 9" wrist circumference



## **Target Wrist Wrap**

- Easy to apply design provides warmth and protection for the wrist
- Indicated for wrist sprains or strains and athletic protection
- Manufactured from premium grade latex-free neoprene

Product Number

UNIV

88-2010

Sizing Guide

#### UNIV

Fits up to 9" wrist circumference





## Ultra Fit Wrist Splint w/Abducted Thumb

Our Ultra Fit Wrist Splint with Abducted Thumb provides consistant clinical results. Designed with dual adjustable aluminum stays that support the palmar surface of the wrist and extensor surface of the thumb. Available in two materials - durable lightweight tricot or breathable cool spacer fabric (see p. 39). May be used for tendonitis, sprains and strains, and post-cast removal.

- Features easy one-hand application and provides comfortable thumb and wrist immobilization while preserving mobility of the other fingers
- Indicated for deQuervain's Syndrome, tendonitis, thumb strains or sprains and post-cast removal
- Manufactured from durable, lightweight material with removable malleable aluminum palmar and thumb stays, flexible dorsal stay and contact closure

Product Number

	XS	S	M	L	XL
8" Right	73-5302	73-5303	73-5304	73-5305	73-5306
8" Left	73-5312	73-5313	73-5314	73-5315	73-5316

Sizing Guide

XS	S	M	L	XL
4"- 5"	5″- 6″	6"-7"	7"- 8"	8"- 9"

Measurement taken from wrist circumference





For additional sizing, laundering or application instructions visit corflex.com/ultrafitwristabthumb

#### **Thumbster Soft**

- Anatomically designed to provide comfortable and secure immobilization of the thumb and wrist while preserving mobility of the other fingers
- Indicated for deQuervain's Syndrome, gamekeeper's thumb, tendonitis, thumb strains or sprains and post-cast removal
- Constructed of vacuum-formed cross linked polyester with stretchable contact closure



Proc	luct	Num	ber
------	------	-----	-----

	UNIV	PEDI	S/M	M/L
Univ	73-5000	73-5011		
Right			73-5001	73-5002
Left			73-5003	73-5004

#### Sizing Guide

PEDI Fits most children & young adolescents S/M Fits most females & smaller framed males
S/M Fits most females & smaller framed males
M/L Fits most males & larger framed females



# Suede Wrist Lacer Splint w/Abducted Thumb

- Features single-pull lace closure for ease of application and provides comfortable thumb and wrist immobilization while preserving mobility of the other fingers
- Indicated for deQuervain's Syndrome, tendonitis, thumb strains or sprains and post-cast removal
- Manufactured from perforated suede-like material with malleable aluminum palmar and dorsal stay

#### Product Number

	XS	S	M	L
8" Right	73-8550	73-8551	73-8552	73-8553
8" Left	73-8560	73-8561	73-8562	73-8563
	XL			
8" Right				

## Sizing Guide

Sizing G	XS	S	М	L	
	4"- 5"	5"- 6"	6"-7"	7"-8"	
	XL				
	8"- 9"				

Measurement taken from wrist circumference



## Ultra Fit Cool Wrist Splint w/Abducted Thumb

- Features easy one-hand application and provides cool and comfortable thumb and wrist immobilization using breathable spacer fabric
- Indicated for deQuervain's Syndrome, tendonitis, thumb strains or sprains and post-cast removal
- Manufactured from durable, lightweight spacer fabric with removable malleable aluminum palmar and thumb stays, flexible dorsal stay and contact closure

#### Product Number

		XS	5	IVI	L	
8"	Right	73-5322	73-5323	73-5324	73-5325	
8"	Left	73-5332	73-5333	73-5334	73-5335	
		XL				
8"	Right	<b>XL</b> 73-5326				

#### Sizing Guide

L
7" 7"-8"

Measurement taken from wrist circumference





Ultra Fit Cool Wrist Splint w/Abducted Thumb has the same great features as the Ultra Fit Wrist Splint w/Abducted Thumb Splint



40



#### Wrist/Hand/Thumb Orthosis (white)

- Anatomically designed to conform to the natural shape of the thumb and wrist providing comfortable and secure immobilization
- Indicated for immobilization of the thumb for post-operative rehabilitation, stabilization of thumb fractures, postcast removal or cast substitute for long-term use
- Constructed of white rigid polyethylene laminated to soft closed cell foam lining; includes two white stockinette undersleeves

Product Number

	XS	S	M	L
White Right	37-0620	37-0621	37-0622	37-0623
White Left	37-0630	37-0631	37-0632	37-0633

XL

White Right 37-0624 White Left 37-0634

Sizing Guide

XS	S	M	L
5"- 6"	6"-7"	7"- 8"	8"- 9"
XL			
9"- 10"			

Measurement taken from wrist circumference



#### Wrist/Hand/Thumb Orthosis (black)

- Anatomically designed to conform to the natural shape of the thumb and wrist providing comfortable and secure immobilization
- Indicated for immobilization of the thumb for post-operative rehabilitation, stabilization of thumb fractures, postcast removal or cast substitute for long-term use
- Constructed of black rigid polyethylene laminated to soft closed cell foam lining; includes two black stockinette undersleeves

#### Product Number

	XS	S	M	L
Black Right	37-0520	37-0521	37-0522	37-0523
Black Left	37-0530	37-0531	37-0532	37-0533

XL

Black Right 37-0524 Black Left 37-0534

Sizing Guide

XS	S	M	L
5"- 6"	6"-7"	7"-8"	8"-9"
XL			
9"-10"			

Measurement taken from wrist circumference



This symbol indicates that pediatric sizes are available for this product.



This symbol indicates that bariatric sizes are available for this product.



This symbol certifies that all products manufactured by Corflex are latex free.



This symbol indicates that certain precautions are needed when using this product.



This symbol identifies CoolTex products that are manufactured from a neoprene-alternative material.

## Low Profile Thumbster

- Anatomically designed to provide comfortable and secure immobilization of the thumb while preserving mobility of the other fingers and wrist
- Indicated for deQuervain's Syndrome, gamekeeper's thumb, tendonitis, thumb strains or sprains and post-cast removal
- Manufactured from ultra-thin, rigid polyethylene with contact closure; includes two stockinette undersleeves

Product Number

	S	M	L	XL
Right	37-5001	37-5002	37-5003	37-5004
Left	37-5011	37-5012	37-5013	37-5014

Sizing Guide

S	M	L	XL	
6"-7"	7"- 8"	8"- 9"	9"- 10"	

Measurement taken from wrist circumference



## **Target Hitchhiker Thumb**

- Designed to provide comfortable immobilization of the thumb and added support for the wrist while preserving mobility of the other fingers
- Indicated for deQuervain's Syndrome, gamekeeper's thumb, tendonitis, and thumb strains or sprains
- Manufactured from premium grade latex-free neoprene with malleable thumb stay and elastic straps

Product Number

# UNIV

88-2080

#### Sizing Guide

#### UNIV

Fits up to 9" wrist circumference



## **Target Sewn Thumb**

- Designed to provide moderate support, warmth, compression and protection for the thumb; available with stays for added support and immobilization
- Indicated for tendonitis and thumb strains or sprains
- Manufactured from premium grade latex-free neoprene; available with optional plastic stays

Product Number

**UNIV** 88-2040

w/Stays 88-2042

Sizing Guide

UNIV

Fits up to 9" wrist circumference



Size

# Upper

## **Digital Splint**

- Unique design allows for use on either palmar or dorsal side of hand to help immobilize any finger, including the thumb (to be used with padded aluminum finger strips; sold separately)
- Indicated for treatment of phalangeal fractures
- Manufactured of durable perforated vinyl with contact closure



# Finger Splints

- Provides comfortable immobilization and protection for the finger; available in packages of 12 or as part of an assortment kit
- Indicated for treatment of finger fractures, dislocations or as protection following injury
- Constructed of malleable aluminum with foam padding

## **Stax Splints**

- Designed to provide immobilization of the distal finger or thumb, while allowing for examination of the finger tip and nailbed; available in packages of 12 or as part of an assortment kit
- Indicated for treatment of finger tip injuries, mallet finger or as protection following injury
- Manufactured from durable plastic

Product Number and Sizing Guide

	UNIV			
	41-3050			
its up to 9"	9" wrist circumference			
<sup>1</sup> / <sub>2</sub> "	3"	41/2"	6"	9"
9-3201	49-3202	49-3205	49-3203	49-3204
1/2"	3"	41/2"	6"	9"
	S	М	L	
	49-3701	49-3702	49-3703	
	2"	21/2"	3"	
	S	М	L	
	49-3501	49-3502	49-3503	
	3"	31/2"	4"	
	S	М	L	
	49-3401	49-3402	49-3403	
	<sup>1</sup> / <sub>2</sub> "x 18"	<sup>3</sup> / <sub>4</sub> "x 18"	1"x 18"	
S	S	М	L	
9-3600	49-3601	49-3602	49-3603	
"	2"	3"	4"	
S	S	М	L	
9-3300	49-3301	49-3302	49-3303	
II .	2"	3"	4"	
1	#2	#3	#4	#5
9-3801	49-3802	49-3803	49-3804	49-3805
"	21/8"	21/4"	21/2"	25/8"
5.5	#6	#7	KIT	
9-3806	49-3807	49-3808	49-3800	
3/4"	27/8"	31/4"	35 pc Assor	tment
it	KIT			
	S 9-3600 "S 9-3800 "1 19-3801 "7-3801 "7-3806	sits up to 9" wrist circumf 1/2" 3" 9-3201 49-3202 1/2" 3"  \$ 49-3701 2" \$ 49-3501 3" \$ 49-3501 3" \$ \$ 49-3401 1/2"x 18" \$ \$ 9-3600 49-3601 2" \$ \$ 9-3800 49-3301 2" \$ \$ 9-3801 49-3802 2" \$ 1 #2 9-3801 49-3802 2" 6 9-3806 49-3807 27/8"	## are selected as a selected	## ## ## ## ## ## ## ## ## ## ## ## ##



45 pc Assortment

## **Forearm Splints**

- Pre-molded design allows for optimal anatomic fit; available in various styles and sizes
- Indicated for treatment of fractures, severe sprains and strains of the wrist and post-operative rehabilitation
- Constructed of malleable aluminum with open cell foam padding

Product Number and Sizing Guide

Forearm Splint XS

Forearm Splint	XS	S	М	L
Right	49-2470	49-2471	49-2472	49-2473
Left	49-2570	49-2571	49-2572	49-2573
Length	6"	8"	10"	12"
Metacarpal Hand Sp	lint	S	M	L
Right		49-2891	49-2892	49-2893
Left		49-2991	49-2992	49-2993
Length		4"	5"	6"
Metacarpal Splint		S	M	L
Right		49-2681	49-2682	49-2683
Left		49-2781	49-2782	49-2783
Length		8"	10"	12"
Cock-up Splint		S	M	L
		49-2161	49-2162	49-2163
Length		8"	10"	12"
Colles Splint		S	М	L
Right		49-2261	49-2262	49-2263
Left		49-2361	49-2362	49-2363
Length		8"	10"	12"





## **Ulnar Splint**

- Design provides for secure fit while allowing full range of wrist and elbow motion
- Indicated for mid shaft and distal ulnar fractures and post-cast removal
- Constructed of rigid polyethylene laminated to soft closed cell foam lining with contact closure

#### Product Number

S	M	L	XL
37-2071	37-2072	37-2073	37-2074

#### Sizing Guide

	S	М	L	XL
Proximal	7"- 9"	9"- 11"	11"- 13"	13"- 15"
Distal	6"-7"	7"- 8"	8"- 9"	9"- 10"

Measurement taken from widest part of forearm and above wrist



## **Target MP Ulnar Drift Support**

- Provides support and maintains neutral alignment of the MP joints while preserving mobility of the fingers
- Indicated for treatment of arthritis and MP joint deformity known as ulnar drift
- Manufactured from latex-free soft loop neoprene with contact closure

#### Product Number

	S/M	M/L
Right	88-8500	88-8501
Left	88-8510	88-8511

#### Sizing Guide

S/M	M/L
21/2"- 31/2"	31/2"- 41/2"

Measurement taken is width at MP joints

## **Extended Length Boxer Splint**

- Anatomically designed to provide comfortable support, protection and immobilization of the 4th and 5th metacarpals; extended length provides added wrist support; average length 95/8"
- Indicated for treatment of acute or sub-acute fractures of the 4th and 5th metacarpals commonly referred to as "boxer fractures"
- Constructed of lightweight lowdensity polyethylene shell with soft foam liner and palmar pad and closure straps

Product Num	iper			
	S	M	L	XL
Right	37-1221	37-1222	37-1223	37-1224
Left	37-1231	37-1232	37-1233	37-1234

9	Sizing	Guide	
			9

	S	M	L	XL
Wrist	4"- 5"	5"- 6"	6"-7"	7"- 8"
4th/5th MCP	23/4"- 31/4"	31/4"- 33/4"	33/4"- 41/4"	41/4"- 43/4"
			1.4.1	0 5.1

Measurement taken from wrist circumference and 4th & 5th



## **Boxer Splint**

- Anatomically designed to provide comfortable support, protection and immobilization of the 4th and 5th metacarpals; average length 65/8"
- Indicated for treatment of acute or sub-acute fractures of the 4th and 5th metacarpals commonly referred to as "boxer fractures"
- Constructed of lightweight lowdensity polyethylene shell with soft foam liner and palmar pad and closure straps

Product Nu	ımber			
	S	M	L	XL
Right	37-1201	37-1202	37-1203	37-1204
Left	37-1211	37-1212	37-1213	37-1214

Sizing Guide

	S	M	L	XL
Wrist	4"- 5"	5"- 6"	6"-7"	7"- 8"
4th/5th MCP	23/4"- 31/4"	31/4"- 33/4"	33/4"- 41/4"	41/4"- 43/4"

Measurement taken from wrist circumference and 4th & 5th



"I designed the Boxer Splint to treat acute fractures of the 4th and 5th metacarpals. It provides comfortable support, protection and immobilization."

Douglas Goumas, MD, Orthopedic Surgeon





## Compression Knee Immobilizer

All styles of our Compression Knee Immobilizer offer optimal immobilization of the knee and maximum patient comfort. Ideal for doctor's offices, emergency rooms and surgery centers because of it's universal sizing. Available with optional pockets to accommodate reusable gel packs, making it popular for post-op use.

- Features sewn-in contoured posterior stays, movable medial/lateral stays and elastic compression straps for optimal immobilization; available with pockets to hold gel packs for cold/hot therapy
- Indicated for post-op immobilization and mild to severe knee injuries where full extension of the leg is desired
- Constructed of breathable foam laminate, aluminum stays, premium quality elastic straps and contact closure

Product N	Number				Sizing Guide
LENGTH	NO	POCKETS	POCKETS	POCKETS	FITS UP TO
ı	POCKETS	0 GELS	2 GELS	4 GELS	CIRC.
6"	75-0651				11"
9"	75-0951				14"
101/2"	75-1051				14"
12"	75-1251				24"
14"	75-1451	75-1450	75-1452	75-1454	28"
16"	75-1651	75-1650	75-1652	75-1654	28"
18"	75-1851	75-1850	75-1852	75-1854	28"
20"	75-2051	75-2050	75-2052	75-2054	29"
22"	75-2251	75-2250	75-2252	75-2254	29"
24"	75-2451	75-2450	75-2452	75-2454	31"



Measurement taken from leg circumference 4" above mid-patella Additional gels sold separately:

72-7001 - Single 5 <sup>1</sup>/<sub>2</sub>"x 11" Gel 72-7000 - Case/10 5 1/2"x 11" Gel





For additional sizing, laundering or application instructions visit corflex.com/compressionkneeimmob



#### Universal Knee Immobilizer

- → Universal design allows for inventory reduction; features contoured posterior stays and movable medial/lateral stay assemblies for optimal positioning and immobilization with 4" elastic patella strap for added compression; available in x-wide width
- Indicated for post-op knee immobilization and mild to severe knee injuries where full extension of the leg is desired
- Constructed of breathable foam laminate to soft loop pile with hook and loop closure straps, premium quality elastic straps and aluminum stays





#### Product Number

LENGTH	6"	9"	101/2"	12"	14"
Univ	51-1506	51-1509	51-5510	51-1512	51-1514
X-Wide				92-5112	92-5114
LENGTH	16"	18"	20"	22"	24"
<b>LENGTH</b> Univ		. •	<b>20"</b> 51-1520		

#### Sizing Guide

	6"	9"	101/2"	12"	14"
Univ	11"	14"	151/2"	20"	20"
X-Wide				25"	25"
	16"	18"	20"	22"	24"
Univ	<b>16"</b> 24"	<b>18"</b> 24"	<b>20"</b> 25"	<b>22"</b> 25"	<b>24"</b> 27"

Measurement taken from leg circumference 4" above mid-patella



## Tri Panel Knee Immobilizer

- Features movable medial/lateral stay panels and sewn-in contoured posterior stays to provide secure immobilization; universal design allows for inventory reduction
- Indicated for post-op knee immobilization and mild to severe knee injuries where full extension of the leg is desired
- Constructed of soft foam material with hook and loop closure straps and aluminum stays

#### Product Number

LENGTH	13"	17"	19"	23"	26"
	51-7113	51-7117	51-7119	51-7123	51-7126

Sizing Guide

#### Fits up to 26" circumference

Measurement taken from leg circumference 4" above mid-patella



## **Ultra Tricot Knee Immobilizer**

- Features sewn-in contoured posterior stays and medial/lateral stays with patella opening for optimal immobilization
- Indicated for post-op knee immobilization and mild to severe knee injuries where full extension of the leg is desired
- Constructed of breathable foam with canvas outer, hook and loop closure straps and aluminum stays



# Product Number

LEINGIL	2/2	VO	3	IVI
13"	51-1529	51-1530	51-1531	51-1532
17"		51-1540	51-1541	51-1542
19"			51-1551	51-1552
23"			51-1561	51-1562
	L	XL		
13"	<b>L</b> 51-1533			
13" 17"	51-1533			
	51-1533	51-1534		
17"	51-1533 51-1543	51-1534 51-1544 51-1554		

#### Sizing Guide

3121119 00	iiac			
	2XS	XS	S	M
	7"- 10"	10"- 13"	13"- 16"	16"- 18"
	L	XL		
	18"- 20"	20"- 22"		

Measurement taken from leg circumference 4" above mid-patella

## **CDS Knee Orthosis**

- Dynamic orthosis features easy-to-set tension spring hinge with multiple extension or flexion adjustments that provides low load prolonged stretch to the connective tissue, stretching patient into knee extension
- Indicated for contractures, R.O.M. limitations, patello-femoral syndrome, post-cast removal and post-surgical applications
- Constructed of lightweight aluminum uprights with double-faced loop adjustable straps with quik-clip closure

Product Number

	RIGHT	LEFT
Extension	89-1001	89-1011
Flexion	89-1002	89-1012

Sizing Guide

Fits up to 29" thigh circumference

Measurement taken 6" above mid-patella



## **Cryotherm Knee Wrap**

- Economical wrap designed to provide compression and concentrated cold/hot therapy for the knee; available with 2 pockets or 3 pockets for circumferential coverage
- Indicated for post-op use, arthritis, strains or sprains and for integration into rehab program
- Manufactured from soft breathable foam laminate with contact closure and pockets to hold reusable gel packs



#### Product Number

	0 GELS	2 GELS	3 GELS	
2 Pockets	72-7020	72-7022		
3 Pockets	72-7080		72-7083	
	4 GELS	6 GELS		
2 Pockets	72-7024			
3 Pockets		72-7086		

Sizing Guide

UNIV

Fits up to 24" circumference

Measurement taken from leg circumference 4" above mid-patella Additional gels sold separately: 72-7001 - Single 5 <sup>1</sup>/<sub>2</sub>"x 11" Gel

72-7001 - Single 5 7/2 x 11" Gel 72-7000 - Case/10 5 1/2"x 11" Gel



# Cryotherm Vari Wrap

- Economical wrap designed to accommodate various parts of the body and provide compression and concentrated cold/hot therapy
- Indicated for post-op use, arthritis, strains or sprains and for integration into rehab program
- Manufactured from soft breathable foam laminate with contact closure and pocket to hold reusable gel pack

Product Number

**0 GELS 1 GEL 2 GELS** 72-7070 72-7071 72-7072

Sizing Guide

#### UNIV

Fits up to 40" circumference

Additional gels sold separately: 72-7001 - Single 5  $^{1}/_{2}$ "x 11" Gel 72-7000 - Case/10 5  $^{1}/_{2}$ "x 11" Gel







## **Cryo Pneumatic Knee Splint**

Uniquely designed to combine the benefits of pneumatic compression and cold therapy. This anatomical wrap is ideal for post-surgical application or as part of a physical therapy or athletic training program to help reduce edema and decrease pain. It features an easy to use detachable pump and contoured reusable gel pack. Available in a universal size that fits either left or right.

- Anatomical knee wrap designed to combine the benefits of pneumatic compression with cold therapy; features pump bulb that detaches once desired compression is achieved
- Indicated for post-op application and as part of a physical therapy or athletic training program; helps diminish swelling/edema, decrease pain and reduce muscle spasms
- Manufactured from plush foam laminate with contact closure straps, removable gel pack with tricot lining and non-toxic gel; additional universal extension straps available

#### Product Number

	0 GELS	1 GEL	2 GELS	
	72-7170	72-7171	72-7172	
			UNIV	
Knee Gel			72-7179	
Extension	n Straps (P	kg/10)	72-7190	

#### Sizing Guide



Fits up to 24" circumference





For additional sizing, laundering or application instructions visit corflex.com/cryopneumaticknee

## Cryo Pneumatic Knee Orthosis w/Hinge

- Anatomical knee wrap that features bi-lateral hinges providing medial/ lateral support; hinges may be contoured and feature a built in stop at full extension; includes two extension straps to facilitate hinge placement
- Indicated for post knee arthroscopy, meniscal tears, osteoarthritis pain and swelling and patellofemoral pain
- Manufactured from plush foam laminate with hook and loop closures; reusable gel packs feature soft fabric lining and non-toxic gel; includes two aluminum hinges and extension straps

Product Number					
	0 GELS	1 GEL	2 GELS		
	72-7270	72-7271	72-7272		
			UNIV		
Knee G	el		72-7179		

72-7275

Sizing Guide

#### UNIV

Hinge Set w/Extension Straps

Fits up to 24" circumference





## Cryo Pneumatic Knee Orthosis w/R.O.M. Hinge

- Features ROM Hinges with malleable aluminum uprights offering full flexion and extension control with settings at 0°,15°, 30°, 60°, and 90°
- Indicated following ACL repair, post total joint surgery, meniscal tears, patella realignment or during rehabilitation
- Manufactured from plush foam laminate with hook and loop closures; reusable gel packs feature soft fabric lining and non-toxic gel; includes two aluminum ROM hinges with plastic covers and extension straps



Product Number

	0 GELS	1 GEL	2 GELS	
	72-7370	72-7371	72-7372	
			UNIV	
Knee Gel	72-7179			
R.O.M. H	72-7375			

Sizing Guide

#### UNIV

Fits up to 24" circumference





#### **Contender Knee Brace**

One of our most original designs, the Contender offers the convenience of a wrap above the knee and a sleeve below the knee. It is recommended for noncontact sports and activities of daily living. The polycentric hinge with contoured condyle shell has aluminum uprights and offers full flexion and extension control with settings at 0°,15°,30°,60° and 90°. Also available in CoolTex, a neoprene alternative material.

- Features R.O.M. hinges with contoured condyle shells that allow for full flexion and extension control along with convenience of a sleeve below the knee and benefits of an adjustable wrap above the knee; available in two lengths
- Indicated for medial/lateral instability, muscle/ligament weakness or injury, ACL/PCL moderate instability, strains or sprains and degenerative joint disease where extension/ flexion control is required
- Manufactured from either premium grade latex-free neoprene with suede hinge cover or CoolTex, a neoprene alternative

#### Product Number


<sup>3</sup> / <sub>16</sub> "	S	M	L	XL	2XL
13"	88-8403	88-8404	88-8405	88-8406	88-8407
16"	88-8453	88-8454	88-8455	88-8456	88-8457
16" CoolTex	86-8453	86-8454	86-8455	86-8456	86-8457
<sup>3</sup> / <sub>16</sub> "	3XL	4XL			
13"	88-8408	88-8409			
16"	88-8458	88-8459			

#### Sizina Guide

5					
	6" ABOVE MID-PATELLA	MID-PATELLA	6" BELOW MID-PATELLA		
Small	16"- 18"	12"- 14"	11"- 13"		
Medium	18"- 19"	14"- 16"	13"- 15"		
Large	19"- 22"	16"- 18"	15"- 17"		
X-Large	22"- 24"	18"- 20"	17"- 19"		
2X-Large	24"- 26"	20"- 22"	19"- 21"		
3X-Large	26"- 28"	22"- 24"	21"- 23"		
4X-Large	28"- 30"	24"- 26"	23"- 24"		









For additional sizing, laundering or application instructions visit corflex.com/contenderknee

## Posterior Adjustable Knee Sleeve w/R.O.M. Hinge

- Features R.O.M. hinges that allow for full flexion and extension control and unique design that offers the convenience of a sleeve below the knee and benefits of an adjustable wrap above the knee; available in two lengths
- Indicated for medial/lateral instability, muscle/ligament weakness or injury, ACL/PCL moderate instability, strains or sprains and degenerative joint disease where extension/flexion control is required
- Manufactured from premium grade latex-free neoprene, or CoolTex , a neoprene alternative

Product Number					
<sup>3</sup> / <sub>16</sub> "	S	M	L	XL	
13"	88-8133	88-8134	88-8135	88-8136	
16"	88-8213	88-8214	88-8215	88-8216	
13" CoolTex	86-8133	86-8134	86-8135	86-8136	
16" CoolTex	86-8213	86-8214	86-8215	86-8216	
<sup>3</sup> / <sub>16</sub> "	2XL	3XL	4XL		
13"	88-8137	88-8138	88-8139		
16"	88-8217	88-8218	88-8219		
13" CoolTex	86-8137	86-8138	86-8139		
16" CoolTex	86-8217	86-8218	86-8219		

Sizing Guide

SAME AS CONTENDER KNEE BRACE SIZING GUIDE PAGE 52





## Posterior Adjustable Knee Sleeve w/Hinge

- Unique design offers the convenience of a sleeve below the knee and benefits of an adjustable wrap above the knee with heavy-duty hinges for optimal medial/lateral support; available in two lengths
- Indicated for medial/lateral instability, muscle/ligament weakness or injury, ACL/PCL moderate instability, strains or sprains and degenerative joint disease
- Manufactured from premium grade latex-free neoprene



Hoducti	TOGUCE INGITIDE					
<sup>3</sup> / <sub>16</sub> "	S	M	L	XL		
13"	88-8083	88-8084	88-8085	88-8086		
16"	88-8203	88-8204	88-8205	88-8206		
<sup>3</sup> / <sub>16</sub> "	2XL	3XL	4XL			
13"	88-8087	88-8088	88-8089			
16"	88-8207	88-8208	88-8209			


Sizing Guide

SAME AS CONTENDER KNEE BRACE SIZING GUIDE PAGE 52



# Posterior Adjustable Knee Sleeve w/Cor-Trak Buttress & Hinge

- Unique design offers the convenience of a sleeve below the knee and benefits of an adjustable wrap above the knee with moveable U-shaped Cor-Trak buttress for patella stabilization and heavy-duty hinges for optimal medial/lateral support
- Indicated for medial/lateral instability, muscle/ligament weakness or injury, ACL/PCL moderate instability, strains or sprains and degenerative joint disease where patella stabilization is required
- Manufactured from premium grade latex-free neoprene and compression molded textured neoprene buttress



Sizing Guide

SAME AS CONTENDER KNEE BRACE SIZING GUIDE PAGE 52







## Anterior Closure Knee Wrap w/R.O.M. Hinge

- Features R.O.M. hinges that allow for full flexion and extension control and wrap-around design for ease of application; available in two lengths
- Indicated for medial/lateral instability, muscle/ligament weakness or injury, and degenerative joint disease where extension/flexion control is required
- Manufactured from premium grade latex-free neoprene



Product Number						
	<sup>3</sup> / <sub>16</sub> "	S	M	L	XL	2XL
	13"	88-1353	88-1354	88-1355	88-1356	88-1357
	13" Op Pop	88-1343	88-1344	88-1345	88-1346	88-1347
	16"	88-1693	88-1694	88-1695	88-1696	88-1697
	16" Op Pop	88-1683	88-1684	88-1685	88-1686	88-1687
	<sup>3</sup> / <sub>16</sub> "	3XL	4XL			
	13"	88-1358	88-1359			
	13" Op Pop	88-1348	88-1349			
	16"	88-1698	88-1699			
	16" Op Pop	88-1688	88-1689			

#### Sizing Guide

Product Number

	6" above Mid-patella	MID-PATELLA	6" BELOW MID-PATELLA
Small	16"- 18"	12"- 14"	11"- 13"
Medium	18"- 19"	14"- 16"	13"- 15"
Large	19"- 21"	16"- 18"	15"- 17"
X-Large	21"- 23"	18"- 20"	17"- 19"
2X-Large	23"- 25"	20"- 22"	19"- 21"
3X-Large	25"- 27"	22"- 24"	21"- 22"
4X-Large	27"- 28"	24"- 26"	22"- 23"



## Anterior Closure Knee Wrap w/Hinge

- Features wrap-around design for ease of application with heavy-duty hinges for optimal medial/lateral support; available in two lengths
- Indicated for medial/lateral instability, muscle/ligament weakness or injury, and degenerative joint disease
- Manufactured from premium grade latex-free neoprene



1/8"	S	M	L	XL	2XL
13"	88-1213	88-1214	88-1215	88-1216	88-1217
13" Op Pop	88-1203	88-1204	88-1205	88-1206	88-1207
1/8"	3XL	4XL			
13"	88-1218	88-1219			

13" Op Pop 88-1208 88-1209 3/<sub>16</sub>" 2XL 13" 88-1313 88-1314 88-1315 88-1316 88-1317 13" Op Pop 88-1303 88-1304 88-1305 88-1306 88-1307 88-1673 88-1674 88-1675 88-1676 88-1677 16" Op Pop 88-1663 88-1664 88-1665 88-1666 88-1667 3/16" 3XL 4XL

13" 88-1318 88-1319 13" Op Pop 88-1308 88-1309 16" 88-1678 88-1679 16" Op Pop 88-1668 88-1669

Sizing Guide

SAME AS ANTERIOR CLOSURE KNEE WRAP W/R.O.M. SIZING GUIDE ON THIS PAGE

## CoolTex Anterior Closure Knee Wrap w/R.O.M. Hinge

- → Features R.O.M. hinges that allow for full flexion and extension control and wrap-around design for ease of application; in CoolTex, a neoprene alternative material
- Indicated for medial/lateral instability, muscle/ligament weakness or injury, and degenerative joint disease where extension/flexion control is required; ideal for patients with sensitivity to neoprene
- Manufactured from lightweight, breathable CoolTex material







CoolTex	S	М	L	XL	2XL
16"	86-1693	86-1694	86-1695	86-1696	86-1697
16" Op Pop	86-1683	86-1684	86-1685	86-1686	86-1687
CoolTex	3XL	4XL			
16"	86-1698	86-1699			

Sizing Guide

SAME AS ANTERIOR CLOSURE KNEE WRAP W/R.O.M. SIZING GUIDE PAGE 54

16" Op Pop 86-1688 86-1689



## CoolTex Anterior Closure Knee Wrap w/Hinge

- Features wrap-around design for ease of application with heavy-duty hinges for optimal medial/lateral support; in CoolTex, a neoprene alternative material
- Indicated for medial/lateral instability, muscle/ligament weakness or injury, and degenerative joint disease; ideal for patients with sensitivity to neoprene
- Manufactured from lightweight, breathable CoolTex material





Product Number

CoolTex	S	M	L	XL	2XL
13"	86-1313	86-1314	86-1315	86-1316	86-1317
13" Op Pop	86-1303	86-1304	86-1305	86-1306	86-1307
CoolTex	3XL	4XL			
CoolTex 13"	<b>3XL</b> 86-1318	.,			

Sizing Guide

SAME AS ANTERIOR CLOSURE KNEE WRAP W/R.O.M. SIZING GUIDE PAGE 54





## Knee Sleeve w/R.O.M. Hinge

- → Features R.O.M. hinges that allow for full flexion and extension control and sleeve design
- Indicated for medial/lateral instability, muscle/ligament weakness or injury and degenerative joint disease where extension/flexion control is required
- Manufactured from premium grade latex-free neoprene



#### Product Number

<sup>3</sup> / <sub>16</sub> "	XS	S	M	L	XL
13"	88-0362	88-0363	88-0364	88-0365	88-0366
13" Op Pop	88-0372	88-0373	88-0374	88-0375	88-0376
<sup>3</sup> / <sub>16</sub> "	2XL	3XL	4XL		
13"	88-0367	88-0368	88-0369		

Sizing Guide

Product Number

XS

SAME AS SIZING GUIDE ON THIS PAGE



## Knee Sleeve w/Hinge

- Features sleeve design with heavy-duty hinges for optimal medial/lateral support
- Indicated for medial/lateral instability, muscle/ligament weakness or injury, and degenerative joint disease
- Manufactured from premium grade latex-free neoprene or CoolTex, a neoprene alternative





13" Std	88-0302	88-0303	88-0304	88-0305	88-0306
13" Op Pop	88-0312	88-0313	88-0314	88-0315	88-0316
<sup>1</sup> / <sub>8</sub> "	2XL	3XL	4XL		
13" Std	88-0307	88-0308	88-0309		
13" Op Pop	88-0317	88-0318	88-0319		
<sup>3</sup> / <sub>16</sub> "	XS	S	М	L	XL
13" Std	88-0322	88-0323	88-0324	88-0325	88-0326
13" Op Pop	88-0332	88-0333	88-0334	88-0335	88-0336
<sup>3</sup> / <sub>16</sub> "	2XL	3XL	4XL		
13" Std	88-0327	88-0328	88-0329		
13" Op Pop	88-0337	88-0338	88-0339		
CoolTex	XS	S	M	L	XL
13" Std	86-0322	86-0323	86-0324	86-0325	86-0326

## Sizing Guide

CoolTex 13" Std 2XL

86-0327

13" Op Pop 86-0337 86-0338 86-0339

3XL

4XL

86-0328 86-0329

	6" ABOVE MID-PATELLA	MID-PATELLA	6" BELOW MID-PATELLA
X-Small	16"- 18"	12"- 13"	11"- 12"
Small	17"- 19"	13"- 14"	12"- 13"
Medium	18"- 20"	14"- 15"	13"- 14"
Large	19"- 21"	15"- 16"	14"- 15"
X-Large	20"- 22"	16"- 18"	15"- 16"
2X-Large	21"- 23"	18"- 20"	16"- 17"
3X-Large	23"- 26"	20"- 22"	17"- 18"
4X-Large	26"- 29"	22"- 24"	18"- 19"

XL

## U-Shaped Patella Stabilizer w/Hinge

- Features sewn-in U-shaped solid tubular buttress for patella stabilization with heavy-duty hinges for optimal medial/lateral support
- Indicated for medial/lateral instability, muscle/ligament weakness or injury, treatment of patella baja and degenerative joint disease where patella stabilization is also required
- Manufactured from premium grade latex-free neoprene or CoolTex, a neoprene alternative





Product Nur	nbei			
<sup>1</sup> / <sub>8</sub> "	XS	S	M	L
13"	88-5112	88-5113	88-5114	88-5115
13" Op Pop	88-5122	88-5123	88-5124	88-5125
1/8"	XL	2XL	3XL	4XL
13"	88-5116	88-5117	88-5118	88-5119
13" Op Pop	88-5126	88-5127	88-5128	88-5129
<sup>3</sup> / <sub>16</sub> "	XS	S	M	L
13"	88-5152	88-5153	88-5154	88-5155
13" Op Pop	88-5162	88-5163	88-5164	88-5165
<sup>3</sup> / <sub>16</sub> "	XL	2X	3XL	4XL
<sup>3</sup> / <sub>16</sub> " 13"	<b>XL</b> 88-5156	<b>2X</b> 88-5157	<b>3XL</b> 88-5158	<b>4XL</b> 88-5159
	88-5156			
13"	88-5156	88-5157	88-5158	88-5159
13" 13" Op Pop	88-5156 88-5166	88-5157 88-5167	88-5158 88-5168	88-5159 88-5169
13" 13" Op Pop <b>CoolTex</b>	88-5156 88-5166 <b>XS</b> 86-5152	88-5157 88-5167 <b>S</b>	88-5158 88-5168 <b>M</b>	88-5159 88-5169 <b>L</b>
13" 13" Op Pop <b>CoolTex</b> 13"	88-5156 88-5166 <b>XS</b> 86-5152	88-5157 88-5167 <b>S</b> 86-5153	88-5158 88-5168 <b>M</b> 86-5154	88-5159 88-5169 <b>L</b> 86-5155
13" 13" Op Pop <b>CoolTex</b> 13" 13"Op Pop	88-5156 88-5166 <b>XS</b> 86-5152 86-5162	88-5157 88-5167 <b>S</b> 86-5153 86-5163	88-5158 88-5168 <b>M</b> 86-5154 86-5164	88-5159 88-5169 <b>L</b> 86-5155 86-5165



#### SAME AS SIZING GUIDE ON THIS PAGE

# CoolTex Hinged Knee Sleeve w/Lo Pro Hinge

- Features sleeve design with low profile hinge for medial/lateral support
- Indicated for medial/lateral instability, muscle/ligament weakness or injury, and degenerative joint disease
- Manufactured from lightweight, breathable CoolTex material







CoolTex	XS	S	M	L
13"	86-0382	86-0383	86-0384	86-0385
13" Op Pop	86-0392	86-0393	86-0394	86-0395
CoolTex	XL	2XL	3XL	4XL
CoolTex 13"	/ <del>-</del>	<b>2XL</b> 86-0387	U/ 11_	.,

	6" above Mid-patella	MID-PATELLA	6" BELOW MID-PATELLA
X-Small	16"- 18"	12"- 13"	11"- 12"
Small	17"- 19"	13"- 14"	12"- 13"
Medium	18"- 20"	14"- 15"	13"- 14"
Large	19"- 21"	15"- 16"	14"- 15"
X-Large	20"- 22"	16"- 18"	15"- 16"
2X-Large	21"- 23"	18"- 20"	16"- 17"
3X-Large	23"- 26"	20"- 22"	17"- 18"
4X-Large	26"- 29"	22"- 24"	18"- 19"






## Patella Stabilizer w/Cor-Trak Buttress & Hinge


- Features movable U-shaped Cor-Trak buttress with textured finish for maximum patella stabilization and heavy-duty hinges for optimal medial/ lateral support
- Indicated for medial/lateral instability, muscle/ligament weakness or injury and degenerative joint disease where patella stabilization is also required
- Manufactured from premium grade latex-free neoprene and compression molded textured neoprene buttress

Product Nur	mber			
<sup>1</sup> / <sub>8</sub> "	XS	S	M	L
13"	88-5032	88-5033	88-5034	88-5035
13" Op Pop	88-5042	88-5043	88-5044	88-5045
<sup>1</sup> / <sub>8</sub> "	XL	2X	3XL	4XL
13"	88-5036	88-5037	88-5038	88-5039
13" Op Pop	88-5046	88-5047	88-5048	88-5049
<sup>3</sup> / <sub>16</sub> "	XS	S	M	L
13"	88-5072	88-5073	88-5074	88-5075
13" Op Pop	88-5082	88-5083	88-5084	88-5085
<sup>3</sup> / <sub>16</sub> "	XL	2XL	3XL	4XL
13"	88-5076	88-5077	88-5078	88-5079
13" Op Pop	88-5086	88-5087	88-5088	88-5089

#### Sizing Guide

9			
	6" ABOVE MID-PATELLA	MID-PATELLA	6" BELOW MID-PATELLA
X-Small	16"- 18"	12"- 13"	11"- 12"
Small	17"- 19"	13"- 14"	12"- 13"
Medium	18"- 20"	14"- 15"	13"- 14"
Large	19"- 21"	15"- 16"	14"- 15"
X-Large	20"- 22"	16"- 18"	15"- 16"
2X-Large	21"- 23"	18"- 20"	16"- 17"
3X-Large	23"- 26"	20"- 22"	17"- 18"
4X-Large	26"- 29"	22"- 24"	18"- 19"





- Indicated for patella subluxation, patello-femoral syndrome and as part of a rehabilitation program for patella tendon injuries where medial/lateral support is also required
- Manufactured from premium grade latex-free neoprene and compression molded textured neoprene buttress



<sup>3</sup> / <sub>16</sub> "	XS	S	M	L	XL
13"	88-7272	88-7273	88-7274	88-7275	88-7276
<sup>3</sup> / <sub>16</sub> "	2XL	3XL	4XL		

	6" above Mid-patella	MID-PATELLA	6" BELOW MID-PATELLA
X-Small	16"- 18"	12"- 13"	11"- 12"
Small	17"- 19"	13"- 14"	12"- 13"
Medium	18"- 20"	14"- 15"	13"- 14"
Large	19"- 21"	15"- 16"	14"- 15"
X-Large	20"- 22"	16"- 18"	15"- 16"
2X-Large	21"- 23"	18"- 20"	16"- 17"
3X-Large	23"- 26"	20"- 22"	17"- 18"
4X-Large	26"- 29"	22"- 24"	18"- 19"



## Kinetic Posterior Adj Knee Sleeve w/Cor-Trak Buttress & Hinge

- → Features dynamic buttress for adjustable tension control with unique design that offers the convenience of a sleeve below the knee and benefits of an adjustable wrap above the knee
- Indicated for patella subluxation, patello-femoral syndrome and as part of a rehab program for patella tendon injuries where medial/lateral support is required
- Manufactured from premium grade latex-free neoprene and compression molded textured neoprene buttress



Product Nu	mber			
<sup>3</sup> / <sub>16</sub> "	S	M	L	XL
13"	88-8153	88-8154	88-8155	88-8156
<sup>3</sup> / <sub>16</sub> "	2XL	3XL	4XL	
13"	88-8157	88-8158	88-8159	

#### Sizing Guide

	6" above Mid-patella	MID-PATELLA	6" BELOW MID-PATELLA
Small	16"- 18"	12"- 14"	11"- 13"
Medium	18"- 19"	14"- 16"	13"- 15"
Large	19"- 22"	16"- 18"	15"- 17"
X-Large	22"- 24"	18"- 20"	17"- 19
2X-Large	24"- 26"	20"- 22"	19"- 21"
3X-Large	26"- 28"	22"- 24"	21"- 23"
4X-Large	28"- 30"	24"- 26"	23"- 24"



- Features sewn-in J-shaped solid tubular buttress with adjustable pressure strap that controls force applied to patella; low profile hinge
- Indicated for medial/lateral instability, patella instability, lateral patella subluxation and dislocation
- Manufactured from lightweight, breathable CoolTex material





1 TOGGCE TAGITIDE				
CoolTex	XS	S	М	L
13" Right	86-6182	86-6183	86-6184	86-6185
13" Left	86-6192	86-6193	86-6194	86-6195
CoolTex	XL	2XL	3XL	4XL
CoolTex 13" Right	<b>XL</b> 86-6186		<b>3XL</b> 86-6188	.,

	6" above Mid-patella	MID-PATELLA	6" BELOW MID-PATELLA
X-Small	16"- 18"	12"- 13"	11"- 12"
Small	17"- 19"	13"- 14"	12"- 13"
Medium	18"- 20"	14"- 15"	13"- 14"
Large	19"- 21"	15"- 16"	14"- 15"
X-Large	20"- 22"	16"- 18"	15"- 16"
2X-Large	21"- 23"	18"- 20"	16"- 17"
3X-Large	23"- 26"	20"- 22"	17"- 18"
4X-Large	26"- 29"	22"- 24"	18"- 19"
		20 22	.,





## 6" Knee-O-Trakker w/Hinge

- Features textured C-shaped buttress for maximum patella control, versatile strapping system that provides dynamic treatment and hinges for medial/lateral support
- Indicated for patella subluxation, patello-femoral syndrome and as part of a rehabilitation for patella tendon injuries where medial/lateral support is required
- Manufactured from premium grade latex-free neoprene, or CoolTex, a neoprene alternative; compression molded textured neoprene buttress and high density plastic hinges



6"	XS	S	M	L
<sup>3</sup> / <sub>16</sub> "	88-7242	88-7243	88-7244	88-7245
CoolTex*	86-7242	86-7243	86-7244	86-7245
6"	XL	2XL	3XL	4XL
<b>6"</b> <sup>3</sup> / <sub>16</sub> "		_,	<b>3XL</b> 88-7248	.,

Sizing Guide

XS	S	M	L
12"- 13"	13"- 14"	14"- 15"	15"- 16"
XL	2XL	3XL	4XL
16"- 18"	18"- 20"	20"- 22"	22"- 24"

Measurement taken at mid-patella

\*Cor-Trak buttress contains neoprene



#### 6" Knee-O-Trakker

- Features textured C-shaped buttress for maximum patella control and versatile strapping system that provides dynamic treatment with adjustable tension control
- Indicated for patella subluxation, patello-femoral syndrome and as part of a rehabilitation program for patella tendon injuries
- Manufactured from premium grade latex-free neoprene or CoolTex, a neoprene alternative; and compression molded textured neoprene buttress



#### Product Number

6"	XS	S	M	L
<sup>3</sup> / <sub>16</sub> "	88-7232	88-7233	88-7234	88-7235
CoolTe	x* 86-7232	86-7233	86-7234	86-7235
6"	XL	2XL	3XL	4XL
<b>6"</b> <sup>3</sup> / <sub>16</sub> "		<b>2XL</b> 88-7237		

#### Sizing Guide

XS	S	M	L
12"- 13"	13"- 14"	14"- 15"	15"- 16"
XL	2XL	3XL	4XL
16"- 18"	18"- 20"	20"- 22"	22"- 24"

Measurement taken at mid-patella

\*Cor-Trak buttress contains neoprene



#### 6" Knee-O-Trakker w/Tilt Strap

- Features textured C-shaped buttress for maximum patella control with tilt strap to assist in proper tracking and versatile strapping system that provides dynamic treatment with adjustable tension control
- Indicated for patella subluxation, patello-femoral syndrome and as part of a rehabilitation program for patella tendon injuries
- Manufactured from premium grade latex-free neoprene and compression molded textured neoprene buttress

Product Number

6"	XS	S	M	L
<sup>3</sup> / <sub>16</sub> "	88-7252	88-7253	88-7254	88-7255
6"	XL	2XL	3XL	4XL

## Sizing Guide

XS	S	M	L
12"- 13"	13"- 14"	14"- 15"	15"- 16"
XL	2XL	3XL	4XL

Measurement taken at mid-patella



## **Knee-O-Trakker w/Stays**

- Features textured C-shaped buttress for maximum patella control and versatile strapping system that provides dynamic treatment with adjustable tension control with spiral stays to prevent roll over
- Indicated for patella subluxation, patello-femoral syndrome and as part of a rehabilitation program for patella tendon injuries
- Manufactured from premium grade latex-free neoprene and compression molded textured neoprene buttress



Product Number				
13"	XS	S	М	L
<sup>3</sup> / <sub>16</sub> "	88-7262	88-7263	88-7264	88-7265
13"	XL	2XL	3XL	4XL
<sup>3</sup> / <sub>16</sub> "	88-7266	88-7267	88-7268	88-7269

#### Sizing Guide

#### SAME AS SIZING GUIDE ON THIS PAGE



## **U-Shaped Patella Stabilizer**

- → Features sewn-in U-shaped solid tubular buttress to provide optimal contact and patella stability with spiral stays to prevent rollover; some styles available with superior buttress placement
- Indicated for patella subluxation, patello-femoral tracking syndrome, treatment of patella alta/baja and as part of a rehabilitation program for patella tendon injuries
- Manufactured from premium grade latex-free neoprene or Cooltex, a neoprene alternative





Product Number				
13"	XS	S	М	L
<sup>1</sup> / <sub>8</sub> " Inferior U	88-5092	88-5093	88-5094	88-5095
1/8" Inferior U Op Pop	88-5102	88-5103	88-5104	88-5105
13"	XL	2XL	3XL	4XL
<sup>1</sup> / <sub>8</sub> " Inferior U	88-5096	88-5097	88-5098	88-5099
<sup>1</sup> / <sub>8</sub> " Inferior U Op Pop	88-5106	88-5107	88-5108	88-5109
13"	XS	S	M	L
<sup>3</sup> / <sub>16</sub> " Inferior U	88-5132	88-5133	88-5134	88-5135
<sup>3</sup> / <sub>16</sub> " Inferior U Op Pop	88-5142	88-5143	88-5144	88-5145
13"	XL	2XL	3XL	4XL
<sup>3</sup> / <sub>16</sub> " Inferior U	88-5136	88-5137	88-5138	88-5139
<sup>3</sup> / <sub>16</sub> " Inferior U Op Pop	88-5146	88-5147	88-5148	88-5149
CoolTex	XS	S	M	L
13" Inferior U	86-5132	86-5133	86-5134	86-5135
Inferior U Op Pop	86-5142	86-5143	86-5144	86-5145
CoolTex	XL	2XL	3XL	4XL
13" Inferior U	86-5136	86-5137	86-5138	86-5139
Inferior U Op Pop	86-5146	86-5147	86-5148	86-5149

	6" ABOVE MID-PATELLA	MID-PATELLA	6" BELOW MID-PATELLA
X-Small	16"- 18"	12"- 13"	11"- 12"
Small	17"- 19"	13"- 15"	12"- 13"
Medium	18"-20"	14"- 15"	13"- 14"
Large	19"- 21"	15"- 16"	14"- 15"
X-Large	20"- 22"	16"- 18"	15"- 16"
2X-Large	21"- 23"	18"- 20"	16"- 17"
3X-Large	23"- 26"	20"- 22"	17"- 18"
4X-Large	26"- 29"	22"- 24"	18"- 19"




88-5015

88-5025

88-5069



- Features movable U-shaped Cor-Trak buttress with textured finish for maximum patella stabilization and spiral stays for moderate medial/ lateral support
- Indicated for patella subluxation, patello-femoral tracking syndrome, and as part of a rehabilitation program for patella tendon injuries
- Manufactured from premium grade latex-free neoprene and compression molded textured neoprene buttress



13"	XL	2XL	3XL	4XL
1/8"	88-5016	88-5017	88-5018	88-5019
<sup>1</sup> / <sub>8</sub> " Op Pop	88-5026	88-5027	88-5028	88-5029
13"	XS	S	М	L
<sup>3</sup> / <sub>16</sub> "	88-5052	88-5053	88-5054	88-5055
<sup>3</sup> / <sub>16</sub> " Op Pop	88-5062	88-5063	88-5064	88-5065
13"	XL	2XL	3XL	4XL
3/ <sub>16</sub> "	88-5056	88-5057	88-5058	88-5059

88-5067

88-5068

#### Sizing Guide


<sup>3</sup>/<sub>16</sub>" Op Pop 88-5066

Product Number

	6" Above Mid-Patella	MID-PATELLA	6" BELOW MID-PATELLA
X-Small	16"- 18"	12"- 13"	11"-12"
Small	17"- 19"	13"- 15"	12"- 13"
Medium	18"- 20"	14"- 15"	13"- 14"
Large	19"- 21"	15"- 16"	14"- 15"
X-Large	20"- 22"	16"- 18"	15"- 16"
2X-Large	21"- 23"	18"- 20"	16"- 17"
3X-Large	23"- 26"	20"- 22"	17"- 18"
4X-Large	26"- 29"	22"- 24"	18"- 19"



- Unique design offers the convenience of a sleeve below the knee and benefits of an adjustable wrap above the knee with moveable U-shaped Cor-Trak buttress and spiral stays
- Indicated for patella subluxation, patello-femoral tracking syndrome and as part of a rehabilitation program for patella tendon injuries
- Manufactured from premium grade latex-free and neoprene compression molded textured neoprene buttress



	6" ABOVE MID-PATELLA	MID-PATELLA	6" BELOW MID-PATELLA
Small	16"- 18"	12"- 14"	11"- 12"
Medium	18"- 19"	14"- 16"	13"- 14"
Large	19"- 22"	16"- 18"	15"- 16"
X-Large	22"- 24"	18"- 20"	17"- 18"
2X-Large	24"- 26"	20"- 22"	19"- 20"
3X-Large	26"- 28"	22"- 24"	21"- 22"
4X-Large	28"- 30"	24"- 26"	23"- 24"

## J-Pull Patella Stabilizer

- Features sewn-in J-shaped solid tubular buttress with adjustable pressure strap that controls force applied to patella; spiral stays
- Indicated for patella instability, lateral patella subluxation and dislocation
- Manufactured from premium grade latex-free neoprene



Product	Number
3/16"	
7/16	

<sup>3</sup> / <sub>16</sub> "	XS	S	М	L
13" Right	88-6052	88-6053	88-6054	88-6055
13" Left	88-6062	88-6063	88-6064	88-6065
<sup>3</sup> / <sub>16</sub> "	XL	2XL	3XL	4XL
13" Right	88-6056	88-6057	88-6058	88-6059
13" Left	88-6066	88-6067	88-6068	88-6069
<sup>3</sup> / <sub>16</sub> "	XS	S	M	L
3/ <sub>16</sub> " 13" Rt Op Pop	<b>XS</b> 88-6072	<b>S</b> 88-6073	<b>M</b> 88-6074	<b>L</b> 88-6075
		-		_
13" Rt Op Pop	88-6072	88-6073	88-6074	88-6075
13" Rt Op Pop 13" Lf Op Pop	88-6072 88-6082	88-6073 88-6083	88-6074 88-6084	88-6075 88-6085



J			
	6" above Mid-patella	MID-PATELLA	6" BELOW MID-PATELLA
X-Small	16"- 18"	12"- 13"	11"- 12"
Small	17"- 19"	13"- 14"	12"- 13"
Medium	18"- 20"	14"- 15"	13"- 14"
Large	19"- 21"	15"- 16"	14"- 15"
X-Large	20"- 22"	16"- 18"	15"- 16"
2X-Large	21"- 23"	18"- 20"	16"- 17"
3X-Large	23"- 26"	20"- 22"	17"- 18"
4X-Large	26"- 29"	22"- 24"	18"- 19"



## **C-Pull Patella Stabilizer**

- Features sewn-in C-shaped solid tubular buttress with adjustable pressure strap that controls pressure applied to patella and spiral stays
- Indicated for patella instability, lateral patella subluxation and dislocation
- Manufactured from premium grade latex-free neoprene



#### Product Number

<sup>3</sup> / <sub>16</sub> "	XS	S	M	L
13" Right	88-7052	88-7053	88-7054	88-7055
13" Left	88-7062	88-7063	88-7064	88-7065
<sup>3</sup> / <sub>16</sub> "	XL	2XL	3XL	4XL
13" Right	88-7056	88-7057	88-7058	88-7059
13" Left	88-7066	88-7067	88-7068	88-7069
<sup>3</sup> / <sub>16</sub> "	XS	S	М	L
³/ <sub>16</sub> " 13" Rt Op Pop	<b>XS</b> 88-7072	<b>S</b> 88-7073	<b>M</b> 88-7074	<b>L</b> 88-7075
0	-	_		_
13" Rt Op Pop	88-7072	88-7073	88-7074	88-7075
13" Rt Op Pop 13" Lf Op Pop	88-7072 88-7082	88-7073 88-7083	88-7074 88-7084	88-7075 88-7085

_			
	6" ABOVE MID-PATELLA	MID-PATELLA	6" BELOW MID-PATELLA
X-Small	16"- 18"	12"- 13"	11"- 12"
Small	17"- 19"	13"- 14"	12"- 13"
Medium	18"- 20"	14"- 15"	13"- 14"
Large	19"- 21"	15"- 16"	14"- 15"
X-Large	20"- 22"	16"- 18"	15"- 16"
2X-Large	21"- 23"	18"- 20"	16"- 17"
3X-Large	23"- 26"	20"- 22"	17"- 18"
4X-Large	26"- 29"	22"- 24"	18"- 19"





## Kinetic Posterior Adjustable Knee Sleeve w/Cor-Trak Buttress

- Features dynamic buttress for adjustable tension control with unique design that offers the convenience of a sleeve below the knee and benefits of an adjustable wrap above the knee
- Indicated for patella subluxation, patello-femoral tracking syndrome and correction of patella tilt
- Manufactured from premium grade latex-free neoprene and compression molded textured neoprene buttress

•	_	•	`
۷	_		
2	n		١

Product	Number			
13"	S	М	L	XL
<sup>3</sup> / <sub>16</sub> "	88-8143	88-8144	88-8145	88-8146
13"	2XL	3XL	4XL	
<sup>3</sup> / <sub>16</sub> "	88-8147	88-8148	88-8149	

#### Sizing Guide

0121119 04	100		
	6" ABOVE MID-PATELLA	MID-PATELLA	6" BELOW MID-PATELLA
Small	16"- 18"	12"- 14"	11"- 13"
Medium	18"- 19"	14"- 16"	13"- 15"
Large	19"- 22"	16"- 18"	15"- 17"
X-Large	22"- 24"	18"- 20"	17"- 19"
2X-Large	24"- 26"	20"- 22"	19"-21"
3X-Large	26"- 28"	22"- 24"	21"- 23"
4X-Large	28"- 30"	24"- 26"	23"- 24"



## Medial-Lateral Patella Stabilizer

- Features sewn-in medial-lateral tubular buttress to provide optimal contact and patella stability with spiral stays to prevent roll over
- Indicated for patella subluxation, patello-femoral tracking syndrome and as part of a rehabilitation program for patella tendon injuries
- Manufactured from premium grade latex-free neoprene



#### Product Number

13"	XS	S	M	L
1/8"	88-8172	88-8173	88-8174	88-8175
3/16"	88-8182	88-8183	88-8184	88-8185
13"	XL	2XL	3XL	4XL
13"  1/8"		<b>2XL</b> 88-8177		<b>4XL</b> 88-8179

## Sizing Guide

	6" ABOVE MID-PATELLA	MID-PATELLA	6" BELOW MID-PATELLA
X-Small	16"- 18"	12"- 13"	11"- 12"
Small	17"- 19"	13"- 14"	12"- 13"
Medium	18"- 20"	14"- 15"	13"- 14"
Large	19"- 21"	15"- 16"	14"- 15"
X-Large	20"- 22"	16"- 18"	15"- 16"
2X-Large	21"- 23"	18"- 20"	16"- 17"
3X-Large	23"- 26"	20"- 22"	17"- 18"
4X-Large	26"- 29"	22"- 24"	18"- 19"



This symbol indicates that pediatric sizes are available for this product.



This symbol indicates that bariatric sizes are available for this product.



This symbol certifies that all products manufactured by Corflex are latex free.



This symbol indicates that certain precautions are needed when using this product.



This symbol identifies CoolTex products that are manufactured from a neoprene-alternative material.

65

# Osgood-Schlatter Knee Sleeve

- Features patella pad with circumferential strap to provide compression over affected area
- Indicated for treatment of Osgood Schlatter disease and Jumper's Knee
- Manufactured from premium grade latex-free neoprene



## Product Number

13"	XS	S	M	L	XL
1/8"	88-0402	88-0403	88-0404	88-0405	88-0406
13"	2XL	3XL	4XL		

#### Sizing Guide

	6" ABOVE MID-PATELLA	MID-PATELLA	6" BELOW MID-PATELLA
X-Small	16"- 18"	12"- 13"	11"- 12"
Small	17"- 19"	13"- 14"	12"- 13"
Medium	18"- 20"	14"- 15"	13"- 14"
Large	19"- 21"	15"- 16"	14"- 15"
X-Large	20"- 22"	16"- 18"	15"-16"
2X-Large	21"- 23"	18"- 20"	16"- 17"
3X-Large	23"- 26"	20"- 22"	17"- 18"
4X-Large	26"- 29"	22"- 24"	18"- 19"



## **Knee-O-Band**

- Features solid tubular neoprene pad that applies pressure to the patella tendon
- Indicated for Osgood Schlatter disease, patello-femoral syndrome and patella tendonitis
- Manufactured from premium grade latex-free soft loop neoprene



#### Product Number

	3XS	2XS	XS	S	M
<sup>3</sup> / <sub>16</sub> "	88-7210	88-7211	88-7212	88-7213	88-7214
	L	XL	2XL		

#### Sizing Guide

3XS	2XS	XS	S	M
10"-11"	11"-12"	12"- 13"	13"- 14"	14"- 15"
L	XL	2XL		

Measurement taken at mid-patella



## **Kneed It**

- Unique v-notch design applies gentle pressure to the patella tendon
- Ideal for relief of minor knee pain associated with arthritis, tendonitis and chrondromalacia
- Constructed of two-piece rigid plastic with neoprene pad and contact closure strap

## Product Number

UNIV
99-560

Sizing Guide

UNIV

13"- 20"

Measurement taken at mid-patella





# Knee Sleeve w/Stays

- → 13" support provides warmth, support and compression for the knee with spiral stays to prevent roll over
- Indicated for moderate knee strains or sprains, knee injuries and arthritis
- Manufactured from premium grade latex-free neoprene



## Product Number

<sup>1</sup> / <sub>8</sub> "	XS	S	M	L
	88-0052	88-0053	88-0054	88-0055
Ор Рор	88-0062	88-0063	88-0064	88-0065
Op Pat	88-0072	88-0073	88-0074	88-0075
Op Pat/Pop	88-0082	88-0083	88-0084	88-0085
1/8"	XL	2XL	3XL	4XL
	88-0056	88-0057	88-0058	88-0059
Ор Рор	88-0066	88-0067	88-0068	88-0069
Op Pat	88-0076	88-0077	88-0078	88-0079
Op Pat/Pop	88-0086	88-0087	88-0088	88-0089
3/16"	XS	S	M	L
	88-0132	88-0133	88-0134	88-0135
Ор Рор	88-0132 88-0142	88-0133 88-0143	88-0134 88-0144	88-0135 88-0145
Op Pop Op Pat				
	88-0142	88-0143	88-0144	88-0145
Op Pat	88-0142 88-0152	88-0143 88-0153	88-0144 88-0154	88-0145 88-0155
Op Pat Op Pat/Pop	88-0142 88-0152 88-0162	88-0143 88-0153 88-0163	88-0144 88-0154 88-0164	88-0145 88-0155 88-0165
Op Pat Op Pat/Pop	88-0142 88-0152 88-0162 <b>XL</b>	88-0143 88-0153 88-0163 <b>2XL</b>	88-0144 88-0154 88-0164 <b>3XL</b>	88-0145 88-0155 88-0165 <b>4XL</b>
Op Pat Op Pat/Pop <sup>3</sup> / <sub>16</sub> "	88-0142 88-0152 88-0162 <b>XL</b> 88-0136	88-0143 88-0153 88-0163 <b>2XL</b> 88-0137	88-0144 88-0154 88-0164 <b>3XL</b> 88-0138	88-0145 88-0155 88-0165 <b>4XL</b> 88-0139

# **Sizing Chart For All Knee Sleeves**

Sizing Garac			
	6" ABOVE MID-PATELLA	MID-PATELLA	6" BELOW MID-PATELLA
X-Small	16"-18"	12"-13"	11"-12"
Small	17"-19"	13"-14"	12"-13"
Medium	18"-20"	14"-15"	13"-14"
Large	19"-21"	15"-16"	14"-15"
X-Large	20"-22"	16"-18"	15"-16"
2X-Large	21"-23"	18"-20"	16"-17"
3X-Large	23"-26"	20"-22"	17"-18"
4X-Large	26"-29"	22"-24"	18"-19"

## **Knee Sleeve**

- → 13" support provides warmth, support and compression for the knee
- Indicated for moderate knee strains or sprains, knee injuries and arthritis
- Manufactured from premium grade latex-free neoprene



#### Product Number

1/,"

1/8"	XS	5	IVI	L
	88-0012	88-0013	88-0014	88-0015
Ор Рор	88-0022	88-0023	88-0024	88-0025
Op Pat	88-0032	88-0033	88-0034	88-0035
Op Pat/Pop	88-0042	88-0043	88-0044	88-0045
1/8"	XL	2XL	3XL	4XL
	88-0016	88-0017	88-0018	88-0019
Ор Рор	88-0026	88-0027	88-0028	88-0029
Op Pat	88-0036	88-0037	88-0038	88-0039
Op Pat/Pop	88-0046	88-0047	88-0048	88-0049
<sup>3</sup> / <sub>16</sub> "	XS	S	М	L
	88-0092	88-0093	88-0094	88-0095
Op Pop	88-0102	88-0103	88-0104	88-0105
Op Pop Op Pat	88-0102 88-0112	88-0103 88-0113	88-0104 88-0114	88-0105 88-0115
Op Pat	88-0112	88-0113	88-0114	88-0115
Op Pat Op Pat/Pop	88-0112 88-0122	88-0113 88-0123	88-0114 88-0124	88-0115 88-0125
Op Pat Op Pat/Pop	88-0112 88-0122 <b>XL</b>	88-0113 88-0123 <b>2XL</b>	88-0114 88-0124 <b>3XL</b>	88-0115 88-0125 <b>4XL</b>
Op Pat Op Pat/Pop <sup>3</sup> / <sub>16</sub> "	88-0112 88-0122 <b>XL</b> 88-0096	88-0113 88-0123 <b>2XL</b> 88-0097	88-0114 88-0124 <b>3XL</b> 88-0098	88-0115 88-0125 <b>4XL</b> 88-0099



## Knee Sleeve w/Anterior Pad

- → 13" support provides warmth, protection and compression for the knee, also features anterior pad over patella region for added protection
- Indicated for moderate knee strains or sprains, knee injuries and arthritis
- Manufactured from premium grade latex-free neoprene



#### Product Number

	Λ3	3	IVI	-
	88-0212	88-0213	88-0214	88-0215
Ор Рор	88-0242	88-0243	88-0244	88-0245
Op Pat	88-0232	88-0233	88-0234	88-0235
Op Pat/Pop	88-0222	88-0223	88-0224	88-0225
1/8"	XL	2XL	3XL	4XL
	88-0216	88-0217	88-0218	88-0219
Ор Рор	88-0246	88-0247	88-0248	88-0249
Op Pat	88-0236	88-0237	88-0238	88-0239
Op Pat/Pop	88-0226	88-0227	88-0228	88-0229
<sup>3</sup> / <sub>16</sub> "	XS	S	М	L
	88-0172	88-0173	88-0174	88-0175
Ор Рор	88-0172 88-0202	88-0173 88-0203	88-0174 88-0204	88-0175 88-0205
Op Pop Op Pat				
	88-0202	88-0203	88-0204	88-0205
Op Pat	88-0202 88-0192	88-0203 88-0193	88-0204 88-0194	88-0205 88-0195
Op Pat Op Pat/Pop	88-0202 88-0192 88-0182	88-0203 88-0193 88-0183	88-0204 88-0194 88-0184	88-0205 88-0195 88-0185
Op Pat Op Pat/Pop	88-0202 88-0192 88-0182 <b>XL</b>	88-0203 88-0193 88-0183 <b>2XL</b>	88-0204 88-0194 88-0184 <b>3XL</b>	88-0205 88-0195 88-0185 <b>4XL</b>
Op Pat Op Pat/Pop 3/ <sub>16</sub> "	88-0202 88-0192 88-0182 <b>XL</b> 88-0176	88-0203 88-0193 88-0183 <b>2XL</b> 88-0177	88-0204 88-0194 88-0184 <b>3XL</b> 88-0178	88-0205 88-0195 88-0185 <b>4XL</b> 88-0179





## **Universal Knee Wrap**

- Provides warmth, support, and compression for the knee, universal design allows for inventory reduction; available with spiral stays, hinges or dynamic buttress
- Indicated for moderate knee strains or sprains, knee injuries and arthritis
- Manufactured from premium grade latex-free soft brushed loop neoprene and compression molded textured neoprene buttress

Product Num
-------------

111/2"		UNIV
¹/8″ Ultra		88-1000
<sup>1</sup> / <sub>8</sub> " Ultra w/Stays		88-1001
<sup>1</sup> / <sub>8</sub> " Wrap		88-1010
1/8" Wrap w/Stays		88-1020
<sup>1</sup> / <sub>8</sub> " Wrap w/Hing	е	88-1030
Sizing Guide		
	ULTRA UNIV	WRAP UNIV

	ULTRA UNIV	WRAP U
Fits up to	25"	22"

Measurement taken 6" above mid-patella



## **Anterior Closure Knee Wrap w/Stays**

- Features wrap-around design for ease of application and provides warmth, support and compression for the knee with spiral stays to prevent roll over
- Indicated for moderate knee strains or sprains, knee injuries and arthritis
- Manufactured from premium grade latex-free neoprene



## Product Number

13"	S	М	L	XL
3/16"	88-1373	88-1374	88-1375	88-1376
<sup>3</sup> / <sub>16</sub> " Op Pop	88-1363	88-1364	88-1365	88-1366
13"	2XL	3XL	4XL	
<b>13"</b> <sup>3</sup> / <sub>16</sub> "		<b>3XL</b> 88-1378		

Sizing Guide

SAME AS SIZING GUIDE BELOW



## **Anterior Closure Knee Wrap**

- Features wrap-around design for ease of application and provides warmth, support and compression for the knee
- Indicated for moderate knee strains or sprains, knee injuries and arthritis
- Manufactured from premium grade latex-free neoprene



# Product Number

13	3	IVI	_	/L
<sup>3</sup> / <sub>16</sub> "	88-1313-799	88-1314-799	88-1315-799	88-1316-799
<sup>3</sup> / <sub>16</sub> " Op Pop	88-1303-799	88-1304-799	88-1305-799	88-1306-799
13"	2XL	3XL	4XL	
<b>13"</b> <sup>3</sup> / <sub>16</sub> "		<b>3XL</b> 88-1318-799	.,	

#### Sizing Guide

	6" ABOVE MID-PATELLA	MID-PATELLA	6" BELOW MID-PATELLA
Small	16"- 18"	12"- 14"	11"- 12"
Medium	18"- 19"	14"- 16"	13"- 14"
Large	19"- 21"	16"- 18"	15"- 16"
X-Large	21"- 23"	18"- 20"	17"- 18"
2X-Large	23"- 25"	20"- 22"	19"- 20"
3X-Large	25"- 27"	22"- 24"	21"- 22"
4X-Large	27"- 28"	24"- 26"	22"- 23"

CORFLEX

800 426 7353

## **Universal Knee Wrap w/Dynamic Buttress**

- Features dynamic Cor-Trak buttress for adjustable tension control; universal design allows for inventory reduction; available with either spiral stays or low profile hinge
- Indicated for patella subluxation, patella-femoral tracking syndrome and correction of patella tilt
- Manufactured from premium grade latex-free soft brushed loop neoprene and compression molded textured neoprene buttress

Product Number

1	1¹/₂"	UNIV
٧	Vrap w/Dynamic Buttress & Stays	88-1040
٧	Vrap w/Dynamic Buttress & Hinge	88-1050

Sizing Guide

CIVIV	
Fits up to	22"

LININ

Measurement taken 6" above mid-patella



## **Thigh Sleeve**

- Provides warmth, support and compression for the thigh and groin area
- Indicated for strains of the groin, hamstring or quadriceps muscles
- Manufactured from premium grade latex-free neoprene

Product Number

11"	XS	S	M	L
1/8"	88-1502	88-1503	88-1504	88-1505
3/ <sub>16</sub> "	88-1512	88-1513	88-1514	88-1515
11"	XL			
<b>11"</b>	<b>XL</b> 88-1506			

Sizing Guide

XS	S	М	L
18"- 19"	19"- 20"	20"- 22"	22"- 25"
XL			
25"- 28"			

Measurement taken from widest part of thigh



#### **Calf Sleeve**

- Provides warmth, support and compression for the calf
- Indicated for strains of the calf or gastroc muscles and shin splints
- Manufactured from premium grade latex-free neoprene

Product Number

12"	XS	S	M	L
1/8"	88-1402	88-1403	88-1404	88-1405
3/ <sub>16</sub> "	88-1412	88-1413	88-1414	88-1415
12"	XL	2XL		
<b>12"</b>	<b>XL</b> 88-1406			

## Sizing Guide

XS	S	M	L
12"- 13"	13"- 14"	14"- 15"	15"- 17"
XL	2XL		

Measurement taken from widest part of calf





## **Pneumatic Walker**

Designed to increase stability and decrease edema, our walker features a built-in pump to achieve desired level of compression. Ideal for stable foot/ankle fractures, severe sprains, soft tissue injuries, post cast removal and post-surgical application. Cushion sole with low profile non-skid bottom promotes a natural gait. Available in two heights; both feature metal uprights.

- Easy to use air pump provides desired compression level and accommodates for swelling
- > Indicated for stable foot or ankle fractures, severe sprains or strains, soft tissue injuries, post-op and post-cast removal
- Constructed of durable malleable metal uprights, easy to use air pump, circumferential bladder and contact closure straps

#### Product Number

	S	M	L	XL
Lower Leg	79-1080	79-1081	79-1082	79-1083
Ankle	79-1090	79-1091	79-1092	79-1093

#### Sizing Guide

	MENS SIZE	WOMENS SIZE
Small	4 - 6	5 - 7
Medium	6 - 10	7 - 11
Large	10 - 12	11 - 13
X-Large	12 +	13 +

Measured by shoe size





For additional sizing, laundering or application instructions visit corflex.com/pneumaticwalker

# **Lower Leg Fixed Walker**

- Provides stability for the lower leg with low natural heel height to help reduce leg length discrepancy
- Indicated for stable foot or ankle fractures, severe sprains or strains, soft tissue injuries, post-op or post-cast removal
- Constructed of rigid malleable metal uprights, padded liner and contact closure straps





Pediatric

### Product Number

XS	S	M	L
79-1049	79-1050	79-1051	79-1052
XL			
79-1053			

### Sizing Guide

	MENS SIZE	WOMENS SIZE
X-Small	up to 5	up to 5
Small	5 - 6	5 - 7
Medium	6 - 10	7 - 11
Large	10 - 12	11 - 131/2
X-Large	12 +	14 +
Manaurad	bu abaa siza	

### Measured by shoe size

# Also Available in Pediatric Sizes

Product Number

S/M	L	XL
79-1201	79-1202	79-1203

SIZIDO	Guide
JIZIIIQ	Julue

Age	1-2 years	2 ½-6 years	6 ½-9 ½ yea
Foot Lengt	h up to 6"	6"-7 ½"	7 ½"-8 ½"



# **Ankle Fixed Walker**

- Provides stability for the ankle with low natural heel height to help reduce leg length discrepancy
- Indicated for stable foot fractures, severe sprains or strains, post-op or post-cast removal
- Constructed of rigid malleable metal uprights, padded liner and contact closure straps

### Product Number

XS	S	M	L
79-1059	79-1060	79-1061	79-1062
XL			
79-1063			

3121119 00	1100	
	MENS SIZE	WOMENS SIZE
X-Small	up to 5	up to 5
Small	5 - 6	5 - 7
Medium	6 - 10	7 - 11
Large	10 - 12	11 - 13 <sup>1</sup> / <sub>2</sub>
X-Large	12 +	14 +
Measured	by shoe size	





# R.O.M. Walker

- Adjustable hinge allows for or limits range-of-motion during healing process
- Indicated for stable foot or ankle fractures, severe sprains or strains, soft tissue injuries, post-op and post-cast removal
- Constructed of durable malleable metal uprights, easy to set hinge, padded liner and contact closure straps

### Product Number

S	M	L
79-1070	79-1071	79-1072

# Sizing Guide

	MENS SIZE	WOMENS SIZE
Small	5 - 6	5 - 7
Medium	6 - 10	7 - 11
Large	10 - 12	11 -131/2
Measured	by shoe size	



# Pneumatic R.O.M.Walker

- Combines the benefits of adjustable compression with the ability to limit range of motion during the healing process
- Indicated for stable foot fractures, severe sprains, soft tissue injuries, post-op and post-cast removal
- Constructed of durable malleable metal uprights with easy to use air pump with circumferential bladder and contact closure straps

### Product Number

S	M	L
79-1075	79-1076	79-1077

### Sizing Guide

	MENS SIZE	WOMENS SIZE	
Small	5 - 6	5 - 7	
Medium	6 - 10	7 - 11	
Large	10 - 12	11 -131/2	
Measured by shoe size			



# Tri Shell Pneumatic Walker

- Features posterior and anterior shells with dual chamber liner for optimal compression and stability
- Indicated for stable fractures, severe sprains, soft tissue injuries, post-cast removal or post-op care management
- Constructed of durable plastic shell, plush liner with medial/lateral air chambers, contact closure straps and synthetic knit socks

### Product Number

	S	М	L
Lower Leg	79-1021	79-1022	79-1023
Ankle	79-1031	79-1032	79-1033

	MENS SIZE	WOMENS SIZE
Small	up to 8	up to 81/2
Med/Lg	81/2 - 111/2	9 - 12
X-Large	12 +	121/2 +
Measured	by shoe size	

# **Ultra Stretch Plantar Fasciitis Splint**

- Features bi-lateral dorsiflexion assist straps that provide a controlled amount of dorsiflexion and removable wedge for added toe extension; available with or without slip resistant tread
- Indicated for treatment of plantar fasciitis, Achilles tendonitis, plantar flexion contractures, heel spurs and other lower leg over-use injuries
- Constructed of polyethylene shell with foam liner and padded straps with quick-release buckles; available with optional textured vinyl tread

roduct	Number
TOGGC	INGILIDO

	S	М	L	XL
	37-4191	37-4192	37-4193	37-4194
w/Tread	37-4291	37-4292	37-4293	37-4294

# Sizing Guide

	MENS SIZE	WOMENS SIZE
Small	51/2 - 61/2	5 - 6
Medium	7 - 8	61/2 - 71/2
Large	81/2 - 91/2	8 - 9
X-Large	10 - 11	91/2 - 101/2
Measured	bv shoe size	



# Padded Plantar Fasciitis Night Splint

- Features polyethylene shell with soft padded cover and straps to provide maximum patient comfort; available in 90° fixed, 5° dorsiflexion and adjustable models
- Indicated for treatment of plantar fasciitis, Achilles tendonitis, plantar flexion contractures, heel spurs and other lower leg over-use injuries
- Constructed of high density rigid shell, replaceable soft padded cover and contact closure

Product Number

	S	M	L	XL	2XL
90° Fixed	37-4161	37-4162	37-4163	37-4164	37-4165
5° DF	37-4171	37-4172	37-4173	37-4174	37-4175
Adjustable	37-4181	37-4182	37-4183	37-4184	37-4185

# Sizing Guide

	MENS SIZE	WOMENS SIZE
Small	51/2 - 61/2	5 - 6
Medium	7 - 8	61/2 - 71/2
Large	81/2 - 91/2	8 - 9
X-Large	10 - 11	91/2 - 101/2
2X-Large	111/2 +	11 +
Measured I	by shoe size	



# **Ultra Plantar Fasciitis Night Splint**

- Features polyethylene shell with soft liner and d-ring closure for ease of application and secure fit
- Indicated for treatment of plantar fasciitis, Achilles tendonitis, plantar flexion contractures, heel spurs and other lower leg over-use injuries
- Constructed of high density rigid shell with soft plush replaceable liner

Product Number

	S	М	L	XL	2XL
90° Fixed	37-4031	37-4032	37-4033	37-4034	37-4035
5° DF	37-4041	37-4042	37-4043	37-4044	37-4045
Adjustable	37-4051	37-4052	37-4053	37-4054	37-4055

0121119 0 01		
	MENS SIZE	WOMENS SIZE
Small	51/2 - 61/2	5 - 6
Medium	7 - 8	61/2 - 71/2
Large	$8^{1}/_{2}$ - $9^{1}/_{2}$	8 - 9
X-Large	10 - 11	91/2 - 101/2
2X-Large	111/2 +	11 +
Measured l	by shoe size	





# **Ultra Dorsal Night Splint**

- Simple design allows for a stable and secure fit with improved patient comfort; features memory foam pad to help eliminate any possible pressure points
- Indicated for treatment of plantar fasciitis, Achilles tendonitis, plantar flexion contractures, drop foot and stroke patients
- Constructed of high density polyethylene shell with straps and sewn-in memory foam support pad

Product Number

S/M	L/XL
37-4111	37-4112

Sizing Guide

	MENS SIZE	WOMENS SIZE	
S/M	5 - 9	6 - 10	
L/XL	91/2 - 14	10 <sup>1</sup> / <sub>2</sub> - 15	

Measured by shoe size



# **Universal Plantar Fasciitis Support**

- Anatomically designed to provide continuous support and lift by means of a durable elastic tension strap; fits left or right
- Indicated for treatment of mild to moderate plantar fasciitis; aids in reducing pain and inflammation
- Constructed of 11/2" woven elastic with high resistance level; breathable 3/8" foam; alligator hook closures

Product Number

	S/M	L/XL
Single	61-1000	61-1001
10/Pack	61-1002	61-1003

Sizing (	Guide		
	MENS SIZE	WOMENS SIZE	
S/M	up to 8	up to 8 <sup>1</sup> / <sub>2</sub>	
L/XL	81/2 +	9 +	
Measure	ed by shoe size		



### Soft AFO

- Features an ankle wrap with two interchangeable forefoot attachments, making it ideal for use with a tie shoe, slip-on shoe or barefoot
- Indicated for treatment of foot-drop, dorsiflexion injury, stroke patients or leg muscle disorders
- Manufactured from breathable foam laminate with nylon straps and center release buckle; contact closure

### Product Number

S/M	L/XL
86-4100	86-4110

# Sizing Guide

S/M 8 - 10" L/XL 9<sup>1</sup>/<sub>2</sub> - 12<sup>1</sup>/<sub>2</sub>"

Measurement taken at base of ankle



# **Target Plantar Fasciitis Soft Boot**

- Features dorsiflexion assist strap that provides continuous stretching of the plantar fascia; low profile design makes it easy to apply; optional rigid foot plate available
- Indicated for treatment of plantar fasciitis and other lower leg over use injuries
- Manufactured from premium grade latex-free neoprene with contact closure strap

### Product Number

	S	M	L	XL
Soft Boot	88-2201	88-2202	88-2203	88-2204
Foot Plate	37-6141	37-6142	37-6143	37-6144
Gauntlet	88-2133	88-2134	88-2135	88-2136

# Sizing Guide

	MENS SIZE	WOMENS SIZE
Small	3 - 5	4 - 6
Medium	6 - 7	7 - 9
Large	8 - 10	10 - 11
X-Large	11 +	11 +

Measured by shoe size

# **Ankle Stirrup**

- Universal design conforms to either ankle to provide protection and prevention of ankle inversion or eversion
- Indicated for chronic ankle instability, ankle sprains or strains and for athletic protection
- Constructed of semi-rigid shells connected with soft adjustable foot-piece

### Product Number

	BLACK	WHITE	
Air	75-1015	75-1115	
Air / Gel	75-1025	75-1125	
Memory Foam	75-1030	75-1130	

### Sizing Guide

# UNIV

Fits up to 14" ankle circumference

Measurement is circumference from base of heel around ankle



# Game On Ankle Stabilizer

- Anatomically designed to support and protect the ankle; utilizes a dual figure-8 strapping system to simulate a professional taping technique
- Indicated for treatment of moderate ankle sprains, strains, ankle instability and weakness, or for prophylactic use
- Constructed of durable ballistic nylon, ventilated spacer fabric, nylon figure-8 straps and wide elastic compression strap



### Product Number

2XS	XS	S	M
75-8020	75-8021	75-8022	75-8023
L	XL	2XL	3XL

### Sizing Guide

2XS	XS	S	M
9"-10"	10"-11"	11"-12"	12"-13"
L	XL	2XL	3XL

Measurement is circumference from base of heel around ankle. Fits either left or right foot.



# Lace-Up Ankle

- Features lace-up design, nylon figure-8 support straps and removable medial/lateral stays; also available with stabilizing strap for additional support
- Indicated for chronic ankle instability, sprains or strains and protection against inversion or eversion ankle injuries
- Constructed from durable black tricot with brushed tricot lining and removable plastic stays



# Product Number

	2/\3	7.5	9	141
	75-8000	75-8001	75-8002	75-8003
w/Strap	75-8010	75-8011	75-8012	75-8013
	L	XL		
	75-8004	75-8005		
w/Strap	75-8014	75-8015		
c· ·				

### Sizing Guide

2XS	XS	S	М
9"- 10"	10"- 11"	11"- 12"	12"- 13"
L	XL		
13"- 14"	14"- 15"		

Measurement is circumference from base of heel around ankle





# **Classic Post-Op Shoe**

- Features flat footbed with rocker-style Product Number bottom to promote a natural gait pattern with tongue closure for protection
- Indicated for post-operative foot protection and support
- Constructed of corrugated semi-rigid rubber sole with water-resistant upper and contact closure

	S	M	L	XL	2XL
Mens	69-1022	69-1023	69-1024	69-1025	69-1029
Womens	69-1026	69-1027	69-1028		

# Sizing Guide

LENGTH	S	M	L	XL	2XL
Mens	10"	11"	12"	13"	14"
Womens	9"	10"	11"		



# **Rigid Post-Op Shoe**

- > Features flat footbed with rigid sole to prevent excessive foot motion while walking
- Indicated for post-operative foot protection and support
- Constructed of wood-like sole with water-resistant upper and contact closure



### Product Number

	PEDI	S	M	L	XL
Mens		69-3122	69-3123	69-3124	69-3125
Womens	5	69-3132	69-3133	69-3134	
Pedi	69-3121				

### Sizing Guide

LENGTH PEDI	S	M	L	XL
Mens	11"	12"	13"	14"
Womens	9"	10"	11"	
Pedi 8"				



# **Open Toe Cast Sandal**

- > Features open toe and heel for comfort and breathability
- Indicated for lower extremity casted
- Constructed of non-skid, ribbed rocker sole with durable upper and contact closure



### Product Number

3XS	2XS	XS	S	М
69-3078	69-3079	69-3080	69-3081	69-3082
L	XL			
69-3083	69-3084			

# Sizing Guide

	3XS	2XS	XS	S	М
Length	4"	5″	7"	9"	11"
	L	XL			
	12"	14"			



# **Enclosed Toe Cast Boot**

- > Features enclosed toe and heel to provide complete protection and warmth for the foot
- Indicated for lower extremity casted fractures
- Constructed of non-skid, ribbed rocker sole with water-resistant vinyl upper and contact closure



# Product Number

XS	S	М	L	XL
69-3100	69-3101	69-3102	69-3103	69-3104

	XS	S	М	L	XL
Length	7"	9"	11"	12"	14"

# **All Purpose Boot**

- Features square toe design for optimal comfort with rocker bottom to promote natural gait and help reduce plantar pressure
- Ideal for post-op use; provides closed toe protection that allows dressing to be kept clean and toes dry
- Constructed of EVA insole with tricot foam upper and contact closure

Product Number

**XS S M L XL** 69-1030 69-1031 69-1032 69-1033 69-1034

Sizing Guide

	XS	S	M	L	XL
Mens	up to $5^{1}/_{2}$	6 - 8	81/2 - 10	101/2 - 12	121/2 - 14
Womens	4 - 7	$7^{1}/_{2}$ -10	$10^{1}/_{2}$ -13	131/2+	

Measured by shoe size



# **Heelwedge Shoe**

- Features square toe design for optimal comfort with modified rocker bottom that off-loads pressure of the heel and encourages smooth toe-off while maintaining stable mid-foot base
- Ideal for post-op use or following trauma involving the heel to help promote healing
- Constructed of EVA insole with tricot foam upper and contact closure


Product Number

XS	S	M	L	XL
69-1050	69-1051	69-1052	69-1053	69-1054

Sizing Guide

	XS	S	M	L	XL
Mens	up to $5^{1}/_{2}$	6 - 8	81/2-10	101/2 - 12	121/2 - 14
Womens	4 - 7	71/2 -10	101/2-13	$13^{1}/_{2}+$	

Measured by shoe size



# **Orthowedge Shoe**

- Features square toe design that acts as protective bumper with 15° dorsi-angle bottom to shift weight back to mid-foot and heel
- Ideal for post-op use or following trauma of the forefoot or treatment of ulcers/infections to help promote healing
- Constructed of EVA insole with tricot foam upper and contact closure

Product Number

XS	S	M	L	XL
69-1040	69-1041	69-1042	69-1043	69-1044

Sizing Guide

	XS	S	М	L	XL
Mens	up to $5^{1}/_{2}$	6 - 8	81/2 - 10	101/2 - 12	12 <sup>1</sup> / <sub>2</sub> - 14
Womens	4 - 7	$7^{1}/_{2}$ -10	$10^{1}/_{2} - 13$	$13^{1}/_{2}+$	

Measured by shoe size



# **Pressure Reduction Insole**

- Insole features removable pegs that help off-load pressure from the plantar surface of the foot to provide optimal shock absorption and moldability; easily fits into walker or post-op shoe
- Indicated to help reduce pressure and promote healing after surgery, trauma or when wounds are present
- Constructed of 18mm Plastazote®, Poron® and multi-foam

Product Number

XS	S	M	L	XL
69-1060	69-1061	69-1062	69-1063	69-1064

Sizing Guide

XS	S	M	L	XL
Mens	up - 4 <sup>1</sup> / <sub>2</sub>	5 - 71/2	8 - 10	101/2+
Womens 4 - $5^{1}/_{2}$	6 - 81/2	9 - 11	$11^{1}/_{2}+$	

Measured by shoe size





# **Med-Surg Post-Op Shoe**


- Features added support under the metatarsal heads, square toe design for added comfort and protection, and rocker bottom to promote a natural gait
- Indicated for post-op use, when K-wires are present or after acute injuries
- Constructed of EVA insole and tricot foam upper and contact closure

### Product Number

	S	M	L	XL	Pedi
Mens	69-1002	69-1003	69-100	69-1005	
Womens	69-1007	69-1008	69-1009		
Child					69-1001

### Sizing Guide

	S	M	L	XL	Pedi
Mens	6-8	81/2-10	101/2-12	121/2-14	
Womens	4-6	61/2-8	81/2-10		
Child					12-1





## Slimline Cast Boot

- Features square toe design for optimal comfort and protection, forefoot closure adjusts to accommodate for casts or bandages and can be switched between left/right closure; ankle strap provides added stability
- Indicated for post-surgical application or use with casts or compression bandaging of the foot and ankle
- Constructed of EVA insole with tricot foam upper and contact closure

### Product Number

	XS	S	M	L	XL
Adult	69-1090	69-1091	69-1092	69-1093	69-1094

### Sizing Guide

	XS	S	M	L	XL
Mens			6-8	81/2-11	11½-14
Womens	up to 5	51/2-8	81/2-10	101/2-13	131/2+

Measured by shoe size See p. 89 Pediatric Slimline Cast Boot



# **Soft Heel Cushions**

- Features soft center comfort dot to help provide effective relief from heel pain; sold by pair
- Indicated for treatment of heel spurs, leg pain, general foot fatigue and for treatment of plantar fasciitis or Achilles tendonitis
- Manufactured from medical grade silicone with soft center comfort dot

### Product Number

	S	M	L/XL
Pkg/PR	99-4301-100	99-4302-100	99-4303-100

### Sizing Guide

	MENS SIZE	WOMENS SIZE
Small	4 - 6	4 - 71/2
Medium	61/2 - 10	8 - 11
Large/XL	$10^{1}/_{2}$ +	111/2 +
Measured	by shoe size	



# **Target Ankle Wrap**

- Features wrap-around design to provide warmth, compression and protection with Velstretch figure-8 strap for compression over malleolus
- Indicated for treatment of ankle strains or sprains and athletic protection
- Manufactured from premium grade latex-free neoprene

### Product Number

TTOGGCTT	Valifica	
	S/M	L/XL
	88-4010	88-4011

### Sizing Guide

S/M	L/XL
10"- 13"	13"- 16"

Measurement is circumference from base of heel around ankle

# **Bunion Splint**

- Designed to provide comfortable immobilization and alignment of the hallux
- Indicated for post-operative alignment following hallux valgus surgery
- Constructed of <sup>1</sup>/<sub>8</sub>" soft foam laminated to flannel lining with movable stay and contact closure

Product Number

		S	M	L
R	ight	61-5237	61-5238	61-5239
Le	eft	61-5337	61-5338	61-5339

Sizing Guide

S	M	L
9"-10"	10"-11"	11"-12"

Measurement taken at widest part of foot



# Cryo Pneumatic Ankle Wrap

- Anatomical ankle wrap designed to combine the benefits of pneumatic compression with cold therapy; features pump bulb that detaches once desired compression is achieved
- Indicated for post-op application, severe ankle strains / sprains, as part of a physical therapy or athletic training program; helps diminish swelling / edema, decrease pain and reduce muscle spasms
- Manufactured from plush foam laminate that is Velcro compatible; removable gel pack with tricot lining; non-toxic gel



### Product Number

	0 GELS	1 GEL	2 GELS
	72-7210	72-7211	72-7212
Ankle Gel	72-7219		

Sizing Guide

### UNIV

Fits up to 16" ankle circumference  $\,$ 



# **Cryotherm Ankle Wrap**

- Economical wrap designed to provide compression and concentrated cold / hot therapy for the ankle
- Indicated for post-op use, arthritis, strains or sprains and as part of rehab program
- Manufactured from a soft breathable foam laminate with contact closure and pockets to hold reusable gel packs

Product Number

0 GELS 2 GELS 4 GELS

72-7140 72-7142 72-7144

Sizing Guide

### UNIV

Fits up to 14" ankle circumference

Additional gels sold separately: 72-7003 - Single 3 <sup>1</sup>/<sub>2</sub>"x 11" Gel 72-7002 - Case/10 3 <sup>1</sup>/<sub>2</sub>"x 11" Gel

Measurement is circumference from base of heel around ankle







# **Next Step Suspension Belt**

- Designed to provide trans-femoral amputee optimal comfort, suspension and rotational control
- Ideal form of primary suspension, ancillary rotational control for lanyard and liner wearers
- Manufactured from latex-free soft-grip neoprene against the skin with durable nylon outer layer and heavy-duty 2" nylon straps with contact closure

Product Number

24"	26"	28"	30"	32"
89-1450	89-1451	89-1452	89-1453	89-1454
34"	36"	38"		

Sizing Guide

24"	26"	28"	30"	32"
2"x 24"	2"x 26"	2"x 28"	2"x 30"	2"x 32"
34"	36"	38"		
2"x 34"	2"x 36"	2"x 38"		

Measurement taken from anteriorly-mounted chafe across opposite hip to posterior chafe, subtract 6" from measurement



# **Hip Abduction Pillow**

- Features concave shape to ensure comfortable positioning and anatomical fit; vacuum packed for easy storage
- Indicated for post-operative hip positioning
- Manufactured from sturdy mediumdensity foam with soft foam straps and contact closure

Product Number

S	M	L
89-1411	89-1412	89-1413

Sizing Guide (Length x Width x Height)

0121119	Garac (Length)	· · · · · · · · · · · · · · · · · · ·	9110
	S	M	L
	17"×12"× 7"	22"x15"x7"	25"x16"x7"





# Inflatable Heel Boot

- Features two adjustable, pre-inflated air chambers to elevate the heel and help protect the Achilles tendon; venting holes promote healthy skin; vinyl exterior slides easily to help reduce sheer forces during bed mobility and soft interior helps protect fragile skin and improves patient comfort
- Indicated for treatment and prevention of pressure sores up to Stage IV
- Manufactured from medical-grade vinyl with contact closure

Product Number

UNIV

64-0101

Sizing Guide

UNI

Fits up to 19" calf circumference



# **Ultra Bucks Traction Splint**

- Features medial and lateral stays for added support and convoluted foam to help prevent migration
- Ideal for post-op treatment of hip fractures, hip dislocation and femur fractures
- Constructed of soft convoluted foam and aluminum spreader bar with contact closure

Product Number

S	М	L	XL
54-2161	54-2162	54-2163	54-2164

Sizing Guide

5	IVI	L	XL
10"- 12"	12"- 14"	14"- 16"	16"- 18"

Measurement taken at widest part of calf



# **Economizer Bucks Traction Splint**

- Universal design features medial and lateral stays for added support
- Ideal for post-op treatment of hip fractures, hip dislocation and femur fractures
- Constructed of breathable foam laminate with contact closure

Product Number

**UNIV** 54-2165

Sizing Guide

UNIV

12"- 18"

Measurement taken at widest part of calf



# **Universal Bucks Traction Splint**

- Universal design features medial and lateral stays for added support and convoluted foam to help prevent migration
- Ideal for post-op treatment of hip fractures, hip dislocation and femur fractures
- Constructed of soft convoluted foam with polyester pile leg and heel rest and aluminum spreader bar with contact closure

Product Number

UNIV

54-2150


Sizing Guide

UNIV

12"- 18"

Measurement taken at widest part of calf





# Lower Extremity

# Ŭ

# Philadelphia Cervical Collar w/Trachea

- Features two-piece design anatomically contoured to provide superior stability and comfort with large trachea opening for easy access
- Indicated for immobilization for both extrication and rehabilitation, post-operative discectomy and spinal fusion
- Constructed of durable plastic laminated to soft hypoallergenic foam with contact closure



Product Number

**PEDI** 2 1/4" 19-4100

Sizing Guide

PEDI

8"- 11"

Measurement taken from circumference of the neck at "Adam's apple"

Height measurement taken from posterior portion of mandible to the clavicle directly below



# **Ambu Perfit Collar**

- Features one-piece flat design with flip out chin piece and trachea opening for easy access
- Indicated for immobilization for extrication and rehabilitation
- Constructed of durable plastic laminated to foam with contact closure



Product Number

INFANT	PEDI	MINI ADJ
19-9001	19-9002	19-9008

Sizing Guide

INFANT	PEDI	MINI ADJ
Purple	Aqua	Univ Pediatric



## **Ultra Cervical Collar**

- Features anatomically contoured foam with soft cotton stockinette covering for maximum patient comfort
- Indicated for cervical facet syndrome, discogenic pain, degenerated disc disease and muscle sprains or strains of cervical region
- Manufactured from medium firm density foam with 2" contact closure



Product Number

	285	XS
2"	15-5000	15-5001

Sizing Guide

**2XS XS** 10"-12" 14"-16"

Measurement taken from circumference of the neck at "Adam's apple"

Height measurement taken from posterior portion of mandible to the clavicle directly below





This symbol indicates that pediatric sizes are available for this product.



This symbol indicates that bariatric sizes are available for this product.



This symbol certifies that all products manufactured by Corflex are latex free.



This symbol indicates that certain precautions are needed when using this product.



This symbol identifies CoolTex products that are manufactured from a neoprene-alternative material.



# **Ultra Shoulder Immobilizer**

- Features 1½" wide waist strap for secure immobilization and thumb loop to help prevent wrist drop
- Indicated for immobilization following rotator cuff repairs, S/P glenohumeral acromioplasty surgery or other shoulder-related injuries
- Manufactured from seamless, lightweight, durable material with added pad for shoulder strap; crayons included

### Product Number

2XS	XS	S	М
23-1869	23-1870	23-1871	23-1872

### Sizing Guide

	2XS	XS	S	M
Width/ Height	10"x 4"	12"x 5"	15"x 6"	17"x 7"

Measurements are sling pouch width x height



# **Economizer Shoulder Immobilizer**

- Features 1" shoulder and waist straps with slide buckle
- Indicated for immobilization following rotator cuff repairs, S/P glenohumeral acromioplasty surgery or other shoulder-related injuries
- Manufactured from seamless lightweight muslin material; crayons included



XS	S	M
23-1740	23-1741	23-1742

### Sizing Guide

	XS	S	M
Width/ Height	12"x 5"	15"x 6"	17"x 7"

Measurements are sling pouch width x height



### Tricot Shoulder Immobilizer

- Features waist strap for secure fit; thumb loop helps prevent wrist drop;
   Velcro closure at top of sling envelope
- Indicated for immobilization following rotator cuff repairs, S/P glenohumeral acromioplasty surgery or other shoulder-related injuries
- Manufactured from soft, durable tricot material

# Product Number

2XS	XS	S	M
23-5059	23-5060	23-5061	23-5062

### Sizing Guide

	2XS	XS	S	M
Width/ Height	10"x 4"	12"x 5"	15"x 6"	17"x 7"

Measurements are sling pouch width x height



# **Ultra Clavicle Strap**

- Features heavy-duty felt center point pad with hook and loop closure for optimal comfort
- Indicated for clavicle fractures, Grade II/III acromiclavicular separations, sternoclavicular dislocations and postural control
- Constructed of breathable, open cell stockinette covered foam with contact closure

Product Number

2XS	XS	S
25-1729	25-1730	25-1731

### Sizing Guide

Jizing Guiac			
	2XS	XS	S
	15"- 21"	17"- 24"	25"- 31"

Measurement taken from center back, over shoulder under arm and back to center point of back



# **Ultra Arm Sling**

- Features double d-ring hook and loop closure for ease of application and thumb loop to help prevent wrist drop
- Indicated for injuries to the arm, shoulder, wrist or hand and cast support to maintain arm in neutral position
- Manufactured from seamless, lightweight, durable material with added shoulder strap pad; crayons included



XS	S	M
21-1870	21-1871	21-1872

Sizing Guide

_			
	XS		
Width/ Height	12"x 5"	15"x 6"	17"x 7

Measurements are sling pouch width x height



# **Economizer Arm Sling**

- Allows for slide buckle or contact closure for ease of application
- Indicated for injuries to the arm, shoulder, wrist or hand and cast support to maintain arm in neutral position
- Manufactured from seamless lightweight muslin material; crayons included



### Product Number

2XS	XS	S	M
21-1969	21-1970	21-1971	21-1972

### Sizing Guide

	2XS	XS	S	M
Width/ Height	10"x 4"	12"x 5"	15"x 6"	17"x 7"

Measurements are sling pouch width x height



# **Tricot Arm Sling**

- Features thumb loop to prevent wrist drop and Velcro closure at top of sling envelope to provide a secure fit
- Indicated for injuries to the arm, shoulder, wrist or hand and cast support to maintain arm in neutral position
- Manufactured from soft yet durable tricot material



### Product Number

2XS	XS	S	M
21-5059	21-5060	21-5061	21-5062

# Sizing Guide

	2XS	XS	S	M
Width/ Height	10"x 4"	12"x 5"	15"x 6"	17"x 7"

Measurements are sling pouch width x height



# Mesh Arm Sling

- Allows for slide buckle or contact closure
- Indicated for injuries to the arm, shoulder, wrist or hand and cast support to maintain arm in neutral position
- Manufactured from seamless breathable mesh material



### Product Number

XS	S	M
21-1730	21-1731	21-1732

# Sizing Guide

	XS	S	M
Width/ Height	12"x 5"	15"x 6"	17"x 7"

Measurements are sling pouch width x height





# Signature and Ultra Vinyl Wrist Splint

- Features v-style elastic loop locks with ultra thin web strap to provide secure immobilization of wrist while preserving finger mobility
- Indicated for Carpal Tunnel Syndrome, tendonitis, wrist strains or sprains and post-cast removal
- Constructed of premium grade perforated vinyl laminated to soft flannel lining with malleable aluminum palmar and dorsal stays and contact closure

Product Number

Signature	2XS	XS	
5" Right	73-5540	73-5542	
5" Left	73-5541	73-5543	
Ultra	2XS	XS	
<b>Ultra</b> 6" Right	<b>2XS</b> 41-1510	<b>XS</b> 41-1513	

Sizina Guide

Jizirig Guide			
	2XS	XS	
	3"_ 1"	4"- 5"	

Measurement taken from wrist circumference



# Signature Foam Flannel Wrist Splint

- Provides support of the wrist while preserving finger mobility
- Indicated for Carpal Tunnel Syndrome, tendonitis, wrist strains or sprains and post-cast removal
- Constructed of cotton with flannel lining; malleable aluminum palmar stay, flexible palmar stay and contact closure

Product Number

	4XS	3XS
5" Right	73-5687	73-5688
5" Left	73-5697	73-5698
6" Right		73-5718
6" Left		73-5728

Sizina Guide

5		
	4XS	3XS
	21/2"- 31/2"	3"- 4"

Measurement taken from wrist circumference



# Target Wrist Stabilizer w/Stay

- Easy to apply wrap-around design provides support, warmth, compression and protection for the wrist
- Indicated for Carpal Tunnel Syndrome, tendonitis and wrist strains or sprains
- Manufactured from latex-free soft loop neoprene with removable malleable palmar stay

Product Number

	UNIV K	OINIV L
Pedi	88-2056	88-2057

- J	
	UNIV
Pedi	Fits up to 7" wrist circumference



# **Thumbster Soft**

- Anatomically designed to provide comfortable and secure immobilization of the thumb and wrist support while preserving mobility of the other fingers
- Indicated for deQuervain's Syndrome, gamekeeper's thumb, tendonitis, thumb strains or sprains and post-cast removal
- Constructed of vacuum-formed cross linked polyester with stretchable contact closure

Product Number

UNIV 7 ½" 73-5011

Sizing Guide		
	UNIV	
Pedi	Fits most children & young adolescents	

# **Compression Knee Immobilizer**

- Features sewn-in contoured posterior stays, movable medial/lateral stays and elastic compression straps for optimal immobilization
- Indicated for post-op immobilization and mild to severe knee injuries where full extension of the leg is desired
- Constructed of breathable foam laminate, aluminum stays, premium quality elastic straps and contact closure

Product Number

LENGTH	6"	9"	101/2"	12"
Univ	75-0651	75-0951	75-1051	75-1251

Sizing Guide

	6"	9"	101/2"	12"
Univ	11"	14"	14"	24"

Measurement taken from leg circumference 4" above mid-patella



# **Universal Knee Immobilizer**

- Universal design allows for inventory reduction; features contoured posterior stays and movable medial/ lateral stay assemblies for optimal positioning and immobilization with 4" elastic patella strap for added compression
- Indicated for post-op knee immobilization and mild to severe knee injuries where full extension of the leg is desired
- Constructed of breathable foam laminate to soft loop pile with hook and loop closure straps, premium quality elastic straps and aluminum stays



# Product Number

LENGTH	6"	9"	101/2"	12"
Univ	51-1506	51-1509	51-5510	51-1512

## Sizing Guide

	6"	9"	101/2"	12"
Univ	11"	14"	15 <sup>1</sup> / <sub>2</sub> "	20"

Measurement taken from leg circumference 4" above mid-patella



# Ultra Tricot Knee Immobilizer

- Features sewn-in contoured posterior stays and medial/lateral stays with patella opening for optimal immobilization
- Indicated for post-op knee immobilization and mild to severe knee injuries where full extension of the leg is desired
- Constructed of breathable foam with canvas outer, hook and loop closure straps and aluminum stays

Product Number

LENGTH	2XS	XS	S	М
13"	51-1529	51-1530	51-1531	51-1532

# Sizing Guide

2XS	XS	S	M
7"- 10"	10"- 13"	13"- 16"	16"- 18"

Measurement taken from leg circumference 4" above mid-patella







# **Open Toe Cast Sandal**

- Features open toe and heel for comfort and breathability
- Indicated for lower extremity casted fractures
- Constructed of non-skid, ribbed rocker sole with durable upper and contact closure

Product Numbe
---------------

	3XS	2XS	XS	S
	69-3078	69-3079	69-3080	69-3081
Sizing Guic	de			
	3XS	2XS	XS	S
Length	4"	5"	7"	9"





# **Enclosed Toe Cast Boot**

- Features enclosed toe to provide complete protection and warmth for the foot
- Indicated for lower extremity casted fractures
- Constructed of ribbed rocker sole with water-resistant vinyl upper and contact closure

Product	Number

XS	S	M
69-3100	69-3101	69-3102

# Sizing Guide

	XS	S	М
Length	7"	9"	11"





# **Slimline Cast Boot**

- Forefoot closure adjusts to accommodate for casts or bandages; ankle strap provides added stability
- Indicated for post-surgical application or use with casts or compression bandaging of the foot and ankle
- Constructed of EVA insole with neoprene upper; contact closure

# Product Number

S	M	L
69-1101	69-1102	69-1103

### Sizing Guide

	S	М	L
Length	5-71/2	8-11½	12-11/2





# **Pediatric Fixed Walker**

- Easy to apply design provides stability for the lower leg with low natural heel height; features a removable washable foam liner and aluminum uprights for adjustability; fits children from 1½ to 9 years
- Indicated for stable foot or ankle fractures, severe sprains or strains, post cast removable or as a comfortable alternative to casting
- Constructed of covered malleable aluminum uprights, padded liner and contact closure straps

### Product Number

S/M

	79-1201	79-1202	79-1203
Sizing Guide			
	S/M	L	XL



# **Med-Surg Post-Op Shoe**

- Features added support under the metatarsal heads, square toe design for added comfort and protection; rocker bottom to promote a natural gait
- Indicated for post-op uses, when K-wires are present or after acute injuries
- Constructed of EVA insole and tricot foam upper and contact closure

### Product Number

	PEDI	S
Pedi Mens Womens	69-100	69-1002 69-1007

# Sizing Guide

	PEDI	S
Pedi Mens Womens	12-1	6-8 4-6
Measured by	/ shoe size	4-0





# **Rigid Post-Op Shoe**

- Features flat footbed with rigid sole to prevent excessive foot motion while walking
- Indicated for post-operative foot protection and support
- Constructed of wood-like sole with water-resistant upper and contact closure

### Product Number

	PEDI	S
Pedi Mens Womens	69-3121	69-3122 69-3132

### Sizing Guide

	PEDI	S
Length Pedi Length Mens Length Womens	8"	11" 9"





# Game On Ankle Stabilizer

- Anatomically designed to support and protect the ankle; utilizes a dual figure-8 strapping system to simulate a professional taping technique
- Indicated for treatment of moderate ankle sprains, strains, ankle instability and weakness, or for prophylactic use
- Constructed of durable ballistic nylon, ventilated spacer fabric, nylon figure-8 straps and wide elastic compression strap

### Product Number

2XS	XS	S
75-8020	75-8021	75-8022

# Sizing Guide

oizing carao			
	2XS	XS	S
	9"-10"	10"-11"	11"-12"

Measurement is circumference from base of heel around. Fits either left or right foot.





# **Knee-O-Band**

- Features solid tubular neoprene pad that applies pressure to the patella tendon
- Indicated for Osgood Schlatter disease, patello-femoral syndrome and patella tendonitis
- Manufactured from premium grade latex-free soft loop neoprene

### Product Number

	3XS	2X	XS	S
3/16"	88-7210	88-7211	88-7212	88-7213

# Sizing Guide

3XS	2X	XS	S
10"-11"	11"-12"	12"- 13"	13"- 14"

Measurement is circumference from base of heel around. Fits either left or right foot.







Α		BINDERS/BELTS	
All Purpose Boot	77	Elastic Rib Belt	10
Aluminum Finger Strip		Panel Elastic Abdominal Binder	11
Ambo Elastic Wrist Splint		Trochanter Belt	
Ambu Perfit Collar		Ultra Abdominal Binder	11
Ankle Fixed Walker		Universal Abdominal Binder	11
		Universal Rib Belt	10
ANKLE IMMOBILIZERS		Boxer Splint	45
Ankle Fixed Walker		Bunion Splint	79
Lower Leg Fixed Walker		С	
Pneumatic R.O.M. Walker			_
Pneumatic Walker		Calf Sleeve	69
R.O.M. Walker		CALF/THIGH SLEEVES	
Tri Shell Pneumatic Walker		Calf Sleeve	69
Ankle Stirrup	/5	Thigh Sleeve	69
ANKLE SUPPORTS		Canvas Wrist Lacer Splint	35
Ankle Stirrup	75	CDS Knee Orthosis	49
Game-On Ankle Stabilizer	75	CERVICAL SUPPORTS	
Lace-Up Ankle	75	Ambu Perfit Collar	12
Anterior Closure Knee Wrap w/Hinge	54	Cryotherm Neck Wrap	
Anterior Closure Knee Wrap w/R.O.M. Hinge	54	Medic-Air Cervical Sleep Pillō	
Anterior Closure Knee Wrap w/Stays	68	Medic-Air Cervical Sleep Fillo	
Anterior Closure Knee Wrap	68	Philadelphia Cervical Collar	
ARM SLINGS		Ultra Cervical Collar	
Collar & Cuff	25		
Cradle Arm Sling		Classic Post-Op Shoe Cock-up Splint	
Economizer Arm Sling		·	
•		Collar & Cuff Sling	
Hanging Cast Sling		Colles Splint	
Mesh Arm Sling		Compression Knee Immobilizer	
Strap Arm Sling		Contender Knee Brace	
Tricot Arm Sling		Contour Back Support	
Ultra Arm Sling	24	CoolTex Anterior Closure Knee Wrap w/Hinge	
В		CoolTex Anterior Closure Knee Wrap w/R.O.M. Hinge	
Back-2-Form Support	5	CoolTex Hinged Knee Sleeve w/Lo Pro Hinge	
		CoolTex J-Pull Patella Stabilizer w/Lo Pro Hinge	
BACK MISCELLANEOUS	0	CoolTex Soft AFO	
Medic-Air Back Pillō		C-Pull Patella Stabilizer	
Medic-Air Lumbar Röll		Cradle Arm Sling	
Medic-Air Seat Cushion	9	Criss-Cross Back Support	
BACK SUPPORTS		Cryo Pneumatic Ankle Wrap	
Back-2-Form Support.	5	Cryo Pneumatic Elbow Wrap	
Contour Back Support	5	Cryo Pneumatic Knee Orthosis w/Hinge	
Criss-Cross Back Support	6	Cryo Pneumatic Knee Orthosis w/ROM Hinge	
Cryotherm Back Wrap	8	Cryo Pneumatic Knee Splint	
E/N Lumbar Support	5	Cryo Pneumatic Shoulder	
Industrial Back Support	8	Cryo Pneumatic Wrist Support	
Low Profile Industrial Back Support	8	Cryotherm Ankle Wrap	
Target Back Wrap	6	Cryotherm Back Wrap	
Thoracolumbar Support	7	Cryotherm Elbow Wrap	
Ultra Lumbo Sacral Support		Cryotherm Knee Wrap	
V-Lock Back Support		Cryotherm Neck Wrap	
V-Lock Lumbar Orthosis w/Rigid Insert		Cryotherm Shoulder Wrap	
Band It Elbow Strap		Cryotherm Vari Wrap	
Baseball Splint		Cryotherm Wrist Wrap	
		Curved Finger Splint	42

D	
Digital Splint	42
Disc Unloader LO	3
Disc Unloader LSO	4
E	
	г
E/N Lumbar Support	
Economizer Arm Sling.	
Economizer Bucks Traction Splint	
Economizer Shoulder Immobilizer	
Elastic Rib Belt	
Elastic Shoulder Immobilizer	
Elastic Wrist Splint	
Elbow Orthosis w/R.O.M. Hinge	21
ELBOW SLEEVES	
Target Elbow Sleeve	
Target Elbow Sleeve w/Pad	
Target Elbow Sleeve w/Strap	28
ELBOW SPLINTS	
Cryo Pneumatic Elbow Wrap	30
Cryotherm Elbow Wrap	30
Elbow Orthosis w/R.O.M. Hinge	27
Padded Elbow Wrap	27
Ultra Clavicle Strap	26
Universal Elbow Immobilizer	27
Ventral Cubital Tunnel Splint	26
XR Cubital Tunnel Splint	26
Enclosed Toe Cast Boot	76
ER Shoulder Abduction Pillow w/Sling	17
Extended Length Boxer Splint	45
Extended Length Humeral Splint	23
F	
Finger Cot Splint	42
Finger Protector	
FINGER SPLINTS	
	10
Aluminum Finger Strip	
Baseball Splint	
Curved Finger Splint	
Digital Splint	
Finger Cot Splint	
Finger Protector	
Frog Splint	
Stax Splint	42
FOOT/ANKLE MISCELLANEOUS	
All Purpose Boot	
Bunion Splint	
Classic Post-Op Shoe	76
Cryo Pneumatic Ankle Wrap	79
Cryotherm Ankle Wrap	79 79
Cryotherm Ankle Wrap  Enclosed Toe Cast Boot	79 79 76
Cryotherm Ankle Wrap  Enclosed Toe Cast Boot  Heelwedge Shoe	79 79 76 77
Cryotherm Ankle Wrap  Enclosed Toe Cast Boot	79 79 76 77 78

Orthowedge Shoe	. 77
Pressure Reduction Insole	. 77
Rigid Post-Op Shoe	
Slimline Cast Boot	
Soft Heel Cushions	
Target Ankle Wrap	
Forearm Splint	
Frog Splint	. 42
G	
Game-On Ankle Stabilizer	. 75
н	
HAND/FOREARM SPLINTS	
Boxer Splint	. 45
Cock-up Splint	. 43
Colles Splint	. 43
Extended Length Boxer Splint	. 45
Forearm Splint	. 43
Metacarpal Hand Splint	. 43
Metacarpal Splint	. 43
Target MP Ulnar Drift Support	. 44
Ulnar Splint	. 44
Hanging Cast Sling	. 25
Heelwedge Shoe	. 77
HINGED KNEE SUPPORTS	
Anterior Closure Knee Wrap w/Hinge	. 54
Anterior Closure Knee Wrap w/R.O.M. Hinge	. 54
Contender Knee Brace	. 52
CoolTex Anterior Closure Knee Wrap w/Hinge	. 55
CoolTex Anterior Closure Knee Wrap w/R.O.M. Hinge	. 55
CoolTex Hinged Knee Sleeve w/Lo Pro Hinge	
CoolTex J-Pull Patella Stabilizer w/Lo Pro Hinge	. 59
Kinetic Posterior Adj Knee Sleeve w/Cor-Trak Buttress & Hinge	.59
Knee Sleeve w/Hinge	. 56
Knee Sleeve w/R.O.M. Hinge	
Knee-O-Trakker w/Hinge	. 58
Knee-O-Trakker w/Hinge 6"	
Patella Stabilizer w/Cor-Trak Buttress & Hinge	. 58
Posterior Adjustable Knee Sleeve w/Cor-Trak Buttress & Hinge	
Posterior Adjustable Knee Sleeve w/Hinge	
Posterior Adjustable Knee Sleeve w/R.O.M. Hinge	
U-Shaped Patella Stabilizer w/Hinge	
Hip Abduction Pillow	. 80
HIP MISCELLANEOUS	
Hip Abduction Pillow	. 80
Next Step Suspension Belt	. 80
HIP TRACTION	
Economizer Bucks Traction Splint	
Inflatable Heel Boot	
Ultra Bucks Traction Splint	
Universal Bucks Traction Splint	
Humeral Splint	. 23

HUMERAL SPLINTS  Extended Length Humeral Splint	23
Humeral Splint Mid Length Humeral Splint	
T I	
Industrial Back Support Industrial Wrist Splint Inflatable Heel Boot	35
J	
J-Pull Patella Stabilizer	63
К	
Kinetic Posterior Adj Knee Sleeve w/Cor-Trak Buttress Kinetic Posterior Adj Knee Sleeve w/Cor-Trak Buttress & Hinge	
KNEE IMMOBILIZERS	
Compression Knee Immobilizer	
Tri Panel Knee Immobilizer	
Ultra Tricot Knee Immobilizer Universal Knee Immobilizer	
	40
KNEE MISCELLANEOUS	<b>-</b> 4
Cryo Pneumatic Knee Orthosis w/Hinge Cryo Pneumatic Knee Orthosis w/ROM Hinge	
Cryo Pneumatic Knee Orthosis W/KOW Filinge Cryo Pneumatic Knee Splint	
Cryotherm Knee Wrap	
Cryotherm Vari Wrap	
KNEE ORTHOSIS	
CDS Knee Orthosis	49
Knee Sleeve	
Knee Sleeve w/Anterior Pad	
Knee Sleeve w/Hinge	56
Knee Sleeve w/R.O.M. Hinge	56
Knee Sleeve w/Stays	66
KNEE SLEEVES	
Knee Sleeve	67
Knee Sleeve w/Anterior Pad	67
Knee Sleeve w/Stays	66
KNEE WRAPS	
Anterior Closure Knee Wrap w/Stays	68
Anterior Closure Wrap	68
Universal Knee Wrap	
Universal Knee Wrap w/Dynamic Buttress	
Kneed It	
Knee-O-Band	
Knee-O-Trakker 6"  Knee-O-Trakker w/Hinge	
Knee-O-Trakker w/Hinge 6"	
Knee-O-Trakker w/Stays	
Knee-O-Trakker w/Tilt Strap 6"	

L	
Lace-Up Ankle Low Profile Industrial Back Support Low Profile Thumbster Lower Leg Fixed Walker	8 41
M	
Med-Surg Post-Op Shoer Medial-Lateral Patella Stabilizer Medic-Air Back Pillō Medic-Air Cervical Sleep Pillō Medic-Air Lumbar Rōll Medic-Air Seat Cushion Medic-Air Snooze Pillō Mesh Arm Sling Metacarpal Hand Splint Metacarpal Splint Mid Length Humeral Splint	64 9 12 9 13 24 43
N	∠∠
Neutral Shoulder Pillow Next Step Suspension Belt	
0	
Open Toe Cast Sandal Orthowedge Shoe Osgood-Schlatter Knee Sleeve	77
Р	
Padded Elbow Wrap Padded Plantar Fasciitis Night Splint Panel Elastic Abdominal Binder	73
PATELLA STABILIZER SUPPORTS	
C-Pull Patella Stabilizer J-Pull Patella Stabilizer Kinetic Posterior Adj Knee Sleeve w/Cor-Trak Buttress Kneed It Knee-O-Band Knee-O-Trakker 6" Knee-O-Trakker w/Stays Knee-O-Trakker w/Tilt Strap 6" Medial-Lateral Patella Stabilizer Osgood-Schlatter Knee Sleeve Patella Stabilizer w/Cor-Trak Buttress Posterior Adjustable Knee Sleeve w/Cor-Trak Buttress	63 64 65 65 60 61 60 62
U-Shaped Patella Stabilizer	61
Patella Stabilizer w/Cor-Trak Buttress	
Patella Stabilizer w/Cor-Trak Buttress & Hinge Pediatric Fixed Walker	

# PEDIATRIC PRODUCTS Ambu Perfit Collar 83 Compression Knee Immobilizer 87 Economizer Arm Sling 85 Economizer Shoulder Immobilizer 84 Enclosed Toe Cast Boot 88 Game-On Ankle Stabilizer 89 Knee-O-Band 89 Med-Surg Post-Op Shoe 89 Mesh Arm Sling 85 Open Toe Cast Sandal 88 Pediatric Fixed Walker 88 Philadelphia Cervical Collar 83 Rigid Post-Op Shoe 89 Signature Vinyl Wrist Splint 86 Slimline Cast Boot 88 Thumbster Soft 86 Tricot Arm Sling 85 Tricot Shoulder Immobilizer 84 Ultra Arm Sling 85 Ultra Cervical Collar 83 Ultra Clavicle Strap 84 Ultra Shoulder Immobilizer 84 Ultra Vinyl Wrist Splint 86 Universal Knee Immobilizer 87 Philadelphia Cervical Collar 12 **PLANTAR FASCIITIS SPLINTS** Soft AFO 74 Ultra Dorsal Night Splint 74 Ultra Stretch Plantar Fasciitis Splint 73 Pneumatic R.O.M. Walker 72 Pneumatic Walker 70 Posterior Adjustable Knee Sleeve w/Cor-Trak Buttress & Hinge .. 53 Posterior Adjustable Knee Sleeve w/Hinge \_\_\_\_\_53 Pressure Reduction Insole 77 R R.O.M. Walker 72 S Shoulder Abduction Pillow w/Sling \_\_\_\_\_\_\_16

Economizer Shoulder Immobilizer Elastic Shoulder Immobilizer ER Shoulder Abduction Pillow w/Sling Neutral Shoulder Pillow w/Sling Shoulder Abduction Pillow w/Firm Fit Sling Shoulder Abduction Pillow w/Sling Sling & Swathe Ultra Shoulder Abduction Pillow w/Sling Ultra Shoulder Abduction Pillow w/Firm Fit Sling Ultra Shoulder Immobilizer Ultra Sling & Swathe Universal Comfort Shoulder Immobilizer Universal Shoulder Immobilizer Universal Sling & Swathe SHOULDER MISCELLANEOUS Cryo Pneumatic Shoulder Cryotherm Shoulder Wrap	18 17 15 16 19 16 16 20 19 18 20 18
Tricot Shoulder Immobilizer Signature Foam Flannel Wrist Splint Signature Vinyl Wrist Splint Slimline Cast Boot Sling & Swathe Soft Heel Cushions	35 34 78
Disc Unloader LO Disc Unloader LSO Stax Splint Strap Arm Sling Suede Wrist Lacer Splint Suede Wrist Lacer Splint w/Abducted Thumb	4 42 25 34
Т	
Target Ankle Wrap  Target Back Wrap  Target Contoured Tennis Elbow Wrap  Target Elbow Sleeve  Target Elbow Sleeve w/Pad  Target Elbow Sleeve w/Strap  Target Hitchhiker Thumb  Target Plantar Fasciitis Soft Boot  Target Sewn Thumb  Target Tennis Elbow Sleeve w/Pad  Target Universal Tennis Elbow Wrap  Target Wrist Splint  Target Wrist Stabilizer  Target Wrist Stabilizer  Target Wrist Stabilizer  Target Wrist Stabilizer	6 29 28 28 41 74 41 29 29 37 37
Farget MP Ulnar Drift Support	44

SHOULDER IMMOBILIZERS

TENNIS ELBOW SUPPORTS	
Band It Elbow Strap	30
Target Contoured Tennis Elbow Wrap	29
Target Tennis Elbow Sleeve w/Pad	29
Target Universal Tennis Elbow Wrap	
Thigh Sleeve	
Thoracolumbar Support	
THUMB SPLINTS	
Low Profile Thumbster	41
Suede Wrist Lacer Splint w/Abducted Thumb	39
Target Hitchhiker Thumb	
Target Sewn Thumb	
Thumbster Soft	
Ultra Fit Cool Wrist Splint w/Abducted Thumb	
Ultra Fit Wrist Splint w/Abducted Thumb	
Wrist/Hand/Thumb Orthosis	
Thumbster Soft	
Tri Panel Knee Immobilizer	
Tri Shell Pneumatic Walker	
Tricot Arm Sling	
Tricot Shoulder Immobilizer	
Trochanter Belt	
	10
U	
Ulnar Splint	44
•	
Ultra Abdominal Binder	.11
Ultra Abdominal Binder Ultra Arm Sling	.11 .24
Ultra Abdominal Binder Ultra Arm Sling Ultra Bucks Traction Splint	. 11 . 24 . 81
Ultra Abdominal Binder Ultra Arm Sling Ultra Bucks Traction Splint Ultra Cervical Collar	.11 .24 .81 .12
Ultra Abdominal Binder Ultra Arm Sling Ultra Bucks Traction Splint Ultra Cervical Collar Ultra Clavicle Strap	. 11 . 24 . 81 . 12 . 26
Ultra Abdominal Binder Ultra Arm Sling Ultra Bucks Traction Splint Ultra Cervical Collar Ultra Clavicle Strap Ultra Dorsal Night Splint	. 11 . 24 . 81 . 12 . 26
Ultra Abdominal Binder Ultra Arm Sling Ultra Bucks Traction Splint Ultra Cervical Collar Ultra Clavicle Strap Ultra Dorsal Night Splint Ultra Fit Cool Wrist Splint	.11 .24 .81 .12 .26 .74
Ultra Abdominal Binder Ultra Arm Sling Ultra Bucks Traction Splint Ultra Cervical Collar Ultra Clavicle Strap Ultra Dorsal Night Splint	.11 .24 .81 .12 .26 .74
Ultra Abdominal Binder Ultra Arm Sling Ultra Bucks Traction Splint Ultra Cervical Collar Ultra Clavicle Strap Ultra Dorsal Night Splint Ultra Fit Cool Wrist Splint Ultra Fit Cool Wrist Splint w/Abducted Thumb Ultra Fit Wrist Splint	.11 .24 .81 .12 .26 .74 .33 .39
Ultra Abdominal Binder Ultra Arm Sling Ultra Bucks Traction Splint Ultra Cervical Collar Ultra Clavicle Strap Ultra Dorsal Night Splint Ultra Fit Cool Wrist Splint Ultra Fit Cool Wrist Splint w/Abducted Thumb	.11 .24 .81 .12 .26 .74 .33 .39
Ultra Abdominal Binder Ultra Arm Sling Ultra Bucks Traction Splint Ultra Cervical Collar Ultra Clavicle Strap Ultra Dorsal Night Splint Ultra Fit Cool Wrist Splint Ultra Fit Cool Wrist Splint w/Abducted Thumb Ultra Fit Wrist Splint	.11 .24 .81 .12 .26 .74 .33 .39 .32
Ultra Abdominal Binder Ultra Arm Sling Ultra Bucks Traction Splint Ultra Cervical Collar Ultra Clavicle Strap Ultra Dorsal Night Splint Ultra Fit Cool Wrist Splint Ultra Fit Wrist Splint Ultra Fit Wrist Splint Ultra Fit Wrist Splint Ultra Fit Wrist Splint	.11 .24 .81 .12 .26 .74 .33 .39 .32 .38
Ultra Abdominal Binder Ultra Arm Sling Ultra Bucks Traction Splint Ultra Cervical Collar Ultra Clavicle Strap Ultra Dorsal Night Splint Ultra Fit Cool Wrist Splint Ultra Fit Wrist Splint w/Abducted Thumb Ultra Fit Wrist Splint w/Abducted Thumb Ultra Lumbo Sacral Support	.11 .24 .81 .12 .26 .74 .33 .39 .32 .38 6
Ultra Abdominal Binder Ultra Arm Sling Ultra Bucks Traction Splint Ultra Cervical Collar Ultra Clavicle Strap Ultra Dorsal Night Splint Ultra Fit Cool Wrist Splint Ultra Fit Wrist Splint w/Abducted Thumb Ultra Fit Wrist Splint w/Abducted Thumb Ultra Lumbo Sacral Support Ultra Plantar Fasciitis Night Splint Ultra Shoulder Abduction Pillow w/Sling Ultra Shoulder Abduction Pillow w/Firm Fit Sling	11 24 81 12 26 74 33 39 32 38 6 73 16
Ultra Abdominal Binder Ultra Arm Sling Ultra Bucks Traction Splint Ultra Cervical Collar Ultra Clavicle Strap Ultra Dorsal Night Splint Ultra Fit Cool Wrist Splint Ultra Fit Wrist Splint w/Abducted Thumb Ultra Fit Wrist Splint w/Abducted Thumb Ultra Lumbo Sacral Support Ultra Plantar Fasciitis Night Splint Ultra Shoulder Abduction Pillow w/Sling	11 24 81 12 26 74 33 39 32 38 6 73 16
Ultra Abdominal Binder Ultra Arm Sling Ultra Bucks Traction Splint Ultra Cervical Collar Ultra Clavicle Strap Ultra Dorsal Night Splint Ultra Fit Cool Wrist Splint Ultra Fit Wrist Splint w/Abducted Thumb Ultra Fit Wrist Splint w/Abducted Thumb Ultra Lumbo Sacral Support Ultra Plantar Fasciitis Night Splint Ultra Shoulder Abduction Pillow w/Sling Ultra Shoulder Abduction Pillow w/Firm Fit Sling	11 24 81 12 26 74 33 39 32 38 6 73 16 16 20
Ultra Abdominal Binder Ultra Arm Sling Ultra Bucks Traction Splint Ultra Cervical Collar Ultra Clavicle Strap Ultra Dorsal Night Splint Ultra Fit Cool Wrist Splint Ultra Fit Cool Wrist Splint w/Abducted Thumb Ultra Fit Wrist Splint Ultra Fit Wrist Splint w/Abducted Thumb Ultra Fit Wrist Splint w/Abducted Thumb Ultra Fit Wrist Splint w/Abducted Thumb Ultra Lumbo Sacral Support Ultra Plantar Fasciitis Night Splint Ultra Shoulder Abduction Pillow w/Sling Ultra Shoulder Immobilizer	11 24 81 12 26 74 33 39 32 38 6 73 16 16 20
Ultra Abdominal Binder Ultra Arm Sling Ultra Bucks Traction Splint Ultra Cervical Collar Ultra Clavicle Strap Ultra Dorsal Night Splint Ultra Fit Cool Wrist Splint Ultra Fit Cool Wrist Splint w/Abducted Thumb Ultra Fit Wrist Splint w/Abducted Thumb Ultra Fit Wrist Splint w/Abducted Thumb Ultra Lumbo Sacral Support Ultra Plantar Fasciitis Night Splint Ultra Shoulder Abduction Pillow w/Sling Ultra Shoulder Immobilizer Ultra Sling & Swathe	11 24 81 12 26 74 33 39 32 38 6 73 16 16 20 19 73
Ultra Abdominal Binder Ultra Arm Sling Ultra Bucks Traction Splint Ultra Cervical Collar Ultra Clavicle Strap Ultra Dorsal Night Splint Ultra Fit Cool Wrist Splint Ultra Fit Wrist Splint w/Abducted Thumb Ultra Fit Wrist Splint w/Abducted Thumb Ultra Lumbo Sacral Support Ultra Plantar Fasciitis Night Splint Ultra Shoulder Abduction Pillow w/Sling Ultra Shoulder Immobilizer Ultra Sling & Swathe Ultra Stretch Plantar Fasciitis Splint	11 24 81 12 26 74 33 39 32 38 6 73 16 16 20 19 73 48
Ultra Abdominal Binder Ultra Arm Sling Ultra Bucks Traction Splint Ultra Cervical Collar Ultra Clavicle Strap Ultra Dorsal Night Splint Ultra Fit Cool Wrist Splint Ultra Fit Wrist Splint w/Abducted Thumb Ultra Lumbo Sacral Support Ultra Plantar Fasciitis Night Splint Ultra Shoulder Abduction Pillow w/Sling Ultra Shoulder Abduction Pillow w/Firm Fit Sling Ultra Shoulder Immobilizer Ultra Sling & Swathe Ultra Tricot Knee Immobilizer	11 24 81 12 26 74 33 39 32 38 6 73 16 16 20 19 73 48 34
Ultra Abdominal Binder Ultra Arm Sling Ultra Bucks Traction Splint Ultra Cervical Collar Ultra Clavicle Strap Ultra Dorsal Night Splint Ultra Fit Cool Wrist Splint Ultra Fit Wrist Splint w/Abducted Thumb Ultra Lumbo Sacral Support Ultra Plantar Fasciitis Night Splint Ultra Shoulder Abduction Pillow w/Sling Ultra Shoulder Abduction Pillow w/Firm Fit Sling Ultra Shoulder Immobilizer Ultra Stretch Plantar Fasciitis Splint Ultra Tricot Knee Immobilizer Ultra Vinyl Wrist Splint	11 24 81 12 26 74 33 39 32 38 6 73 16 16 20 19 73 48 34
Ultra Abdominal Binder Ultra Arm Sling Ultra Bucks Traction Splint Ultra Cervical Collar Ultra Clavicle Strap Ultra Dorsal Night Splint Ultra Fit Cool Wrist Splint Ultra Fit Wrist Splint w/Abducted Thumb Ultra Fit Wrist Splint w/Abducted Thumb Ultra Lumbo Sacral Support Ultra Plantar Fasciitis Night Splint Ultra Shoulder Abduction Pillow w/Sling Ultra Shoulder Abduction Pillow w/Firm Fit Sling Ultra Shoulder Immobilizer Ultra Sling & Swathe Ultra Stretch Plantar Fasciitis Splint Ultra Tricot Knee Immobilizer Ultra Vinyl Wrist Splint Universal Abdominal Binder	11 24 81 12 26 74 33 39 32 38 6 73 16 16 20 19 73 48 34 11 81
Ultra Abdominal Binder Ultra Arm Sling Ultra Bucks Traction Splint Ultra Cervical Collar Ultra Clavicle Strap Ultra Dorsal Night Splint Ultra Fit Cool Wrist Splint Ultra Fit Wrist Splint w/Abducted Thumb Ultra Fit Wrist Splint w/Abducted Thumb Ultra Lumbo Sacral Support Ultra Plantar Fasciitis Night Splint Ultra Shoulder Abduction Pillow w/Sling Ultra Shoulder Abduction Pillow w/Firm Fit Sling Ultra Shoulder Immobilizer Ultra Sling & Swathe Ultra Stretch Plantar Fasciitis Splint Ultra Tricot Knee Immobilizer Ultra Vinyl Wrist Splint Universal Abdominal Binder Universal Bucks Traction Splint	11 24 81 12 26 74 33 39 32 38 6 73 16 16 20 19 73 48 34 11 81 18

Universal Knee Wrap	68
Universal Knee Wrap w/Dynamic Buttress	69
Universal Plantar Fasciitis Support	74
Universal Rib Belt	
Universal Shoulder Immobilizer	20
Universal Sling & Swathe	18
U-Shaped Patella Stabilizer	61
U-Shaped Patella Stabilizer w/Hinge	57
V	
Ventral Cubital Tunnel Splint	26
V-Lock Back Support	
V-Lock Lumbar Orthosis w/Rigid Insert	
W	
Wrist & Forearm Immobilizer	36
Wrist/Hand Orthosis	
Wrist/Hand/Thumb Orthosis	
	10
WRIST SPLINTS	٠,
Ambo Elastic Wrist Splint	
Canvas Wrist Lacer Splint	
Cryo Pneumatic Wrist Support	
Cryotherm Wrist Wrap  Elastic Wrist Splint	
•	
Industrial Wrist Splint Signature Foam Flannel Wrist Splint	
Signature Vinyl Wrist Splint	
Suede Wrist Lacer Splint	
Target Wrist Splint	
Target Wrist Stabilizer	
Target Wrist Stabilizer w/Stay	
Target Wrist Wrap	
Ultra Fit Cool Wrist Splint	
Ultra Fit Wrist Splint	
Ultra Vinyl Wrist Splint.	
Wrist & Forearm Immobilizer	
Wrist/Hand Orthosis	
X	
	01
XR Cubital Tunnel Splint	26

Back Support Comparison Chart							
Back Support	Anterior Panel Height	Posterior Panel Height	Page #				
Disc Unloader LO	6"	12"	3				
Disc Unloader LO w/10" Anterior Panel	10"	12"	3				
Disc Unloader LSO	6"	15"	4				
Disc Unloader LSO w/10" Anterior Panel	10"	15"	4				
E/N Lumbar Support	6"	9"	5				
Back-2-Form Support	6"	9"	5				
Contour Back Support	6"	9"	5				
Ultra Lumbo Sacral Support	6"	11"	6				
Criss-Cross Back Support	6"	10"	6				
Target Back Wrap	6"	10"	6				
Thorocolumbar Support	9"	20"	7				
V-Lock Back Support	6"	11"	7				
Industrial Back Support	6"	9"	8				
Low Profile Industrial Back Support	6"	6"	8				
Cryotherm Back Wrap	6"	12"	8				

Metric Conversion Chart							
IN	CM	IN	CM	IN	CM		
1	2.5	11	27.9	21	53.3		
2	5.1	12	30.5	22	55.9		
3	7.6	13	33.0	23	58.4		
4	10.2	14	35.6	24	61.0		
5	12.7	15	38.1	25	63.5		
6	15.2	16	40.6	26	66.0		
7	17.8	17	43.2	27	68.6		
8	20.3	18	45.7	28	71.1		
9	22.9	19	48.3	29	73.7		
10	25.4	20	50.8	30	76.2		

Formula - INCH  $\times$  2.54 = CM

# **CUSTOMER SERVICE**

800.426.7353, Ext.1 or 603.623.3344, Ext.1 or service@corflex.com







Designed for life.

669 East Industrial Park Drive Manchester, NH 03109 USA

800.426.7353

**CORFLEX.COM**