



CORFLEX®

Designed for life.

800 426 7353

Pediatric Solutions



Corflex offers a comprehensive line of pediatric orthopedic rehabilitation products.

CORFLEX.COM



Pediatric Fixed Walker

- > Easy to apply design provides stability for the lower leg with low natural heel height; features a removable washable foam liner and aluminum uprights for adjustability; fits children from 1 ½ to 9 years
- > Indicated for stable foot or ankle fractures, severe sprains or strains, post cast removable or as a comfortable alternative to casting
- > Constructed of covered malleable aluminum uprights, padded liner and contact closure straps

Product Number

	S/M	L	XL
	79-1201	79-1202	79-1203

Sizing Guide

Age	1-2 years	2 ½-6 years	6 ½-9 ½ years
Foot Length	up to 6"	6"-7 ½"	7 ½"-8 ½"

Suggested Billing Code L-4386

All products manufactured by Corflex are latex free.



Please feel free to contact us for our full catalog or visit our website at www.corflex.com



Ultra Shoulder Immobilizer

- > Features 1½" wide waist strap for secure immobilization and thumb loop to help prevent wrist drop
- > Indicated for immobilization following rotator cuff repairs, S/P glenohumeral acromioplasty surgery or other shoulder-related injuries
- > Manufactured from seamless, lightweight, durable material with added pad for shoulder strap; crayons included

Product Number

2XS	XS	S	M
23-1869	23-1870	23-1871	23-1872

Sizing Guide

Width/Height	10" x 4"	12" x 5"	15" x 6"	17" x 7"
--------------	----------	----------	----------	----------

Measurements are sling pouch width x height

Suggested Billing Code L-3670



Economizer Shoulder Immobilizer

- > Features 1" shoulder and waist straps with slide buckle
- > Indicated for immobilization following rotator cuff repairs, S/P glenohumeral acromioplasty surgery or other shoulder-related injuries
- > Manufactured from seamless lightweight muslin material; crayons included

Product Number

XS	S	M
23-1740	23-1741	23-1742

Sizing Guide

Width/Height	12" x 5"	15" x 6"	17" x 7"
--------------	----------	----------	----------

Measurements are sling pouch width x height

Suggested Billing Code L-3670

2



Tricot Shoulder Immobilizer

- > Features waist strap for secure fit; thumb loop helps prevent wrist drop; Velcro closure at top of sling envelope
- > Indicated for immobilization following rotator cuff repairs, S/P glenohumeral acromioplasty surgery or other shoulder-related injuries
- > Manufactured from soft, durable tricot material

Product Number

2XS	XS	S	M
23-5059	23-5060	23-5061	23-5062

Sizing Guide

Width/Height	10" x 4"	12" x 5"	15" x 6"	17" x 7"
--------------	----------	----------	----------	----------

Measurements are sling pouch width x height

Suggested Billing Code L-3670



Ultra Clavicle Strap

- > Features heavy-duty felt center point pad with hook and loop closure for optimal comfort
- > Indicated for clavicle fractures, Grade II/III acromioclavicular separations, sternoclavicular dislocations and postural control
- > Constructed of breathable, open cell stockinette covered foam with contact closure

Product Number

2XS	XS	S
25-1729	25-1730	25-1731

Sizing Guide

15" - 21"	17" - 24"	25" - 31"
-----------	-----------	-----------

Measurement taken from center back, over shoulder under arm and back to center point of back

Suggested Billing Code L-3650



Ultra Arm Sling

- > Features double d-ring hook and loop closure for ease of application and thumb loop to help prevent wrist drop
- > Indicated for injuries to the arm, shoulder, wrist or hand and cast support to maintain arm in neutral position
- > Manufactured from seamless, lightweight, durable material with added shoulder strap pad; crayons included

Product Number

XS	S	M
21-1870	21-1871	21-1872

Sizing Guide

Width/Height	12"x 5"	15"x 6"	17"x 7"
--------------	---------	---------	---------

Measurements are sling pouch width x height



Economizer Arm Sling

- > Allows for slide buckle or contact closure for ease of application
- > Indicated for injuries to the arm, shoulder, wrist or hand and cast support to maintain arm in neutral position
- > Manufactured from seamless lightweight muslin material; crayons included

Product Number

2XS	XS	S	M
21-1969	21-1970	21-1971	21-1972

Sizing Guide

Width/Height	10"x 4"	12"x 5"	15"x 6"	17"x 7"
--------------	---------	---------	---------	---------

Measurements are sling pouch width x height



Tricot Arm Sling

- > Features thumb loop to prevent wrist drop and Velcro closure at top of sling envelope to provide a secure fit
- > Indicated for injuries to the arm, shoulder, wrist or hand and cast support to maintain arm in neutral position
- > Manufactured from soft yet durable tricot material

Product Number

2XS	XS	S	M
21-5059	21-5060	21-5061	21-5062

Sizing Guide

Width/Height	10"x 4"	12"x 5"	15"x 6"	17"x 7"
--------------	---------	---------	---------	---------

Measurements are sling pouch width x height



Mesh Arm Sling

- > Allows for slide buckle or contact closure
- > Indicated for injuries to the arm, shoulder, wrist or hand and cast support to maintain arm in neutral position
- > Manufactured from seamless breathable mesh material

Product Number

XS	S	M
21-1730	21-1731	21-1732

Sizing Guide

Width/Height	12"x 5"	15"x 6"	17"x 7"
--------------	---------	---------	---------

Measurements are sling pouch width x height



Signature and Ultra Vinyl Wrist Splint

- > Features v-style elastic loop locks with ultra thin web strap to provide secure immobilization of wrist while preserving finger mobility
- > Indicated for Carpal Tunnel Syndrome, tendonitis, wrist strains or sprains and post-cast removal
- > Constructed of premium grade perforated vinyl laminated to soft flannel lining with malleable aluminum palmar and dorsal stays and contact closure

Product Number

Signature	2XS	XS
5" Right	73-5540	73-5542
5" Left	73-5541	73-5543
Ultra		
6" Right	41-1510	41-1513
6" Left	41-1511	41-1514

Sizing Guide

3" - 4"	4" - 5"
---------	---------

Measurement taken from wrist circumference

Suggested Billing Code L-3908



Signature Foam Flannel Wrist Splint

- > Provides support of the wrist while preserving finger mobility
- > Indicated for Carpal Tunnel Syndrome, tendonitis, wrist strains or sprains and post-cast removal
- > Constructed of cotton with flannel lining; malleable aluminum palmar stay, flexible palmar stay and contact closure

Product Number

	4XS	3XS
5" Right	73-5687	73-5688
5" Left	73-5697	73-5698
6" Right		73-5718
6" Left		73-5728

Sizing Guide

2 1/2" - 3 1/2"	3" - 4"
-----------------	---------

Measurement taken from wrist circumference

Suggested Billing Code L-3908



Target Wrist Stabilizer w/Stay

- > Easy to apply wrap-around design provides support, warmth, compression and protection for the wrist
- > Indicated for Carpal Tunnel Syndrome, tendonitis and wrist strains or sprains
- > Manufactured from latex-free soft loop neoprene with removable malleable palmar stay

Product Number

UNIV R	UNIV L
88-2056	88-2057

Sizing Guide

Fits up to 7" wrist circumference

Suggested Billing Code L-3908



Thumbster Soft

- > Anatomically designed to provide comfortable and secure immobilization of the thumb and wrist while preserving mobility of the other fingers
- > Indicated for deQuervain's Syndrome, gamekeeper's thumb, tendonitis, thumb strains or sprains and post-cast removal
- > Constructed of vacuum-formed cross linked polyester with stretchable contact closure

Product Number

UNIV
7 1/2" L 73-5011

Sizing Guide

Fits most children & young adolescents

Suggested Billing Code L-3807



Philadelphia Cervical Collar w/Trachea

- > Features two-piece design anatomically contoured to provide superior stability and comfort with large trachea opening for easy access
- > Indicated for immobilization for both extrication and rehabilitation, post-operative discectomy and spinal fusion
- > Constructed of durable plastic laminated to soft hypoallergenic foam with contact closure

Product Number

	PEDI
2 1/4" H	19-4100

Sizing Guide

8" - 11"

Measurement taken from circumference of the neck at "Adam's apple". Height measurement taken from posterior portion of mandible to the clavicle directly below.

Suggested Billing Code L-0172



Ambu Perfit Collar

- > Features one-piece flat design with flip out chin piece and trachea opening for easy access
- > Indicated for immobilization for extrication and rehabilitation
- > Constructed of durable plastic laminated to foam with contact closure

Product Number

INFANT	PEDI	MINI ADJ
19-9001	19-9002	19-9008

Sizing Guide

Purple	Aqua	Univ Pediatric
--------	------	----------------



Mini Adjustable

Suggested Billing Code L-0172



Ultra Cervical Collar

- > Features anatomically contoured foam with soft cotton stockinette covering for maximum patient comfort
- > Indicated for cervical facet syndrome, discogenic pain, degenerated disc disease and muscle sprains or strains of cervical region
- > Manufactured from medium firm density foam with contact closure

Product Number

	2XS	XS	S
2" H	15-5000	15-5001	
2 1/2" H			15-5021

Sizing Guide

10"-12"	14"-16"	16"-20"
---------	---------	---------

Measurement taken from circumference of the neck at "Adam's apple". Height measurement taken from posterior portion of mandible to the clavicle directly below.

Suggested Billing Code L-0120



Open Toe Cast Sandal

- > Features open toe and heel for comfort and breathability
- > Indicated for lower extremity casted fractures
- > Constructed of non-skid, ribbed rocker sole with durable upper and contact closure

Product Number

	3XS	2XS	XS	S
	69-3078	69-3079	69-3080	69-3081

Sizing Guide

Length	4"	5"	7"	9"
--------	----	----	----	----

Suggested Billing Code L-3260



Enclosed Toe Cast Boot

- > Features enclosed toe and heel to provide complete protection and warmth for the foot
- > Indicated for lower extremity casted fractures
- > Constructed of non-skid, ribbed rocker sole with water-resistant vinyl upper and contact closure

Product Number

	XS	S	M
	69-3100	69-3101	69-3102

Sizing Guide

Length	7"	9"	11"
--------	----	----	-----

Suggested Billing Code L-3260

6



Slimline Cast Boot

- > Forefoot closure adjusts to accommodate for casts or bandages; ankle strap provides added stability
- > Indicated for post-surgical application or use with casts or compression bandaging of the foot and ankle
- > Constructed of EVA insole with neoprene upper; contact closure

Product Number

	S	M	L
	69-1101	69-1102	69-1103

Sizing Guide

	5-7½	8-11½	12-1½
--	------	-------	-------

Measured by shoe size

Suggested Billing Code L-3260



Pediatric Fixed Walker

- > Easy to apply design provides stability for the lower leg with low natural heel height; features a removable washable foam liner and aluminum uprights for adjustability; fits children from 1 ½ to 9 years
- > Indicated for stable foot or ankle fractures, severe sprains or strains, post cast removable or as a comfortable alternative to casting
- > Constructed of covered malleable aluminum uprights, padded liner and contact closure straps

Product Number

	S/M	L	XL
	79-1201	79-1202	79-1203

Sizing Guide

Age	1-2 years	2 ½-6 years	6 ½-9 ½ years
Foot Length	up to 6"	6"-7 ½"	7 ½"-8 ½"

Suggested Billing Code L-4386



Med-Surg Post-Op Shoe

- > Features added support under the metatarsal heads, square toe design for added comfort and protection; rocker bottom to promote a natural gait
- > Indicated for post-op uses, when K-wires are present or after acute injuries
- > Constructed of EVA insole and tricot foam upper and contact closure

Product Number

	PEDI	S
Pedi	69-1001	
Mens		69-1002
Womens		69-1007

Sizing Guide

Pedi	12-1	
Mens		6-8
Womens		4-6

Measured by shoe size

Suggested Billing Code L-3260



Rigid Post-Op Shoe

- > Features flat footbed with rigid sole to prevent excessive foot motion while walking
- > Indicated for post-operative foot protection and support
- > Constructed of wood-like sole with water-resistant upper and contact closure

Product Number

	PEDI	S
Pedi	69-3121	
Mens		69-3122
Womens		69-3132

Sizing Guide

Length Pedi	8"	
Length Mens		11"
Length Womens		9"

Suggested Billing Code L-3260

7



Game On Ankle Stabilizer

- > Anatomically designed to support and protect the ankle; utilizes a dual figure-8 strapping system to simulate a professional taping technique
- > Indicated for treatment of moderate ankle sprains, strains, ankle instability and weakness, or for prophylactic use
- > Constructed of durable ballistic nylon, ventilated spacer fabric, nylon figure-8 straps and wide elastic compression strap

Product Number

	2XS	XS	S
	75-8020	75-8021	75-8022

Sizing Guide

	9"-10"	10"-11"	11"-12"
--	--------	---------	---------

Measurement is circumference from base of heel around. Fits either left or right foot.

Suggested Billing Code L-1902



Knee-O-Band

- > Features solid tubular neoprene pad that applies pressure to the patella tendon
- > Indicated for Osgood Schlatter disease, patello-femoral syndrome and patella tendonitis
- > Manufactured from premium grade latex-free soft loop neoprene

Product Number

	3XS	2XS	XS	S
3/16"	88-7210	88-7211	88-7212	88-7213

Sizing Guide

	10"-11"	11"-12"	12"-13"	13"-14"
--	---------	---------	---------	---------

Measurement taken at mid-patella



Compression Knee Immobilizer

- > Features sewn-in contoured posterior stays, movable medial/lateral stays and elastic compression straps for optimal immobilization
- > Indicated for post-op immobilization and mild to severe knee injuries where full extension of the leg is desired
- > Constructed of breathable foam laminate, aluminum stays, premium quality elastic straps and contact closure

Product Number

LENGTH	6"	9"	10½"	12"
	75-0651	75-0951	75-1051	75-1251

Sizing Guide

11"	14"	14"	24"
-----	-----	-----	-----

Measurement taken from leg circumference 4" above mid-patella

Suggested Billing Code L-1830



Universal Knee Immobilizer

- > Universal design allows for inventory reduction; features contoured posterior stays and movable medial/lateral stay assemblies for optimal positioning and immobilization with 4" elastic patella strap for added compression
- > Indicated for post-op knee immobilization and mild to severe knee injuries where full extension of the leg is desired
- > Constructed of breathable foam laminate to soft loop pile with hook and loop closure straps, premium quality elastic straps and aluminum stays

Product Number

LENGTH	6"	9"	10½"	12"
	51-1506	51-1509	51-5510	51-1512

Sizing Guide

11"	14"	15½"	20"
-----	-----	------	-----

Measurement taken from leg circumference 4" above mid-patella

Suggested Billing Code L-1830



Ultra Tricot Knee Immobilizer

- > Features sewn-in contoured posterior stays and medial/lateral stays with patella opening for optimal immobilization
- > Indicated for post-op knee immobilization and mild to severe knee injuries where full extension of the leg is desired
- > Constructed of breathable foam with canvas outer, hook and loop closure straps and aluminum stays

Product Number

	2XS	XS	S	M
13" L	51-1529	51-1530	51-1531	51-1532

Sizing Guide

7"- 10"	10"- 13"	13"- 16"	16"- 18"
---------	----------	----------	----------

Measurement taken from leg circumference 4" above mid-patella

Suggested Billing Code L-1830



CORFLEX®
Designed for life.

To place your order, contact our Customer Service Team

PHONE 800 426 7353, EXT 1

FAX 603 623 4111

EMAIL SERVICE@CORFLEX.COM

669 EAST INDUSTRIAL PARK DRIVE / MANCHESTER, NH 03109

CORFLEX.COM